

RTI Group v1.2 Ocean Next User Guide

Home » RTI Group » RTI Group v1.2 Ocean Next User Guide 🖺

Contents

- 1 RTI Group v1.2 Ocean
- Next
- 2 Product Usage
- Instructions
- 3 FAQs
- 4 Documents / Resources
 - **4.1 References**



RTI Group v1.2 Ocean Next



Specifications

• Product Name: RTI Scatter Probe

Release Date: 2024-09-09Firmware Version: v1.2

• Dll Version: v1.1

• Compatibility: Requires Ocean Next v1.1.0.0 or higher

Supports HW version 1.0 with added buzzer sound amplifier

Product Usage Instructions

Compatibility

This product is compatible with Ocean Next v1.1.0.0 or higher. It does not work with older versions of Ocean 2014 or other RTI software tools.

Installation

- 1. Ensure you have the required Ocean Next version installed on your system.
- 2. Connect the RTI Scatter Probe to your device using the appropriate cables.
- 3. Power on the device and the probe.
- 4. Follow the on-screen instructions to complete the installation process.

Operation

- 1. Launch the Ocean Next software.
- 2. Select the RTI Scatter Probe from the available devices list.
- 3. Follow the software prompts to configure and use the probe for your specific application.
- 4. Refer to the user manual for detailed instructions on utilizing the probe's features.

World-leading software for X-ray QA & testing

- With Ocean Next[™] software you receive all QA data including waveforms directly from your RTI meter or probe. You can perform anything from simple measurements to advanced analyses of X-ray equipment performance, save and retrieve your measurement values, and generate reports for documentation and traceability. Check out the release highlights and download the latest version here!
- When you register your free myRTI account, you get access to the our customer portal (meter info & logs), plus
 access to RTI Support and online training. With myRTI you can track your devices used with Ocean Next™, get
 reminders on calibration due dates, and access your calibration certificates.
- Remember, the software is compliant with all our meters Mako, Piranha and Cobia as well as the RTI Scatter Probe.

Ocean Next requires MS Windows 10 or later.

Plug n Play

Start measuring within seconds! Ocean Next™ detects what instrument and probes you have connected to assist you in the best way possible – just Plug n Play.

Customization

Plan your tests in advance, create checklists, add pop-up windows with user instructions, to simplify the work for you and your co-workers using a streamlined user interface.



Ocean Next™ Quick

Ocean Next™ Quick detects what instrument and probes you have connected to assist you in the best way possible – just Plug n Play. The interface adapts, and all the measured parameters are displayed on one screen,

including waveforms. The results can be retained in the database for later review and compiled in an autogenerated report.

Ocean Next™ Advantage

For more advanced, streamlined X-ray testing and QA measurements, Ocean Next™ Advantage enables the customization of templates that exactly meet your needs. In Studio View, you can design single-page templates, including analysis and checklists, and the result is also shown in reports with your own logo and layout.

Ocean Next™ Professional

A complete X-ray testing and QA system for efficiency and compliance. Ocean Next™ Professional offers trend analysis and full traceability of your measurements, including location and scheduling. Build a holistic solution for your X-ray QA, by using multi-page templates in a searchable database structured for your preventive maintenance.

| Compare Ocean Next™ functions | Quick | Advantage | Professional |
|---|-------|-----------|--------------|
| myRTI/myBox (Ocean Next 4.0) | | | |
| Built-in display function | | | |
| Save function | | | |
| Printing | | | |
| Basic Excel® connection | | | |
| Advanced Excel® connection | | | |
| User-defined tests and checklists | | | |
| Analysis with pass/fail limits | | | |
| Reporting | | | |
| CT Dose Profiler support | | | |
| Single-page templates | | | |
| Multi-page templates | | | |
| Site database (including location, room, equipment, etc.) | | | |
| Scheduling | | | , |

Best test equipment I've used in Biomed/Imaging

- "I was looking into the template more last night and realized this is the best test equipment I've used since I started working in Biomed/Imaging. Being able to create a PM checklist, record measurements and create a report from one device is amazing.
- · With how organized I like to be this will help me tremendously. "

Product Note

- RTI Group Headquarters Phone: +46 (0) 31 746 36 00
- Flöjelbergsgatan 8C info@rtigroup.com Company registration number: 556230-2462
- SE-431 37 Mölndal, SWEDEN www.rtigroup.com VAT number: SE556230246201 1(1)
- RTI Scatter Probe
- News in Ocean Next Release v4.2.1.0
- · Release information
- Release date 2024-09-09
- Release name Ocean Next v4.2.1.0
- Firmware version v1.2
- DII Version v1.1

Note: This product does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. Ocean Next v1.1.0.0 or higher is required. This release includes various bugfixes and support for HW version 1.0 (added buzzer sound amplifier).

Known Limitations

None.

Senior Imaging Engineer, Midwest Hospital Group

RTI Group Headquarters

Flöjelbergsgatan 8C SE-431 37 Mölndal SWEDEN

• E-mail: <u>sales@rtigroup.com</u>

• Phone: +46 31 746 36 00

RTI Group North America

33 Jacksonville Road, Bldg. 1 Towaco, NJ 07082 USA

• E-mail: sales.us@rtigroup.com

• Phone: +1 800-222-7537 (Sales)

RTI Group APAC 3791 Jalan Bukit Merah #08-10 E-Centre @ Redhill SINGAPORE 159471

• E-mail: sales@rtigroup.com

• Phone: +65 8608 2075 (Sales)

RTI Group Europe Ltd.

- 1a Kingsbury's Lane Ringwood BH24 1EL UNITED KINGDOM
- E-mail: sales@rtigroup.com

FAQs

Q: What versions of Ocean are compatible with the RTI Scatter Probe?

A: The RTI Scatter Probe is compatible with Ocean Next v1.1.0.0 or higher. It does not work with older versions of

Ocean 2014 or other RTI software tools.

Q: Is there a known limitation for the RTI Scatter Probe?

A: No, there are currently no known limitations for the RTI Scatter Probe.

Documents / Resources



References

- ♠ RTI Group forty years of leading X-ray QA solutions
- User Manual

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.