

starlab Polycarbonate Erlenmeyer Flasks User Guide

Contents

[1 Polycarbonate Erlenmeyer Flasks. Autoclaving instructions.](#)

[1.1 STEAM STERILISATION INSTRUCTIONS:](#)

[1.1.1 Filled Flasks](#)

[1.1.2 Empty](#)

[2 Documents / Resources](#)

[2.1 References](#)

[3 Related Posts](#)

Polycarbonate Erlenmeyer Flasks. Autoclaving instructions.



starlab.click/erlenmeyers



STEAM STERILISATION INSTRUCTIONS:

Starlab`s Erlenmeyer Flasks are designed to withstand at least 10 steam sterilisation cycles. Recommended sterilisation temperature is **121°C at 15 psi, up to 60 minutes.**

Filled Flasks

Loosen the cap so it's barely engaged, and sterilise according to your lab's usual guidelines or protocol. Before securing the cap, allow the flasks and contents to cool to at least 55 °C. (Tightening the cap on a hot flask will distort the flask as the contents cool down.)

Empty


Flasks Remove the caps from the flasks and sterilise in autoclave bags.

As with glass Erlenmeyers, cover the neck and tops of the flasks with aluminum foil. Crimp the foil enough so it does not fall off during autoclaving.

Use a small piece of autoclave tape to secure the foil to the flask, and to indicate that the flask has been exposed to steam.

Product Data Sheet.

Documents / Resources

	starlab Polycarbonate Erlenmeyer Flasks [pdf] User Guide Polycarbonate Erlenmeyer Flasks, Erlenmeyer Flasks, Flasks
---	--

References

- [User Manual](#)

Manuals+. [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.