



TANDEM t:connect Mobile App User Guide


[Home](#) » [Tandem](#) » TANDEM t:connect Mobile App User Guide 

Contents

- 1 TANDEM t:connect Mobile App
- 2 Product Information
- 3 Product Usage Instructions
- 4 Two Simple Steps
- 5 What You Need
- 6 Install the Uploader
- 7 Software and Drivers
- 8 Create an Account
- 9 Important Safety Information
- 10 PRECAUTION
- 11 Documents / Resources
 - 11.1 References



TANDEM t:connect Mobile App



Product Information

The t:connect web application is a tool designed to help manage diabetes. It allows users to view data from their insulin pump, generate reports, and share them with healthcare providers. The application provides colorful reports displaying blood glucose levels, insulin boluses, insulin on board (IOB), and more, making it easier for users to spot blood glucose trends and make informed decisions. The t:connect web application user guide provides comprehensive information on how to use the application. It can be accessed at tconnect.tandemdiabetes.com/help.

Product Usage Instructions

1. Install Software & Drivers:

1. Visit tconnect.tandemdiabetes.com.
2. Click on “GET STARTED” to download the t:connect uploader.
3. Follow the online steps to install the uploader on your computer.

2. Create an Account:


1. Run the t:connect uploader software. Note: The software should automatically start after installation. If not, double click the program’s icon on your desktop.
2. Connect the t:slim X2 pump to your computer using the Tandem USB cable. Note: You do not need to disconnect your infusion set or stop insulin delivery.
3. Start the upload process by clicking the “Start Upload” button on the t:connect uploader. The uploader will detect your pump and display its serial number on the computer screen.

If you need further assistance, visit support.tandemdiabetes.com.

Welcome to a positively different way to actively manage your diabetes.


The t:connect web application allows you to easily view data from your insulin pump and generate reports that you can share with your healthcare providers, helping you uncover meaningful insights and trends. This guide provides the basic steps to help you get started. For more comprehensive information, please refer to the t:connect web application user guide available at tconnect.tandemdiabetes.com/help

Two Simple Steps



1. Install Software & Drivers
2. Create Account

What You Need



Suggested Items




Visit tandemdiabetes.com/tconnect for a list of compatible browsers and operating systems.
Need help? We are here for you.
Visit support.tandemdiabetes.com

Install the Uploader

Welcome to t:connect

The t:connect™ Diabetes Management Application is a simple tool to help manage your diabetes. Colorful reports displaying blood glucose levels, insulin boluses, insulin on board (IOB), and more make it easy to help you spot blood glucose trends and make informed decisions.




- 1 To download the t:connect uploader, click **GET STARTED** on tconnect.tandemdiabetes.com

tconnect.tandemdiabetes.com

Note: The software download experience may vary depending on the operating system and web browser you are using.

Software and Drivers



- 2 Click on the green button to download the t:connect uploader, and follow the online steps to install it on your computer.



Create an Account

1. Run the t:connect uploader software

Note: The uploader software should automatically start after you complete the installation. If it does not

automatically start, double click the program's icon on your desktop.

2. Connect the pump to your computer

Connect the micro-USB end of the Tandem USB cable to your insulin pump and the other end to a USB port on your computer.



Note: You do not need to disconnect your infusion set or stop insulin delivery

3. Start upload

The t:connect uploader will detect your pump and display its serial number on the computer screen. Click Start Upload button.

4. Get started

The uploader will prompt you to create a new account. Click the Get Started button

5. Select account type

☐ **Personal**

☐ **Parent, Guardian**

Note: If you already have an account, just sign in with your username (email address) and password to add your new pump to your existing t:connect account.

6. Fill in the necessary details to create your t:connect account*

My First Name	My Last Name
<input type="text"/>	<input type="text"/>
Email	Confirm Email
<input type="text"/>	<input type="text"/>
My Password (8 characters minimum)	Confirm Password
<input type="text"/>	<input type="text"/>

7. Define how often your cartridge and infusion set should be changed*

<input checked="" type="radio"/> Target Glucose Range (mg/dL)	<input type="text" value="70"/> - <input type="text" value="180"/>
<input checked="" type="radio"/> Low Glucose Limit (mg/dL)	<input type="text" value="69"/>
<hr/>	
My pump cartridge should be changed:	<input type="text"/>
My infusion set site should be changed:	<input type="text"/>

Congratulations! You have created your account.

Note: Check your email inbox for a welcome message from Tandem Diabetes Care, and click the link to confirm your email address

* Interface may vary slightly based on your operating system.



Now that you have a t:connect account, download the free t:connect mobile app!


Give yourself the freedom to more easily manage your diabetes with the t:slim X2™ insulin pump and t:connect mobile app:



- Bolus from your smartphone without touching your pump†
- Use your phone as a secondary display for your insulin pump
- Wirelessly upload your t:slim X2 pump and therapy data.‡



Note: We recommend using the USB cable that came with your pump to download your pump data to your computer before using the t:connect mobile app. This will expedite wireless uploads. Visit support.tandemdiabetes.com to see a list of supported devices and frequently asked questions.



† Bolus delivery from the t:connect mobile app requires a compatible smartphone model and operating system, an app update, a remote software update on the customer's t:slim X2 insulin pump and additional training. Only available to customers who reside in the United States and who are in warranty at the time they update their pump. ‡ Uploads to the t:connect web application do not take place in real time and should not be relied upon by healthcare providers or caregivers for remote patient monitoring. Standard carrier data rates may apply. Smartphone sold separately.

Important Safety Information

RX ONLY. The t:slim X2 pump and Control-IQ™ technology are intended for single patient use. The t:slim X2 pump and Control-IQ technology are indicated for use with NovoLog or Humalog U-100 insulin. t:slim X2 insulin pump: The t:slim X2 insulin pump with interoperable technology is an alternate controller enabled (ACE) pump that is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in people requiring insulin. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The t:slim X2 pump is indicated for use in individuals six years of age and greater. Control-IQ technology: Control-IQ technology is intended for use with a compatible integrated continuous glucose monitor (iCGM, sold separately) and ACE pump to automatically increase, decrease, and suspend delivery of basal insulin based on iCGM readings and predicted glucose values. It can also deliver correction boluses when the glucose value is predicted to exceed a predefined threshold. Control-IQ technology is intended for the management of Type 1 diabetes mellitus in persons six years of age and greater.

WARNING: Control-IQ technology should not be used by anyone under the age of six years old. It should also not be used in patients who require less than 10 units of insulin per day or who weigh less than 55 pounds.

Control-IQ technology is not indicated for use in pregnant women, people on dialysis, or critically ill patients. Do not use Control-IQ technology if using hydroxyurea. Users of the t:slim X2 pump and Control-IQ technology must: use the insulin pump, CGM, and all other system components in accordance with their respective instructions for use; test blood glucose levels as recommended by their healthcare provider; demonstrate adequate carb-counting skills; maintain sufficient diabetes self-care skills; see healthcare provider(s) regularly; and have adequate vision and/or hearing to recognize all functions of the pump, including alerts, alarms, and reminders. The t:slim X2

pump, and the CGM transmitter and sensor must be removed before MRI, CT, or diathermy treatment. Visit tandemdiabetes.com/safetyinfo for additional important safety information. t:connect mobile app: The feature set available within the t:connect mobile app is dependent on the pump software version and the compatible smartphone's model and operating system (OS).

- The Display and Data Upload feature set provides a secondary display of pump and continuous glucose monitoring (CGM) information, including display of your pump alerts and alarms, and enables wireless uploading of pump and CGM data to the Tandem cloud through an internet or wireless data connection. Standard carrier data rates may apply.
- The Bolus Delivery plus Display and Data Upload feature set additionally allows users to request, stop, and cancel a bolus on the pump from the t:connect mobile app. At least one smartphone security setting must be enabled to utilize the Bolus Delivery feature of the t:connect mobile app.
- **WARNING:** Always rely on your pump to make therapy decisions when using a smartphone that is incompatible with the Bolus Delivery feature.


PRECAUTION

Always rely on your pump to make therapy decisions if using a smartphone that is incompatible, the smartphone is lost or damaged, or the smartphone loses Bluetooth® Connectivity with your pump. Important pump alerts and alarms are only received as app notifications when enabled and the app is either active or running in the background.

© 2022 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, Control-IQ, t:connect, and t:slim X2 are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. The Bluetooth® wordmark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Tandem Diabetes Care, Inc. is under license. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google Play and the Google Play logo are trademarks of Google LLC. All third-party marks are the property of their respective owners

11075 Roselle Street • San Diego, CA 92121
tandemdiabetes.com

Documents / Resources

	<p>TANDEM t:connect Mobile App [pdf] User Guide t connect Mobile App, t connect, Mobile App, App</p>
---	---

References

- [Min.in](#)
- [Diabetes Management | Insulin Pumps | Tandem Diabetes Care](#)
- [Important Safety Information - Tandem Insulin Pumps](#)

- [t:connect Web Application | Tandem Diabetes Care](#)
- [Tandem Diabetes Care, Inc | t:connect](#)
- [Tandem Diabetes Care, Inc | t:connect](#)

Manuals+.