



# TOMLOV DM201 SE LCD Digital Microscope User Manual

[Home](#) » [TOMLOV](#) » TOMLOV DM201 SE LCD Digital Microscope User Manual 

# TOMLOV

LCD Digital Microscope  
Model: DM201 SE



### Contents

- 1 DM201 SE LCD Digital Microscope
- 2 Read Me First
- 3 Packing List
- 4 Installation
- 5 Quick Start
- 6 Parts Description
- 7 Detailed Operation
- 8 Menu Settings
- 9 Specifications
- 10 FAQ
- 11 Get Support
- 12 Documents / Resources
  - 12.1 References
- 13 Related Posts

## DM201 SE LCD Digital Microscope

Please read this user manual carefully and keep it for future reference.  
Should you encounter any problem, please contact us by following e-mail.  
Most technical questions can be answered through our knowledge base.  
Customer service: [support@tomlov.com](mailto:support@tomlov.com)

### Read Me First

1. The product contains tiny parts and is inedible, please avoid accidental ingestion by children. Also, this product is not intended for medical use.
2. Clean the protective UV filter under the lens if you notice a deterioration in image quality. Carefully wipe the microscope with soft and clean cotton cloth (such as an eyeglass cloth) moistened with alcohol.
3. Monitor on the main unit is not shockproof, please use it with care. Do not drop it, press it hard or touch it with sharp object. The monitor is not waterproof, please clean it with dry and soft cloth. Do not press the monitor when adjusting angle.
4. TOMLOV offers one year manufacture warranty and life-long technical support for digital microscope. Please feel free to get in touch with us if you need any help. It is recommended to attach photos or videos showing the problems in the e-mail/message for quick solutions. TOMLOV customer service usually responds in 12 hours on business day.



E-mail: [support@tomlov.com](mailto:support@tomlov.com)



Facebook: @Tomlov




WhatsApp: +86 18123716479



Website: [www.tomlov.com](http://www.tomlov.com)

### Packing List



Please note: Accessories maybe updated without prior notice, please refer to the actual receipt of the items.

## Installation

Please make sure the following position A B C D are installed tightly to avoid shaking.

### 1. Position A

It has been adjusted properly upon delivery. If not, please adjust it by the wrench.



### 2. Position B

Step 1: Rotate the black ring clockwise one or two turns.



Step 2: Screw the rod clockwise.



Step 3: Screw the black ring clockwise to fix the bracket tightly.



### 3. **Position C**

Insert the monitor and tighten the two little screws.



### 4. **Position D**


Adjust the monitor up and down, then tighten this screw.





## Quick Start

<p>1. Connect the USB-C slot on the base and 5V output slot on the monitor by the USB cable to provide power for the side lights.</p>	<p>2. Insert Micro SD card.</p>	<p>3.</p>
<p>4. Adjust to get suitable distance and clear image.</p>	<p>5. Hold Left/Right button for 2 seconds for digital magnification from 1x to 4x.</p>	<p>6. I ph or</p>

## Parts Description





**!** Press OK Button twice to enter menu. Press and hold OK Button to turn off the microscope.  
( The unit can't be turned on by remote control)

## Detailed Operation

### 1. Charge the Microscope Correctly

Connect the USB cable to the USB-C slot of the monitor and plug the the AC charger to get power.





**Tips:**

- (1) There is one rechargeable battery inside the monitor.
- (2) It takes 2-3 hours to charge the microscope. The power indicator is red during charging and is off after charged fully.
- (3) Short press Power button, current battery level is displayed at the right corner of monitor.



**2. Install Anti-reflective Hood**

The anti-reflective hood is used to reduce the reflection of light. Take off the lens protector, then screw the hood on the lens.





### 3. How to make LED lights work

There are main LED lights (around the camera) and 2 side LED lights ( on the base ),they are powered and controlled separately.



#### (1) How to adjust the brightness of the main LED lights

Press the Left button on the monitor or remote control during Photo/ Video mode to dim the brightness of the main LED light.



Press the Right button on the monitor or remote control during Photo/ Video mode to increase the brightness of the main LED light.



## (2) How to make two side lights work

### Step 1: Provide power for the two side lights.

**Method 1:** Connect the AC charger as the photo below to get power from wall charger.



**Method 2:** Connect the USB cable as the photo below to get power from the battery built-in the monitor.



**Note :** Because the two side lights also gets power to work from the battery, which consumes quickly for method 1, please charge the microscope fully before use.

**Step 2:** Turn on the microscope and roll the scroll wheel on the base to adjust 2 side LED lights.



#### **4. Save Photos/Videos with Micro SD Card**

You can save pictures and videos on a Micro SD card to review later.  
Insert the Micro SD card into slot on back of the monitor.



### 5. Switch Different Mode

There are 3 modes with the microscope: Photo, Video and Playback.

Default setting is photo mode. Press Menu button for 3 seconds to circle the Photo, Video and Playback mode.



Photo, Video, File Mode: Press the MENU button for 3 seconds.

### 6. Capture Pictures

Long Press Menu button switch to photo mode.

Press OK button on the monitor or remote control to capture an image, The photo icon on the screen turns yellow when the image is captured.



## 7. Record Videos

Long Press Menu button switch to Video mode.

Press OK button on the monitor or remote control to record a video, Press OK button on the monitor or remote control again to stop recording, Duration of the recording time is displayed at the top right corner.



## 8. Photos/Videos Playback & Files Management

### Photo/Video Playback

Press MENU button to playback mode.

Press Left and Right button to view photo.


press OK button to view video.



### Files Management

Long press Menu button to Playback mode.

Press Menu button to files management.



Pictures/Videos can be deleted,protected.

Delete: Delete current/Delete all

Protect: Lock the current/Unlock current/Lock all/Unlock all

Press Up/Down button to choose.

Press OK button to confirm.

Press M button to return.

Files locked can't be deleted in the microscope.

(Please note the locked files still can be deleted when deleting files in computer.)

### 9. How to Reset the Microscope

Use a ballpoint pen or similar to push the reset hold on the back of the microscope, keep pressing for 1-2 seconds to reset the device.




## 10. Connect to PC

1. Connect the microscope to the PC with the USB-C cable as the photo below.



2. Turn on the microscope and the screen shows Mass storage, PC Camera and REC\_mode (Note: The REC\_mode enables microscope to charge only).

Press OK button to choose PC Camera to enter PC mode.






## 11. Work with Windows

Please make sure your PC installed “Windows Camera”. If not, get the software in Windows/Microsoft Store or search in on Google.


1. Connect the microscope with computer by the USB cable and turn on the microscope.
2. Press OK button to switch UVC mode.
3. Run “Windows Camera”, it shows the image as below,



## 12. Work with Mac/iMac

Please make sure your MacBook installed “Photo Booth”. If not, get the software from APP store.

1. Connect the microscope with computer by the USB cable and turn it on.
2. Run “Photo Booth”(Finder—Application—Photo Booth) and choose “USB PHY 2.0 (default)”on the menu bar.




#### Please note:

1. Left button/Right/Menu button on the monitor and remote control can't work when the microscope connected with computer.
2. The photos/videos only can be saved in computer (can't not be saved in the Micros SD card in the microscope) when connected with computer.

### Menu Settings

#### 1. Photo Setting

Press Menu button for 3 seconds to enter photo mode, press Menu button once again to enter video setting as the photo below.



**Capture mode**

Set photo capture mode : Single shot / 2s delay/5s delay/10s delay 2s delay refers that the microscope takes photo 2 seconds later after pressing OK button.

### **Resolution**

Set photo resolution : 1.3M/2M/3M/5M/8M/10M/12M

### **Continuous shooting**

Turn on/off continuous shooting .

Turn on this function ,the device take 3 photos continuously after pressing OK button.

Turn off this function ,only 1 photo taken after pressing OK button.

### **Quality**

Set photo quality : Fine/Normal/Economy.

### **Sharpness**

Set photo sharpness : Strong/Normal/Soft.

### **Color**

Set photo color : Normal / Black and white / Sepia

### **Iso**

Set photo ISO : Auto /100/200/400

### **Exposure**

Set photo exposure :+3/+2/+1/+0/-1/-2/-3.

### **Anti-shaking**


Turn on/off anti-shaking function.

### **Date Stamp**

Turn on/off date stamp

## **2. Video Setting**

Press Menu button for 3 seconds to enter video mode, press Menu button once again to enter video setting as the photo below.



### **Resolution**

Set video resolution: 1080FHD/720P/VGA.

### **Cyclic Record**

Set video cycle recording : off/ 3 minutes/5 minutes/10 minutes.

### **High Dynamic Range**

Turn on/off high dynamic range.

### **Exposure**


Set video exposure: +3/+2/+1/+0/-1/-2/-3.

### **Date Stamp**

Turn on/off date stamp.

### 3. System Setting

Press Menu button for 3 seconds to enter video/photo mode, press Menu button twice to enter system setting as the photo below.



#### **Auxiliary to the line**

Set auxiliary line : off/centre/cross. Center refers to center of the monitor.  
Cross refers to the cross line.

#### **Fill light**

Set brightness of the main light : off/1/2/3/4/5.

#### **Screen light off**

Set screen light off :off/ 30 second /1 minute /2 minutes.

#### **Auto power off**

Set auto power off : off/3 minutes/10 minutes.

#### **Frequency**

Set Frequency : 50HZ/60HZ

#### **Language**

Set system language

Chinese/Japanese/English/Korean/French/Germany/Italian/Dutch/  
Portuguese/Spanish/Czech/Polish/Russian/Turkey/Hebrew/Thai

#### **Date/Time**

Set system date and time.

#### **Format**

Format memory card.

#### **Default Setting**

Restore the device to factory setting.

#### **Version**

Display current software version.

### Specifications

Screen Size	7 Inches
Screen Resolution	1024×600
Photo Resolution	12MP 4032×3024 IOMP 3648×2736 8MP 3264×2448 5MP 2592×1944 3MP 2048×1536 2MP 1920×1080 1.3MP 1280×960
Video Resolution	1080FHD 720P VGA(AVI)
Battery	Rechargeable 2000mAh Battery
Focus Range	1cm-21cm
Storage	Micro-SD card, 32GB(recommended) or 64GB
PC View	Supported
IR Remote	Max 5 meters
Language	Chinese/Japanese/English/Korean/French/ Germany/Italian/Dutch/Portuguese/Spanish/ Czech/Polish/Russian/Turkey/H ebrew/Thai
Light Source	6 white-light LEDs integrated in the lens and 2 Side LEDs
Warranty	1 Year Warranty

## FAQ

Why I can't turn on the microscope properly?

Please make sure that the cable is connected correctly and that the product is powered on. If the product still cannot be turned on normally after confirming the correct operation, please reset the microscope.

Please check the mode status of the device .You can't take photos by press OK button during video or playback mode. Please press OK button for 3 seconds to switch to photo mode to take photos.

Why I can't get clear pictures?

Make sure that the LED lights work properly. Adjust the light to illuminate the object, adjust the object distance by using the knob on the holder, and then rotate the focus controller. If you still cannot obtain a clear image, gently wipe the lens with a clean, dry cotton cloth or paper towel to ensure that the lens is clean.

Why does the microscope restart repeatedly?

Please check the power adapter and use a high-quality 5V 2A power adapter.

Why does

After inserting the Micro SD card, please format it first. Enter menu, then choose "SETTING", find "FORMAT", then choose "OK"and press the OK button to confirm the formatting. If it still can't work ,please connect the microscope with computer and format the Micro SD card in computer.

How to reduce image reflection of light ?

Solution 1: Please try to turn low or turn off the white LED. Solution 2: Install anti-reflective hood. Solution 3: Adjust angle between the object and lens to avoid the reflection.

Why are the 2 side lights can't work ?

The 2 side lights and the main LED light are controlled separately, Please refer to page 11 and charge the battery of side lights fully before use.

## Get Support

You are very welcomed to get in touch with us if need any help.



E-mail: [support@tomlov.com](mailto:support@tomlov.com)



Facebook: @Tomlov



WhatsApp: +86 18123716479



Website: [www.tomlov.com](http://www.tomlov.com)

This user's manual has been subjected to thorough reviews and verification to ensure accuracy of its content and contains information that's pertinent to the model and series it describes at the time of publication. However, content of this manual and/or product specifications are subject to change without prior notice and with no direct or indirect liability on tomlov Electronics' part concerning such change. Copyright and registered trademark tomlov, tomlov are registered trademarks.


E-mail: [support@tomlov.com](mailto:support@tomlov.com)

# TOMLOV

**More Than A Microscope**

---

## Documents / Resources

	<p><b><a href="#">TOMLOV DM201 SE LCD Digital Microscope</a></b> [pdf] User Manual DM201 SE, DM201 SE LCD Digital Microscope, LCD Digital Microscope, Digital Microscope, Microscope</p>
---	--

## References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)