




**Contents** [ [hide](#) ]

- [1 VivoTex Scaffolds Suspended Cells Seeding](#)
- [2 Intended use](#)
- [3 Content and storage](#)
- [4 OVERVIEW](#)
- [5 PREPARATION](#)
- [6 PROTOCOL](#)
- [7 TROUBLESHOOTING](#)
- [8 FAQ](#)
- [9 Documents / Resources](#)
  - [9.1 References](#)



## **VivoTex Scaffolds Suspended Cells Seeding**



## Scaffold Designs: Box, Aligned, and Isotropic Patterns

### Intended use

VivoTex™ scaffolds are custom-designed to augment your research, enhancing results that advance your scientific discoveries. VivoTex™ scaffolds are shipped non-sterile, placed inside low-attachment tissue culture plates, and can be used directly for cell culture applications. The scaffold kit does not contain any living cells or materials. VivoTex™ scaffolds are suited for research use only; they are not intended for clinical or household use.

### Content and storage

Catalog number	Content	Storage
VBX.0001.24.24.S	24x <b>Box Pattern</b> MEW Scaffolds 10 µm fibers, 150µm spacing, 200µm thick, PCL, 24 Well, 24 Scaffolds (non-sterile)	Store at cool, dry, well-ventilated area.
	1x 24-well plate, Ultra Low-Attachment (non-sterile)	
VAL.0001.24.24.S	24x <b>Aligned Pattern</b> MEW Scaffolds 10 µm fibers, 100µm spacing, 200µm thick, PCL, 24 Well, 24 Scaffolds (non-sterile)	Store at cool, dry, well-ventilated area.
	1x 24-well plate, Ultra Low-Attachment (non-sterile)	
VIS.0001.24.24.S	24x <b>Isotropic Pattern</b> MEW Scaffolds 10 µm fibers, 150µm spacing, 200µm thick, PCL, 24 Well, 24 Scaffolds (non-sterile)	Store at cool, dry, well-ventilated area.
	1x 24-well plate, Ultra Low-Attachment (non-sterile)	

See kit label for expiration dates. See bottle labels for expiry dates of individual components.

See kit label for expiration dates. See bottle labels for expiry dates of individual components.

## **OVERVIEW**

### **Required materials and equipment (not supplied)**

- Cell suspension: Prepare cells in suspension at a recommended concentration of  $1 \times 10^6$  cells/mL in culture medium.
- Culture Medium: It is recommended to use cell culture medium supplemented with 10% (v/v) serum (e.g. Fetal Bovine Serum) and 1% (v/v) antimicrobial agents (e.g., Penicillin-Streptomycin).
- Water bath or heat block
- Sterile Pipettes and pipette tips
- Biosafety cabinet
- Incubator (37°C, 5% CO<sub>2</sub>)
- Ethanol for work area sterilization
- Sterile forceps (optional for scaffold handling)
- Sterile 1.5 mL microcentrifuge tubes (optional for scaffold handling)

### **Procedural guidelines**

- All steps should be performed using biological aseptic techniques using a Biological Safety cabinet to maintain sterility of the product.
- All kit components can be used at room temperature (18-25°C)
- Use disposable, individually wrapped, sterile plasticware and sterile tubes, pipettes tips, and forceps.
- The biological performance is dependent on the cells used. It is recommended to always include a 2D culture control (without VivoTex™ scaffolds) to confirm cell attachment capacity using cell-culture treated well plates.
- Wear appropriate Personal Protective Equipment in accordance with the biological hazards introduced in your culture protocol (cell type and culture medium)

## **PREPARATION**

## **Sterilization of non-sterile VivoTex™ MEW scaffolds**

- One day before starting your experiment, immerse VivoTex™ MEW scaffolds in 70% ethanol for 30 minutes under sterile conditions.
- Wash the scaffolds with Phosphate Buffer Solution (Solution B) three times (5min per wash) under sterile conditions.

## **Pre-treatment 1: Etching VivoTex™ Scaffolds**

- One day before of your experiment, treat the VivoTex™ scaffolds with the supplied etching solution (solution A) for 10-30 minutes under sterile conditions.
- Wash the scaffolds with Phosphate Buffer Solution (Solution B) three times (5min per wash) under sterile conditions.
- Remove solution B from the scaffolds and allow the scaffolds to air dry inside the biosafety cabinet, leaving the lid off the well plate (1-2h).
- Proceed with the steps for protein adsorption or cell seeding, as outlined below.

## **Pre-treatment 2: Protein adsorption of VivoTex™ Scaffolds**

- Before starting your experiment, treat the VivoTex™ scaffolds with serum-containing media (e.g. 10% FBS) for 2-24 hours under sterile conditions.
- Remove the culture media and proceed with the steps for seeding cells outlined below.
- Remove solution B from the scaffolds and allow the scaffolds to air dry inside the biosafety cabinet, leaving the lid off the well plate (1-2h).

## **PROTOCOL**

### **Seeding of cells onto VivoTex™ MEW Scaffolds**

1. Prepare a cell suspension at  $1 \times 10^6$  cells/mL.
2. Carefully pipette 100  $\mu$ L of cell suspension (containing 100,000 cells) directly onto the center of each VivoTex™ MEW Scaffolds.
3. Allow cells to attach for 1-3 hours at 37°C in the incubator without additional media.  
Note: Do not let the samples dry out. As needed, carefully replenish with 10-40  $\mu$ L of

cell culture medium ~each hour without overflowing the scaffold or dissociating the cells from the sample. The incubation time is dependent on the cell type used. Include a 2D control to confirm cell attachment of the specific cell type used. Confirm initial seeding and attachment under an inverted microscope before adding more media.

4. After incubation, gently add 900  $\mu$ L of cell culture medium to each well to bring the total volume to 1 mL per well.

## **Culturing VivoTex™ MEW Scaffolds**

1. Incubate the well plate with the cell-seeded samples at 37°C and 5% CO<sub>2</sub>.  
Note: The scaffolds are compatible with a wide range of standard CO<sub>2</sub> incubators, tri-gas incubators, and hypoxia chambers designed for cell culture applications. Select the appropriate incubation environment based on the experimental requirements.
2. Replace the culture media according to your standard cell culture protocol (e.g., every 2–3 days).
  - Carefully aspirate the spent media from the side of each well to avoid disturbing the scaffolds.
  - Gently add 1 mL of pre-warmed complete DMEM to each well.  
Note: Ensure all media handling is performed under sterile conditions to maintain culture integrity.
3. Determine culture duration based on the intended downstream application:
  - Short-term culture (24–72 hours): Recommended for studies focused on cell attachment and proliferation.
  - Long-term culture (4–56 days): Suitable for applications involving extracellular matrix production or cellular differentiation.

Note: Optimal culture duration may vary depending on cell type and experimental design. Always refer to your specific assay requirements

## **TROUBLESHOOTING**

- If scaffold flotation is observed during seeding or early culture, consider the following approaches to promote stability:
  - Pre-wet the scaffolds with complete culture media for an extended period prior to cell seeding to enhance saturation and reduce buoyancy.

- Use sterile weighting devices, such as stainless steel or Teflon rings, to gently anchor the scaffolds without compromising sterility or cell viability.


Note: Ensure that any materials used for weighting are biocompatible and have been sterilized according to your laboratory's protocols.

- To promote uniform cell distribution across the scaffold, minimize movement of the culture vessel during the first few hours post-seeding. Furthermore, avoid adding excess culture media too early, as this may dislodge cells from the scaffold surface, particularly at the edges. Smaller media volumes can be used by preparing a more concentrated cell suspension ( $1-10 \times 10^6$  cells/mL) to ensure a droplet based seeding method.
- If poor cell attachment is observed, verify the incubation time required for initial cell adhesion by performing a comparative optimization using standard 2D tissue culture conditions. Note: Scaffolds typically require longer incubation times for effective cell attachment compared to standard tissue culture plastic controls.

## FAQ

- **Q: Can VivoTex™ scaffolds be used for clinical applications?**
  - A: No, VivoTex™ scaffolds are intended for research use only and are not suitable for clinical applications.
- **Q: How should I store the VivoTex™ scaffolds?**
  - A: Store the scaffolds in a cool, dry, and well-ventilated area. Check the kit label for expiration dates.

## Documents / Resources

	<p><a href="#">VivoTex Scaffolds Suspended Cells Seeding [pdf]</a> User Guide</p> <p>VBX.0001.24.24.S, VAL.0001.24.24.S, VIS.0001.24.24.S, Scaffolds Suspended Cells Seeding, Scaffolds, Suspended Cells Seeding, Cells Seeding, Seeding</p>
---	--

## References

- [User Manual](#)

---

## Leave a comment

Your email address will not be published. Required fields are marked \*

Comment \*

Name

Email

Website

☐ Save my name, email, and website in this browser for the next time I comment.

**Post Comment**

**Search:**

e.g. whirlpool wrf535swhz

**Search**

[Manuals+](#) | [Upload](#) | [Deep Search](#) | [Privacy Policy](#) | [@manuals.plus](#) | [YouTube](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.