


WHALETEQ
TRF Generator
Software



WHALETEQ TRF Generator Software User Guide


[Home](#) » [WHALETEQ](#) » WHALETEQ TRF Generator Software User Guide 

Contents

- [1 WHALETEQ TRF Generator Software](#)
- [2 Introduction](#)
- [3 Operation Mode](#)
- [4 Import Database Comparison Results](#)
- [5 Troubleshooting](#)
- [6 Ordering Information](#)
- [7 Documents / Resources](#)
 - [7.1 References](#)



WHALETEQ TRF Generator Software



No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language or computer language, in any form, or by any means, electronic, mechanical, magnetic, optical, chemical, manual or otherwise, without the prior written permission of WhaleTeq Co. LTD.

Disclaimer

WhaleTeq Co. LTD. provides this document and the programs “as is” without warranty of any kind, either expressed or implied, including, but not limited to, the implied warranties of merchantability or fitness for a particular purpose.

This document could contain technical inaccuracies or typographical errors. Changes are periodically made to the information herein; these changes will be incorporated in future revisions of this document. WhaleTeq Co. LTD. is under no obligation to notify any person of the changes.

The following trademarks are used in this document:

WHALETEQ is a registered trademark of WhaleTeq Co. LTD

All other trademarks or trade names are property of their respective holders.

Introduction

Basics

WhaleTeq Universal Medical Standard TRF Generator (“TRF Generator”) is designed to improve the efficiency of performance tests for medical devices. Both DUT testing and report generation can be done in TRF Generator with control over test equipment, medical standards assistant, and IEC Test Report Form unified.

Note: For now, the TRF Generator cannot receive data returned from the DUT. Therefore, the test operator will need to observe the results and fill them in the report record of TRF Generator manually.

System Requirements

A PC with CPU frequency at 1.5GHz or higher and 4G RAM or higher. Recommended OS is Windows 10, with Microsoft Word 2016 or higher installed.

Software Setup

Visit WhaleTeq’s website to download and install TRF Generator on PC.

Software Auto Updates

When TRF Generator launches, it automatically checks whether there’re any updates for TRF generator software or medical standard project files via internet access. Update if needed.

Operation Mode

Login

After purchasing TRF Generator, you will receive an activation mail with a link from WhaleTeq. Click the link to activate the account.




Figure 1: TRF Generator Account Activation Mail

Once activation is done, you will receive another email with login details.

Software Account Login Information - trf-2-25

WhaleTeq <whaleteq@whaleteq.com>

To [REDACTED]

Hello!

Thank you for purchasing WhaleTeq software - trf-2-25

Your account is now created successfully.

The below account name and password can be used to login the software.

Account: [REDACTED]

Password: [REDACTED]

You can also visit [here](#) to download user manual, software and related documents.

Please feel free to contact WhaleTeq [Sales](#) if you have any questions.

WHALETEQ Co., LTD | E whaleteq@whaleteq.com | T +886-2-2596-0701 | www.whaleteq.com

Figure 2: TRF Generator Account Login Information

Launch TRF and click any test item. When [Activate TRF License] window shows, enter the Account email from the login information email and click [Login].




Figure 3: Activate TRF License Window

Then [Authentication] window shows. Enter the Password from the login information email and you will be able to access the full functionality of TRF Generator.

If you have any questions or need help with account activation, feel free to contact WhaleTeq for assistance.




Figure 4: Authentication Window

Download Project File for Different Medical Standards

1. Click [Cloud Library] then [Management].




Figure 5: Download Project File for Different Medical Standards (Step 1)

2. In [Cloud TRF], select the Standard TRF profile or Standard Assistant profile to be used. Once selected, click [Download] to download it to the PC.

Standard TRF contains Questionnaire, IEC Test Report Form, and Standard Assistant. Standard Assistant contains test projects and the control capability for WhaleTeq test equipment.

We will update the profiles once there are any IEC Test Report Form updates available. If you have any questions, email us at service@whaleteq.com. We will get back to you as soon as possible.




Figure 6: Download Project File for Different Medical Standards (Step 2)

3. When the download completes, the profiles will show in [Local TRF].
If an item in [Cloud TRF] is in RED, it indicates that the versions are not matching between Local TRF and Cloud RTF. Please update if needed.




Figure 7: Download Project File for Different Medical Standards (Step 3)

Getting Started

Select the preferred medical standard. The instruction below uses IEC 60601-2-47 as an example.




Figure 8: Medical Standard Options

Click [IEC60601-2-47:2012]. You will see the details that consist of 4 parts:

The screenshot displays the 'Medical Standard Compliance Test Platform' interface. It is divided into four main sections, each highlighted with a red box and a letter:

- A. Questionnaire:** Located on the left, it contains fields for 'Product type: Holter ECG', 'Model name: ?', 'Model difference: N/A', 'Intended use: Recording ECG worn or carried by', 'Operating principles' (Design: Digital Ambulatory recorders, Recording Time: 24 hour, Channel: Channel 1, Channel 2, Channel 3, Patient Electrode: Channel 1+ (RA), Channel 2, Gain setting: 5 mm/mV, 10 mm/mV, 20 mm/mV, Speed: 25 mm/s, 50 mm/s), 'Other capability' (ST segment measurement detection, For infant weighing < 10kg, Capable of recording ECG in the presence of a magnetic field, Capable of recording the activity of an implant), 'Filter: 0.67 ~ 40Hz', 'Algorithm capability' (Accuracy: Algorithm: QRS detection, Heart rate, Supraventricular tachycardia, Physiological report: Heart rate, Supraventricular tachycardia, Test data source: Data from DUT), 'Data pre-processing' (Re-sampling, Reforming, Rescaling, Filtering, Conversion from digital to analogue signal, Other: ?), 'Analysis multi-channel data', and 'Can deactivated automatic analysis mode'.
- B. Test Item:** Located in the top center, it lists various test clauses under '201.12.1.101.2.3.3.1 - Heart rate measurement', '201.12.1.101.2.3.3.2 - Heart rate variability or RR interval variability measurement test', 'Analogue Test', 'Digital Test', '201.12.1.101.2.4 - Run-by-run comparison', '201.12.1.101.2.4.3 - Run sensitivity summary matrix', '201.12.1.101.2.4.4 - Run positive sensitivity summary matrix', '201.12.1.101.2.5 - VF and AF comparison', '201.12.4 - PROTECTION AGAINST HAZARDOUS OUTPUT', '201.12.4.4.101 - Linearity and dynamic range...(100% - Pass)', 'Analogue AMBULATORY RECORDERS', 'Digital AMBULATORY RECORDERS(Alternative 4Hz Sine)', 'Digital AMBULATORY RECORDERS(Alternative 4Hz Sine)', and '201.12.4.4.102 - Input impedance...(33%)'.
- C. Test Record:** Located in the bottom center, it displays a table for 'Channel 1 Channel 2 Channel 3' with columns for '(mV)' and 'DC Offset'. The table has rows for 'Amplitude' (0 mV, 300 mV, -300 mV) and '0.5 mV', '1 mV', '2 mV', '10 mV'.
- D. Test Report Form:** Located on the right, it contains fields for 'Report Reference No.', 'Date of issue', 'Total number of pages', 'CB Testing Laboratory', 'Address', 'Applicant's name', 'Test specification' (Standard: IEC 60601-2-47: 2012 (Second Edition) for use in conjunction with IEC 60601-1:2005 (Third Edition)), 'Test procedure' (CB Scheme), 'Non-standard test method' (N/A), 'Test Report Form No.' (IEC60601_2_47D), and 'Test Report Form Originator' (UL/US).

Figure 9: Medical Standard Compliance Test Screen

- A. Questionnaire: As we know, there are some questions in the body of the standards that need to be answered beforehand. Please answer these questions based on the attributes and the usage of the DUT. Some test projects will determine whether the tests need execution or parameter adjustment.
- B. Test Item: It consists of all the test clauses for the standards. Choose the test clause that you want to test.
- C. Test Record: Displays Report, Description, Procedure, Pass Criterion, and Diagram based on your selection for the test clause.
- D. Test Report Form: IEC Test Report Form with the test results and lets you see how the results display in the report form in real-time.

Select to start a test and fill in the records.

By selecting one of the test clauses, you shall see:

➤ Report

Channel 1 Channel 2		DC Offset		
(mV)		0 mV	300 mV	-300 mV
Amplitude				
0.5 mV				
1 mV				
2 mV				
10 mV				

Report Description Procedure Pass Criterion Diagram

Figure 10: Report Tab

➤ Description

Digital AMBULATORY RECORDERS shall be capable of responding to and displaying input signal of 10 mV peak-to-valley (p-v) in amplitude (when set to the 5 mm/mV GAIN setting) and varying at a rate of 125 mV/s in the presence of a direct current (d.c.) offset voltage of ± 300 mV. The indicated time-varying output signal amplitude referred to input shall not change by more than 10 % or 50 μ V, whichever is greater.
--

Report Description Procedure Pass Criterion Diagram

Figure 11: Description Tab



Procedure

1) Setup ECG to 5mm/mV and 25mm/s.
2) Connect input electrodes to Channel 1+ (RA), Channel 1- (LA)/
Channel 2+ (V1), Channel 2- (V2)/Channel 3+ (V3), Channel 3- (V4);
3) Connect the reference electrode (if used) to RL.
4) Select RA as the output lead wire.
5) Output Triangular wave (0.5 mV, 6.25 Hz) and measure Vpp with Pass Criterion.
6) Set DC Offset to 300 mV and wait for 30s; Measure Vpp with Pass Criterion.
7) Set DC Offset to -300 mV and wait for 30s; Measure Vpp with Pass Criterion.
8) Change 0.5 mV in 2) to 1 / 2 / 10 mV and repeat 5) ~ 7).
9) Change output lead wire to V1/V3 and repeat 4) ~ 8).

ReportDescriptionProcedurePass CriterionDiagram

Figure 12: Procedure Tab



Pass Criterion

The output signal amplitude referred to input shall not change by more than 10% or 50 μ V, whichever is greater - PASS

ReportDescriptionProcedurePass CriterionDiagram

Figure 13: Pass Criterion Tab

➤ Diagram





Figure 14: Diagram Tab


Select any fields in the report. Right-click to see [Start Test] and [Edit Result].

Channel 1 Channel 2


(mV)	DC Offset		
Amplitude	0 mV	300 mV	-300 mV
0.5 mV			
1 mV			
2 mV			
10 mV			




Start Test




Edit Result




Report




Description



Procedure



Pass Criterion



Diagram

Figure 15: Start Test Function

Select [Start Test]. A window for controlling the test equipment and filling test results shows.
The window for controlling the test equipment




Figure 16: The Window for Controlling the Test Equipment

- A. Test Parameter: Shows test parameters required by the test clause.
- B. Executed SDK Routines: Controls the test equipment based on the selected test clause and outputs the corresponding waveform.
- C. Prompt Message: Prompts if manual intervention in the process is needed based on the requirements of the test clause. During the test, please pay attention to the messages and act on them accordingly.
- D. Waveform: Displays the waveform based on the test clause.

Test Result window

The window shows the Pass Criterion of the test clause. The text will be in Blue if the test results are within the scope of the criterion; Red if they are out of the range.




Figure 17: Test Result Window

Fill in the test report with the results according to the criterion of IEC Test Report Form to preview. The text will be in Green if the test results are within the scope of the criterion; Red if they are out of the scope.




Figure 18: Test Report Preview




Figure 19: Test Progress

Actions in [User Data]

- Clean Questionnaire: Clean all the content in the questionnaire.
- Clean All Tables: Clean all the data in test table.
- Export User Data: Stores the test results in the selected folder.
- Import User Data: Imports the results from the last run to continue testing.
- Clean Result Table: Clean the test results of the selected test item.
- Export Table in Word: Exports the results of the selected test table to Word.




Figure 20: Functions in the User Data Tab

Generate Test Report

Select [TRF] and click [Generate Word TRF] to export the test report as a Word document to the folder you select. You can still edit the test report after it is generated.




Figure 21: Export the Test Report as a Word Document

Supported Test Equipment

Table 1: Supported Test Equipment

Item	Equipment Name	Manufacturer	Medical Standards	IEC TRF Form
1	SECG 5.0 AIO SECG 4.0	WhaleTeq	IEC 60601-2-25 IEC 60601-2-27 IEC 60601-2-47 YY078 2-2010 YY0885-2013 Y Y1079-2008 YY1139-2013	IEC 60601-2-25 IEC 60601-2-27 IEC 60601-2-47
2	CMRR 3.0+ CMRR 3.0	WhaleTeq	IEC 60601-2-25 IEC 60601-2-27 IEC 60601-2-47 YY078 2-2010 YY0885-2013 Y Y1079-2008 YY1139-2 013	IEC 60601-2-25 IEC 60601-2-27 IEC 60601-2-47

Import Database Comparison Results

Import Database Comparison Results from WhaleTeq RDCA

When using WhaleTeq RDCA for IEC60601-2-47 database compliance analysis, please import the test files of the 5 ECG databases to WhaleTeq RDCA by databases. It is recommended to name the folders with the name of the ECG databases, as shown in the image below.

Clean	Test Annotation File	Format	QRS +P	VEB Se	VEB +P	VEB FPR	RMS HR Err.
1	209w4.txt	TXT-AAMI/AHA	99.88	0.00	0.00	0.08	2.62
2	228w4.txt	TXT-AAMI/AHA	95.69	83.77	84.33	3.21	32.93
3	210w4.txt	TXT-AAMI/AHA	99.82	92.12	89.41	0.89	4.48
4	111w4.txt	TXT-AAMI/AHA	99.83	100.00	4.17	1.30	5.69
5	112w4.txt	TXT-AAMI/AHA	99.81	-	0.00	0.09	3.79
6	212w4.txt	TXT-AAMI/AHA	99.96	-	-	0.00	1.86
7	113w4.txt	TXT-AAMI/AHA	100.00	-	-	0.00	0.40
8	230w4.txt	TXT-AAMI/AHA	99.89	100.00	50.00	0.05	3.10
9	213w4.txt	TXT-AAMI/AHA	99.74	94.87	95.36	0.38	4.20
10	114w4.txt	TXT-AAMI/AHA	100.00	60.00	100.00	0.00	1.27
11	231w4.txt	TXT-AAMI/AHA	100.00	-	-	0.00	0.77
12	214w4.txt	TXT-AAMI/AHA	100.00	37.74	95.24	0.24	1.66
13	115w4.txt	TXT-AAMI/AHA	100.00	-	0.00	0.06	0.13
14	232w4.txt	TXT-AAMI/AHA	99.60	-	0.00	0.34	4.99
15	116w4.txt	TXT-AAMI/AHA	95.49	88.78	92.55	0.35	33.21
16	215w4.txt	TXT-AAMI/AHA	100.00	83.97	98.21	0.08	0.56
17	233w4.txt	TXT-AAMI/AHA	99.92	59.39	58.22	15.79	3.24
18	117w4.txt	TXT-AAMI/AHA	99.69	-	0.00	0.39	14.74
19	234w4.txt	TXT-AAMI/AHA	100.00	100.00	100.00	0.00	0.13

Figure 22: Import Database Comparison Results from WhaleTeq RDCA (Step 1)

Select to export the reports of 5 ECG databases to 5 folders respectively, as shown in the image below.




Figure 23: Import Database Comparison Results from WhaleTeq RDCA (Step 2)

Under [RDCA Reports] below the questionnaire of TRF Generator, click the database name and select the folder with database comparison results to import. Once imported, the reports of relevant test items will display, as shown in the image below.





Figure 24: Import Database Comparison Results from WhaleTeq RDCA (Step 3)

Import Database Comparison Results from WhaleTeq CDCA

Use WhaleTeq CDCA for IEC 60601-2-25 database comparison and generates 3 different kinds of comparison results. Select [Export Result Report] to export them to the local storage, as shown in the image below.




Figure 25: Import Database Comparison Results from WhaleTeq CDCA (Step 1)

Export the 3 reports to 3 RTF files individually, as shown in the image below.

名稱	修改日期	類型
WhaleTeqIEC 60601-2-25 201.12.1.101.2 (CTS Amplitude) (Digital)201106165822	2020/11/6 下午 04:58	RTF 文件
WhaleTeqIEC 60601-2-25 201.12.1.101.2 (CTS Interval & Duration - Absolute) (Digital)201106165931	2020/11/6 下午 04:59	RTF 文件
WhaleTeqIEC 60601-2-25 201.12.1.101.3.2 (CSE Interval - Biological)201106165951	2020/11/6 下午 04:59	RTF 文件

Figure 26: Import Database Comparison Results from WhaleTeq CDCA (Step 2)

Import the 3 RTF files into TRF Generator. Once imported, the reports of relevant test items shall display, as shown in the image below.




Figure 27: Import Database Comparison Results from WhaleTeq CDCA (Step 3)

Troubleshooting

Q1: What if the test equipment cannot be controlled and output test signals?

A: Please use the test equipment supported by TRF Generator. In case the equipment cannot output signals as expected, reconnect the equipment. Close and reopen TRF Generator software. Do not worry, the data you filled in will be retained when you open it next time. If you would like to save it to local storage, click [User Data] and select [Export User Data].




Figure 28: Save Test Results to a Local Computer

Q2: The test is completed, but I cannot export the report to Word.

A: Please check if Microsoft Word 2016 or higher is correctly installed on the PC. If yes, close any Word instances and try to generate the report again.

Q3: What to do when a message shows that the license is invalid?

A: It means your account has expired and become invalid. Please contact us if you want to keep using it.

Q4: The version of IEC Test Report Form is found not to be the latest one.

A: Please let WhaleTeq know your software version and the environment under which it runs. We will help you with the issue as soon as possible.

Ordering Information

Table 2: Ordering Information

Part No.	Description
HB0-TR00001	TRF Report Generator for IEC 60601-2-25:2011
HB0-TR00002	TRF Report Generator for IEC 60601-2-27:2011
HB0-TR00003	TRF Report Generator for IEC 60601-2-47:2012


Revision History
Table 3: Revision History

Version	Modified Contents	Issued Date
2020-12-23	First release	2020-12-23
2024-05-08	<ul style="list-style-type: none">Update Figure 3 and 4Add 5 Ordering Information 6 Revision History	2024-05-21

Contact WhaleTeq

WHALETEQ Co., LTD
service@whaleteq.com (O)+886 2 2517 6255
8F., No. 125, Songjiang Rd., Zhongshan Dist., Taipei City 104474, Taiwan

Documents / Resources

	WHALETEQ TRF Generator Software [pdf] User Guide TRF Generator Software, Software
---	--

References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.