





## ZEBRA CS6080 Baxter Scanner User Guide

[Home](#) » [ZEBRA](#) » ZEBRA CS6080 Baxter Scanner User Guide 



Baxter 7000-CS6080  
Scanner  
Quick Start Guide



## Contents

- 1 Copyright
- 2 Terms of Use
- 3 Introduction
- 4 Features
- 5 Useful Barcodes
- 6 Installing the Corded USB Converter Cable
- 7 Button Actions
- 8 Assembly Instruction for 7000-APM or 7000-MWS Stands
- 9 Cleaning and Disinfection
- 10 Troubleshooting
- 11 Regulatory Information
- 12 United States and Canada Regulatory
- 13 Service Information
- 14 Documents / Resources
  - 14.1 References

## Copyright

ZEBRA and the stylized Zebra head are trademarks of Zebra Technologies Corporation, registered in many jurisdictions worldwide. All other trademarks are the property of their respective owners. ©2023 Zebra Technologies Corporation and/or its affiliates. All rights reserved.

For further information regarding legal and proprietary statements, please go to:

SOFTWARE: [zebra.com/linkoslegal](https://zebra.com/linkoslegal).

COPYRIGHTS: [zebra.com/copyright](https://zebra.com/copyright).

PATENTS: [ip.zebra.com](https://ip.zebra.com).

WARRANTY: [zebra.com/warranty](https://zebra.com/warranty)

END USER LICENSE AGREEMENT: [zebra.com/eula](https://zebra.com/eula).

## Terms of Use

### Proprietary Statement

This manual contains proprietary information of Zebra Technologies Corporation and its subsidiaries ("Zebra Technologies"). It is intended solely for the information and use of parties operating and maintaining the equipment described herein. Such proprietary information may not be used, reproduced, or disclosed to any other parties for any other purpose without the express, written permission of Zebra Technologies.

### Product Improvements

Continuous improvement of products is a policy of Zebra Technologies. All specifications and designs are subject to change without notice.

### Liability Disclaimer

Zebra Technologies takes steps to ensure that its published Engineering specifications and manuals are correct; however, errors do occur. Zebra Technologies reserves the right to correct any such errors and disclaims liability resulting therefrom.

### Limitation of Liability

In no event shall Zebra Technologies or anyone else involved in the creation, production, or delivery of the accompanying product (including hardware and software) be liable for any damages whatsoever (including, without limitation, consequential damages including loss of business profits, business interruption, or loss of business information) arising out of the use of, the results of use of, or inability to use such product, even if Zebra Technologies has been advised of the possibility of such damages. Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

## Introduction


This guide provides basic instructions for setting up, programming, and using the 7000-CS6080 scanner.

## Firmware

Welch Allyn Connex® Spot Monitors -Firmware 1.40.10 A0002 or greater is required for use with 7000-CS6080.


## Features

### Corded Front View



**NOTE:** The Secondary Button\* is not functional with Baxter devices.

### Bottom View



## Useful Barcodes

### Set Defaults



Return to Factory Defaults

### Night Mode with Vibrate

The night mode disables \*Beep\* sound.



Enable Night Mode Trigger



Disable Night Mode Trigger



Toggle Night Mode


**NOTE:** For additional configurations, please contact your Baxter service representative.

### Installing the Corded USB Converter Cable

Use the USB Converter Cable to connect the corded scanner to a device or stand.


#### Inserting the USB Cable

1. Insert the corded USB Converter Cable into the scanner.
2. Connect the other end of the corded USB Converter Cable to a USB port on the device or stand.




## Removing the USB Cable

To remove the corded USB Converter Cable, press the release and pull the cable away from the scanner.




## To scan a barcode:

1. Aim the scanner at the barcode.
2. Press the Scan button (...).



3. Ensure the aiming dot is centered on the barcode.



The scanner beeps and the LED turns green to indicate a successful \*Good Read\*.

## Button Actions

To actuate the buttons, press the button recess.


For subsequent button actuations, remove force, and then re-apply force.

| Action  | Description       |
|---|-------------------|
| Scan Button (...)                                       |                   |
| Single press  | Trigger           |
| Hold until beam goes off plus 5 seconds longer          | Toggle Night mode |
| Secondary Button (-) Not functional with Baxter devices |                   |

## LED Indicators

| Action               | LED   |
|----------------------|-------|
| <b>Indicator LED</b> |       |
| Good read            | Green |
| No read              | Red   |


## Assembly Instruction for 7000-APM or 7000-MWS Stands



### Below are the steps to assemble

7000-CS6800 Scanner with the 7000-APM or 7000-MWS Stands:

1. Peel off the adhesive sticker from the holster.
2. Attach the holster to the stand, as indicated by the arrow in image 2.
3. After attaching the holster to the stand, you can hang the 7000-CS6800 Scanner using that holster.
4. When the scanner is hung to the holster, please insert the USB converter cable in the hole as indicated by the arrow in image 4.



Below are the steps to assemble 7000-CS6800 Scanner with the 7000-MS3 Stand:

1. Please recycle the holster which is only used with the 7000-APM/7000-MWS stands.
2. The 7000-CS6800 scanner can be placed in the wire basket.

## Cleaning and Disinfection

The 7000-CS6800 is manufactured to withstand the most popular cleaners used in the health care market.

- Isopropyl alcohol 70% (including wipes)
- 91% Isopropyl Alcohol and 9% Water Solution
- 10% Bleach (Sodium Hypochlorite 0.55%) and 90% Water Solution
- 3% Hydrogen Peroxide and 97% water solution
- Mild dish soap
- Azowipe
- Clorox Dispatch Hospital Cleaner Disinfectant Towels with Bleach
- Clorox Formula 409 Glass and Surface Cleaner
- Clorox Healthcare Bleach Germicidal Wipes
- Clorox Healthcare Hydrogen Peroxide Wipes
- Clorox Healthcare Multi-Surface Quat Alcohol Wipes
- Diversey D10 Concentrate Detergent Sanitizer
- Diversey Dimension 256 Neutral Disinfectant Cleaner
- Diversey Oxivir Tb Wipes
- Diversey Virex II 256 One-Step Disinfectant Cleaner
- Metrex CaviCide
- Metrex CaviCide1
- Medipal Alcohol Wipes
- Metrex CaviWipes
- Metrex CaviWipes1

- PDI Easy Screen® Cleaning Wipe
- PDI Sani-Cloth AF3 Germicidal Disposable Wipe
- PDI Sani-Cloth Bleach Germicidal Disposable Wipe
- PDI Sani-Cloth HB Sani-Germicidal Disposable Wipe
- PDI Sani-Cloth Plus Germicidal Disposable Cloth
- PDI Super Sani-Cloth Germicidal Wipe
- Progressive Products Wipes Plus
- Sani Professional Disinfecting Multi-Surface Wipes
- Sani-Hands® Instant Hand Sanitizing Wipes
- SC Johnson Windex Original Glass Cleaner with Ammonia-D
- Sterets Alcowipe
- Steris Coverage Plus Germicidal Surface Wipes
- Veridien Viraguard

## Troubleshooting

| Problem  | Possible Solutions  |
|--|---|
| Aiming pattern does not appear when pressing the scan button.                      | Connect/re-connect the USB cable.   |
| Imager comes on, but scanner does not successfully read the barcode.               | Ensure the scanner is programmed to read the type of barcode being scanned.     |
|  | Ensure the symbol is not defaced. Scan other barcodes of the same barcode type. |
|  | Move scanner closer to or further from barcode.                                 |
| Scanner gets a successful barcode read but does not transmit the data to the host. | Scan the appropriate host type barcode.   |
|  | Re-connect the cable.   |
| Host displays scanned data incorrectly.  | Scan the appropriate host type programming barcode.                             |

## Regulatory Information

This device is approved under Zebra Technologies Corporation.

This guide applies to the following model number CS6080.

All Zebra devices are designed to be compliant with the rules and regulations in the locations they are sold and will be labeled as required.

Any changes or modifications to Zebra equipment not expressly approved by Zebra could void the user's authority to operate the equipment.

Declared maximum operating temperature: 50°C.

**CAUTION:** All components must be dry before connecting to an external power source.

## Health and Safety Recommendations Ergonomic Recommendations

In order to avoid or minimize the potential risk of ergonomic injury, always follow good ergonomic workplace practices. Consult with your local Health & Safety Manager to ensure that you are adhering to your company's safety programs to prevent employee injury.

### Safety in Hospitals and Aircraft

Wireless devices transmit radio frequency energy that may affect medical electrical equipment and aircraft's



operation.

Wireless devices should be switched off wherever you are requested to do so in hospitals, clinics, healthcare facilities or by airline staff. These requests are designed to prevent possible interference with sensitive equipment's.

**WARNING:** When using this product basic safety precautions should always be followed, including the following: The area in which the units are charged should be clear of debris and combustible materials or chemicals. Particular care should be taken where the device is charged in a non-commercial environment.

- Read all the instructions before using the product.

To charge the mobile device, the charger temperatures must be between 0°C and +40°C (+32°F and +104°F).

Do not use incompatible chargers. Use of an incompatible charger may present a risk of fire, explosion, leakage, or other hazard. If you have any questions about the compatibility of a charger, contact Baxter service support.

Do not leave or store the equipment in or near areas that might get very hot, such as in a parked vehicle or near a radiator or other heat source. Do not place the device into a microwave oven or dryer.

To reduce the risk of injury, close supervision is necessary when used near children.

Please follow local regulations to promptly dispose of used re-chargeable charger.

If you suspect damage to your equipment, contact Baxter service support to arrange for inspection.

Do not use a power pack in excess of its output rating.

Overload outputs may result in risk of fire or injury to persons.

Have service performed by a qualified Baxter service representative. This will ensure safety of the product is maintained.

## **RF Exposure Guidelines**

### **Safety Information**

#### **Reducing RF Exposure – Use Properly**

Only operate the device in accordance with the instructions supplied.

The device complies with internationally recognized standards covering human exposure to electromagnetic fields. For information on international human exposure to electromagnetic fields, refer to the Zebra Declaration of Conformity (DoC) at [zebra.com/doc](http://zebra.com/doc).

Use only Zebra tested and approved headset, belt-clips, holsters, and similar accessories to ensure RF exposure compliance. If applicable, follow the instructions for use as detailed in the accessory guide.

The use of third-party belt clips, holsters, and similar accessories may not comply with RF exposure compliance requirements and should be avoided.

For further information on the safety of RF energy from wireless devices, refer to RF exposure and assessment standards section at: [zebra.com/responsibility](http://zebra.com/responsibility).

## **LED Devices**

Classified as “EXEMPT RISK GROUP” according to IEC 62471:2006 and EN 62471:2008.

Pulse Duration: 22.826 ms (CS6080 with scanner SE4100)

## **Marking and European Economic Area (EEA)**

### **Statement of Compliance**

Zebra hereby declares that this equipment is in compliance with Directives 2014/53/EU and 2011/65/EU.

Any radio operation limitations within EEA countries are identified in Appendix A of EU Declaration of Conformity.

The full text of the EU Declaration of Conformities is available at: [zebra.com/doc](http://zebra.com/doc).

EU Importer : Zebra Technologies B.V

Address: Mercurius 12, 8448 GX Heerenveen, Netherlands

## **Waste Electrical and Electronic Equipment (WEEE)**



For EU and UK Customers: For products at the end of their life, please refer to recycling/disposal advice at: [zebra.com/weee](http://zebra.com/weee).

For EU and UK Customers: All products at the end of their life must be returned to Zebra for recycling. For information on how to return product, please go to: [www.zebra.com/weee](http://www.zebra.com/weee).

## United States and Canada Regulatory

### Radio Frequency Interference Notices

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.



**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help.

### Radio Frequency Interference

#### Requirements- Canada

Innovation, Science and Economic Development Canada ICES-003 Compliance Label: CAN ICES-3 ([B])/NMB-3([B])

This device complies with Innovation, Science and Economic Development Canada's licence-exempt RSSs. Operation is subject to the following two conditions: (1) This device may not cause interference; and (2) This device must accept any interference, including interference that may cause undesired operation of the device.

### RF Exposure Requirements – FCC and ISED

The FCC has granted an Equipment Authorization for this device with all reported SAR levels evaluated in compliance with the FCC RF emission guidelines. SAR information on this device is on file with the FCC and can be found under the Display Grant section of [fcc.gov/oet/ea/fccid](http://fcc.gov/oet/ea/fccid).

### India

The sending and receiving of texts is supported in the following Indian languages:

Assamese, Bengali, Gujarati, Hindi, Kannada, Kashmiri, Konkani, Malayalam, Manipuri, Marathi, Nepali, Oriya, Punjabi, Sanskrit, Sindhi, Tamil, Telugu, Urdu, Bodo, Santhali, Maithili, and Dogri.

## Service Information

Before you use the unit, it must be configured to operate in your facility's network and run your applications.

If you have a problem using your equipment, contact your facility's Technical or Systems Support. If there is a problem with the equipment, please go to CSM website to contact the Baxter Customer Support at [hillrom.com](http://hillrom.com)

For the latest version of this guide go to: [hillrom.com](http://hillrom.com).

### Warranty

Barcode scanner has limited 1 year warranty from date of the purchase.


#### For Australia Only:

This warranty is given by Zebra Technologies Asia Pacific Pte.Ltd., 71 Robinson Road, #05-02/03, Singapore 068895, Singapore. Our goods come with guarantees that cannot be excluded under the Australia Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably

foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. Zebra Technologies Corporation Australia’s limited warranty above is in addition to any rights and remedies you may have under the Australian Consumer Law. If you have any queries, please call Zebra Technologies Corporation at +65 6858 0722. You may also visit our website: [zebra.com](http://zebra.com) for the most updated warranty terms.

**Manufactured by:**  
Zebra Technologies Corporation  
Lincolnshire, IL U.S.A.  
[www.zebra.com](http://www.zebra.com)

©2023 Zebra Technologies Corporation and/or its affiliates. All rights reserved. ZEBRA and the stylized Zebra head are trademarks of Zebra Technologies Corporation, registered in many jurisdictions worldwide. All other trademarks are the property of their respective owners.



MN-004769-01EN-P Revision A – May 2023

Documents / Resources

|   |  |
|---|--|
| <div>Baxter<br/>7000-CS6080<br/>Scanner<br/>Quick Start Guide</div> <div></div> | <div><a href="#">ZEBRA CS6080 Baxter Scanner</a> [pdf] User Guide<br/>CS6080 Baxter Scanner, CS6080, Baxter Scanner, Scanner</div> |
|---|--|

References

- [Advancing Connected Care | Hillrom](#)
- [Patents | Zebra](#)
- [Patents | Zebra](#)
- [Advancing Connected Care | Hillrom](#)
- [Patents | Zebra](#)
- [Zebra Technologies | Visible. Connected. Optimized.](#)
- [Contact Zebra | How to Buy, Support & Repair, Questions](#)
- [Copyright and Trademarks | Zebra](#)

- [!\[\]\(effbd7993c63c039a58fd3395789cf3f\_img.jpg\) Declarations of Conformity | Zebra Technologies](#)
- [!\[\]\(144980d038f2541d7b588a8a9132bd70\_img.jpg\) Software Product Categories and End User License Agreements | Zebra](#)
- [!\[\]\(c4ce2d477989700c971cf3d240ad9283\_img.jpg\) RF Exposure and Assessment Standards | Zebra](#)
- [!\[\]\(5013555a72072875cb154b597e002a46\_img.jpg\) Support and Downloads | Zebra](#)
- [!\[\]\(bf2038c114ec21ea58ad011774351c98\_img.jpg\) Product Warranty Information | Zebra](#)
- [!\[\]\(1ad0c3425edfa4762c2f20e33e3e5bbf\_img.jpg\) WEEE and EU Battery Directive FAQs | Zebra](#)
- [!\[\]\(74d2fc5645add84f8511beb934060048\_img.jpg\) Zebra Technologies | Visible. Connected. Optimized.](#)
- [!\[\]\(ddc5533caa0187e636e3d3234e0983a3\_img.jpg\) Copyright and Trademarks | Zebra](#)
- [!\[\]\(7c207f8f59385c6dd11f9d9bdc7a0d1d\_img.jpg\) Declarations of Conformity | Zebra Technologies](#)
- [!\[\]\(57a1bf910af99362b80b3ac4f2eecbac\_img.jpg\) Software Product Categories and End User License Agreements | Zebra](#)
- [!\[\]\(6b114000ab07dda576e2920e2dc838fa\_img.jpg\) RF Exposure and Assessment Standards | Zebra](#)
- [!\[\]\(9129a6a4a4b11facb5cf665660eef788\_img.jpg\) Support and Downloads | Zebra](#)
- [!\[\]\(7cbd4924d31aa4c46df484d5c5bf4696\_img.jpg\) Product Warranty Information | Zebra](#)
- [!\[\]\(b3461332979d340882fda628798720af\_img.jpg\) WEEE and EU Battery Directive FAQs | Zebra](#)
- [!\[\]\(5dc380853f1779b9d7a8d6fbbcdbb357\_img.jpg\) Declarations of Conformity | Zebra Technologies](#)
- [!\[\]\(2894a5ead456ff0145f35b01b4a6c7df\_img.jpg\) Support and Downloads | Zebra](#)
- [User Manual](#)