



Manuals.plus /

› Drive Medical /

› User Manual for Drive Medical 15033BV-PKG-T Electric Bed

## Drive Medical 15033BV-PKG-T

# Drive Medical 15033BV-PKG-T Electric Bed

## USER MANUAL

### Introduction

---

This manual provides essential information for the safe and effective setup, operation, and maintenance of your Drive Medical 15033BV-PKG-T Electric Bed. Please read this manual thoroughly before using the product to ensure proper function and to prevent injury.

The Drive Medical 15033BV-PKG-T Electric Bed is designed for ease of use and assembly. It features a lightweight, interchangeable design with durable head and foot boards. The bed includes a quiet, self-contained motor for height and section adjustments, and a battery backup for emergency lowering during power outages.

### Safety Information

---

Always follow basic safety precautions when using this electric bed to reduce the risk of injury or damage.

- Ensure all components are securely assembled before use.
- Do not exceed the maximum recommended load of 162 pounds.
- Keep hands and feet clear of moving parts during operation.
- Do not allow children or unauthorized persons to operate the bed.
- Disconnect power before cleaning or performing maintenance.
- Use only manufacturer-approved replacement parts.

### Setup and Assembly

---

The Drive Medical Electric Bed is designed for tool-free assembly. Follow these steps for initial setup:

1. **Unpack Components:** Carefully remove all parts from the packaging. Verify all components are present according to the packing list.
2. **Assemble Bed Frame:** Connect the head and foot sections of the bed frame. Ensure all locking mechanisms are engaged.
3. **Install Head and Foot Boards:** The universal style head and foot boards are interchangeable. Slide them into the designated slots on the bed frame until they are secure.

4. **Motor Installation:** The height adjustment motor arrives pre-installed on the foot section. If it needs to be re-installed or removed, it can be done tool-free.
5. **Position Bed:** Place the bed in its desired location, ensuring there is adequate space around it for operation and access.
6. **Connect Power:** Plug the bed's power cord into a grounded electrical outlet.
7. **Install Mattress:** Place the mattress (if included or purchased separately) onto the bed frame.



*Image: The bed frame assembly, illustrating the motor and manual crank for height adjustment.*



*Image: A fully assembled Drive Medical Electric Bed, complete with mattress and optional side rails.*

## Operating Instructions

The bed is operated using a hand pendant with large, easy-to-use controls.

- **Adjusting Bed Height:** Use the designated buttons on the hand pendant to raise or lower the entire bed frame. The bed height can be adjusted from a minimum of 12 inches to a maximum of 19.5 inches (measured from the floor to the bed deck).
- **Adjusting Head Section:** Use the appropriate buttons to raise or lower the head section of the bed.
- **Adjusting Foot Section:** Use the appropriate buttons to raise or lower the foot section of the bed.
- **Simultaneous Adjustments:** Some pendants may allow for simultaneous adjustment of head and foot sections.

## Power Failure Operation

In the event of a power failure, the bed is equipped with a 9-volt battery backup system. This battery allows for emergency lowering of the head and foot sections up to nine times. To activate, locate the 9-volt battery in the motor housing and ensure it is properly connected. Use the hand pendant as normal to lower the sections.

## Maintenance

Regular maintenance ensures the longevity and safe operation of your electric bed.

- **Cleaning:** Clean the bed frame and head/foot boards with a damp cloth and mild detergent. Avoid abrasive cleaners. Ensure the bed is unplugged before cleaning.
- **Inspection:** Periodically inspect all cables, connections, and moving parts for wear or damage. Ensure all bolts and fasteners are tight.
- **Battery Check:** Check the 9-volt battery in the motor housing periodically and replace it as needed to ensure the power failure lowering function is operational.
- **Lubrication:** Moving parts are generally self-lubricating, but consult a qualified technician if any squeaking or

stiffness occurs.

## Troubleshooting

If you encounter issues with your electric bed, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Bed does not operate (no power).	Power cord unplugged; circuit breaker tripped; faulty outlet.	Check power cord connection; reset circuit breaker; test outlet with another device.
Motor runs but bed does not move.	Mechanism jammed; loose connection to motor.	Inspect for obstructions; ensure motor connections are secure. Do not force movement.
Bed only lowers during power outage.	9-volt battery engaged; main power not connected.	Ensure main power cord is securely plugged in. The 9V battery is for emergency lowering only.
Unusual noises during operation.	Loose components; lack of lubrication; motor issue.	Check for loose parts and tighten. If noise persists, contact customer support.

If the problem persists after attempting these solutions, please contact Drive Medical customer support.

## Specifications

<b>Model Number</b>	15033BV-PKG-T
<b>Brand</b>	Drive Medical
<b>Material</b>	Metal
<b>Color</b>	Brown
<b>Product Dimensions (L x W x H)</b>	88.5 x 35.75 x 35.75 inches
<b>Item Weight</b>	169 pounds
<b>Maximum Recommended Load</b>	162 pounds
<b>Bed Height (Min)</b>	12 inches
<b>Bed Height (Max)</b>	19.5 inches
<b>Batteries Required</b>	No (1 x 9V battery for emergency lowering only)
<b>Special Feature</b>	Lightweight, Interchangeable Head/Foot Boards
<b>UPC</b>	822383211848

## Warranty and Returns

This product is subject to the return policy of the retailer from which it was purchased. Please note the following general return guidelines:

- Returns may not be accepted after 14 calendar days from the purchase date without a valid Return Merchandise Authorization (RMA) number.
- Returns may be subject to a restocking fee (e.g., 15%) plus the cost of outbound and inbound shipping.

- Return requests beyond 14 days may be considered on a case-by-case basis and could incur a higher restocking fee (up to 60%).

For specific warranty information, please refer to the documentation provided with your purchase or contact Drive Medical directly. This section outlines return policy, not a manufacturer's warranty.

## **Customer Support**

---

For further assistance, technical support, or inquiries regarding your Drive Medical Electric Bed, please contact Drive Medical customer support through their official channels or visit their online store.

You can visit the official Drive Medical store for more information:[Drive Medical Store](#)