

## Aidapt VB539

# Aidapt VB539 Solo Deluxe Shower Chair User Manual

Model: VB539

## 1. INTRODUCTION

---

Thank you for choosing the Aidapt VB539 Solo Deluxe Shower Chair. This manual provides essential information for the safe and effective use, assembly, and maintenance of your shower chair. Please read these instructions carefully before use and retain them for future reference.

The Solo Deluxe Shower Chair is designed to provide a safe and comfortable seating solution within a shower environment. It features adjustable height, independently lifting armrests for easy transfers, and a durable, rust-resistant construction.

## 2. SAFETY INFORMATION

---

- **Maximum User Weight:** Do not exceed the maximum user weight of 190 kg (30 stones).
- **Stability:** Ensure all four legs are evenly adjusted and firmly on the shower floor before use. Check for any wobbling.
- **Assembly:** Ensure all components are correctly assembled and securely fastened before first use.
- **Inspection:** Regularly inspect the chair for any signs of wear, damage, or loose parts. Discontinue use immediately if any damage is found and contact Aidapt for assistance.
- **Transfers:** Use caution when transferring on and off the chair. Utilize the independently lifting armrests for support as needed.
- **Cleaning:** Clean the chair regularly as per the maintenance instructions to prevent build-up and ensure hygiene.
- **Children:** Keep children away from the shower chair unless supervised by an adult.

## 3. PACKAGE CONTENTS

---

Please check the contents of your package to ensure all parts are present:

- Shower Chair Frame (lightweight, rust-resistant aluminum)
- Clip-on Seat (injection-molded with drainage holes)
- Backrest (with drainage holes)
- Adjustable Legs with non-slip ferrules

- Integrated Armrests (pre-attached to frame, independently lifting)

## 4. ASSEMBLY INSTRUCTIONS

---

The Aidapt VB539 Solo Deluxe Shower Chair is designed for simple assembly. No tools are typically required.

1. **Unpack Components:** Carefully remove all parts from the packaging.
2. **Attach Seat:** Align the clip-on seat with the frame and press down firmly until it securely clicks into place.
3. **Attach Backrest:** Slide the backrest into the designated slots on the frame until it locks.
4. **Adjust Leg Height:** Press the spring-loaded pins on each leg and slide the inner leg tube up or down to the desired height. Ensure all four legs are adjusted to the same height for stability. The legs can be precisely adjusted to compensate for uneven shower trays. Release the pin to lock the leg in position.
5. **Final Check:** Before first use, ensure all connections are secure and the chair is stable on a flat surface.



Image 1: The Aidapt VB539 Solo Deluxe Shower Chair fully assembled, showing its adjustable legs and armrests.

## 5. OPERATING INSTRUCTIONS

---

1. **Placement:** Position the shower chair securely within your shower area, ensuring all four non-slip ferrules are in full contact with the shower floor.
2. **Height Adjustment:** If further height adjustment is needed, refer to step 4 of the Assembly Instructions. Ensure the chair is empty and stable before adjusting.

3. **Using Armrests:** The integrated armrests can be lifted independently. To lift an armrest, gently pull it upwards. This facilitates easier side transfers onto and off the chair. Push it down firmly to lock it back into position for support during showering.
4. **Sitting:** Approach the chair carefully. Use the armrests for support as you slowly lower yourself onto the seat. Ensure you are fully seated before beginning your shower.
5. **Standing:** Use the armrests for support as you push up to a standing position. Take your time and ensure you have full balance before moving away from the chair.

## 6. MAINTENANCE AND CLEANING

---

Regular cleaning and maintenance will ensure the longevity and hygiene of your shower chair.

- **Cleaning:** Clean the chair regularly with a mild, non-abrasive household cleaner and warm water. Wipe dry immediately after cleaning to prevent water spots and maintain the finish.
- **Removable Parts:** The seat and backrest are removable for thorough cleaning. Simply unclip them from the frame, clean, and reattach securely.
- **Rust-Resistant Frame:** The aluminum frame is designed to be rust-resistant. However, prolonged exposure to harsh chemicals or standing water may affect its finish.
- **Inspection:** Periodically check all bolts, screws, and connections to ensure they are tight. Inspect the non-slip ferrules on the legs for wear and tear; replace if necessary.
- **Storage:** Store the chair in a dry place when not in use for extended periods.

## 7. TROUBLESHOOTING

---

- **Chair is Unstable/Wobbles:**
  - Ensure all four legs are adjusted to the same height.
  - Verify that the spring-loaded pins on each leg are fully engaged in the adjustment holes.
  - Check if the shower floor is uneven. Adjust leg heights slightly to compensate if possible, ensuring stability.
- **Armrests do not stay up/down:**
  - Ensure the armrests are pushed firmly into their locked position when down, or fully lifted when up.
  - Inspect the armrest hinge mechanism for any obstructions or damage.
- **Seat/Backrest feels loose:**
  - Ensure the clip-on seat and backrest are fully engaged and clicked into their respective positions on the frame.

If you encounter any issues not listed here, or if troubleshooting steps do not resolve the problem, please contact Aidapt customer support.

## 8. SPECIFICATIONS

---

Feature	Specification
Model Number	VB539
Brand	Aidapt
Material	Rust-resistant Aluminum frame, Injection-molded polymer seat/backrest

Feature	Specification
Color	White
Overall Dimensions (L x W x H)	590 x 610 x 860-970 mm (approx. 59 x 61 x 86-97 cm)
Overall Depth	535 mm
Seat Depth	380 mm (38 cm)
Seat Height (Adjustable)	465 - 615 mm
Seat Width	490 mm
Width between Arms	465 mm
Height from Seat to Top of Back	420 mm
Item Weight	5.5 kg (12.13 lbs)
Maximum User Weight	190 kg (30 stones)
Batteries Required	No
Country of Origin	China
EAN	05021196594970



Image 2: Detailed dimensional and weight specifications for the shower chair.

## 9. WARRANTY AND SUPPORT

---

For warranty information, please refer to the documentation provided at the time of purchase or contact your retailer. Aidapt is committed to providing quality products and customer satisfaction.

If you have any questions, require technical support, or need replacement parts, please contact Aidapt customer service through their official website or the contact information provided with your purchase.