



[Manuals.plus](#) /

> [Drive Medical](#) /

> Drive Medical Shower Chair with Back - Adjustable Bathroom Bench - Instruction Manual

Drive Medical RTL12202KDR

Drive Medical Shower Chair with Back - Instruction Manual

Model: RTL12202KDR

1. INTRODUCTION

This manual provides instructions for the safe and proper assembly, operation, and maintenance of your Drive Medical Shower Chair with Back. This adjustable bathroom bench is designed to offer stability, comfort, and safety during bathing for individuals requiring additional support.

2. SAFETY INFORMATION

WARNING: This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm. (Proposition 65 Warning)

- Always ensure all components are securely fastened before each use.
- The chair features non-slip suction feet designed to maintain grip on wet surfaces. Verify their secure placement before use.
- Drainage holes in the seat prevent water accumulation, reducing slip hazards.
- Do not exceed the maximum weight capacity of 300 lbs.
- Consult a healthcare professional for guidance on proper use if you have specific mobility needs or are recovering from surgery.

3. PACKAGE CONTENTS

Before assembly, ensure all the following components are present:

- Contoured Seat with Drainage Holes
- Backrest
- Two (2) U-shaped Leg Frames
- Four (4) Adjustable Legs with Non-Slip Suction Feet

- Four (4) Short Thumb Screws
- Two (2) Long Thumb Screws



Image: Components for tool-free assembly.

4. ASSEMBLY INSTRUCTIONS

The Drive Medical Shower Chair is designed for quick and easy tool-free assembly. Follow these steps:

1. **Attach Leg Frames to Seat:** Locate the bottom of the seat. Insert the upper parts of the two U-shaped leg frames into the designated grooves on the underside of the seat.
2. **Secure Leg Frames:** Use the four (4) short thumb screws to secure the leg frames to the seat. Tighten them sufficiently to hold the parts together, but do not overtighten yet.
3. **Attach Backrest Supports:** Take the two L-shaped metal rods (backrest supports). Insert them into the corresponding slots on the backrest.
4. **Attach Backrest to Seat:** Align the backrest supports with the rear holes on the seat's leg frames. Use the two (2) long thumb screws to secure the backrest to the seat.
5. **Install Adjustable Legs:** Pinch the pin on each adjustable leg and slide it into the bottom of each leg frame. Adjust to your desired height (see Section 5.1 for details). Ensure the pin snaps securely into a hole.
6. **Final Tightening:** Once all components are in place and aligned, double-check and fully tighten all thumb screws to ensure the chair is stable and secure. Do not overtighten to avoid damaging the plastic components.

Video: Step-by-step assembly guide for the Drive Medical Shower Chair.

This video demonstrates the process of unboxing and assembling the shower chair, including attaching the seat, backrest, and adjustable legs.



Image: Visual guide to key features and tool-free assembly.

5. OPERATING INSTRUCTIONS

5.1. Height Adjustment

The shower chair's legs are adjustable in 1-inch increments to customize the seat height. To adjust:

1. Locate the push-button pins on each leg.

2. Press the pin and slide the inner leg tube up or down to the desired height.
3. Ensure the pin fully protrudes through the corresponding hole, locking the leg securely in place.
4. Repeat for all four legs, ensuring they are set to the same height for stability.

Video: Demonstration of height adjustment and other features.

This video highlights the adjustable legs, non-slip feet, and the overall functionality of the shower chair.

Safe & Helpful Features



Removable Backrest

Customize the chair for your needs and make transfers easier



In-Seat Handles

Help you raise and lower yourself while getting in and out of the bathtub



Durable Construction

Blow-molded plastic seat and corrosion-proof aluminum frame offer dependable support



Tool-Free Assembly

Simple to assemble in minutes

Image: Close-up of the height adjustment mechanism.

5.2. Removable Backrest

The backrest can be removed to convert the chair into a stool, offering versatility for different bathing and transfer needs.

1. To remove, loosen the two long thumb screws securing the backrest to the seat frame.
2. Carefully slide the backrest upwards and off the support rods.
3. To reattach, align the backrest supports with the holes and secure with the long thumb screws.

Video: Demonstration of backrest removal and portability.

This video illustrates how to easily remove the backrest and highlights the chair's lightweight design for portability.



Image: Removable backrest feature.

5.3. In-Seat Handles and Drainage

The contoured seat includes integrated handles for assistance when sitting down or standing up. The multiple drainage holes prevent water from pooling, ensuring a cleaner and more stable seating surface.



Image: In-seat handles for support.



Image: Drainage holes on the seat.

6. MAINTENANCE AND CARE

To ensure the longevity and hygiene of your shower chair:

- Clean the chair regularly with mild soap and water. A scrub brush can be used for textured surfaces.
- Rinse thoroughly after cleaning to remove all soap residue.
- The aluminum frame is corrosion-proof; however, regular drying after use can help maintain its appearance.
- Periodically check all screws and connections to ensure they remain tight. Retighten as necessary to maintain stability.
- Inspect the non-slip suction feet for wear or damage. Replace if they show signs of deterioration to ensure continued safety.

7. TROUBLESHOOTING

7.1. Chair feels unstable or wobbly

- Ensure all four adjustable legs are set to the same height.
- Verify that all thumb screws (both short and long) are fully tightened.
- Check that the non-slip suction feet are making full contact with the shower or tub surface.

7.2. Screws appear to be pushing through the seat

- This may occur if the incorrect length screws are used in certain positions. Ensure that the shorter thumb screws are used for attaching the main leg frames to the seat, and the longer thumb screws are used for the backrest attachment.
- Avoid overtightening the screws, as this can stress the plastic material.

8. SPECIFICATIONS

Feature	Detail
Model Number	RTL12202KDR
Material	Aluminum frame, Plastic seat
Color	Grey
Product Dimensions	11.5"D x 19.5"W x 28"H (Adjustable height)
Adjustable Height Range	15.5" - 19.5"
Item Weight	7 Pounds
Weight Capacity	Up to 300 lbs
Features	Removable backrest, In-seat handles, Non-slip feet, Drainage holes, Tool-free assembly



Image: Key features and dimensions of the shower chair.

9. WARRANTY AND SUPPORT

For warranty information or product support, please refer to the documentation included with your purchase or contact Drive Medical customer service directly. Please note that return policies are set by the

manufacturers. Returns may be subject to specific conditions, including restocking fees and time limits, as outlined in the purchase agreement.