

- › [Serene Innovations](#) /
- › [Serene Innovations Renew Dryer & Freshener \(Model DB-100\) Instruction Manual](#)

Serene Innovations DB-100

Serene Innovations Renew Dryer & Freshener (Model DB-100) Instruction Manual

Model: DB-100 | Brand: Serene Innovations

1. INTRODUCTION

The Serene Innovations Renew Dryer & Freshener is designed to maintain the cleanliness and longevity of your hearing aids, cochlear implants, and wireless earbuds. This electronic dehumidifier effectively removes moisture and helps eliminate harmful microscopic elements, ensuring optimal performance and extending the lifespan of your delicate devices.

Regular use of this device is crucial for preventing moisture-related damage, which can lead to costly repairs and reduced sound quality. This manual provides detailed instructions for setting up, operating, and maintaining your Renew Dryer & Freshener.



Image 1: The Serene Innovations Renew Dryer & Freshener with its lid open, revealing the blue drying tray and control panel.

2. SETUP

1. **Unpack the Device:** Carefully remove the Renew Dryer & Freshener from its packaging. Ensure all components are present: the main unit, the power adapter, and the removable tray.
2. **Placement:** Place the unit on a flat, stable surface, such as a bedside table or countertop, away from direct sunlight or heat sources.
3. **Connect Power:** Insert the small end of the power adapter into the DC input jack on the side of the Renew Dryer & Freshener. Plug the other end of the power adapter into a standard electrical outlet. The power indicator light on the control panel should illuminate.



Image 2: The Renew Dryer & Freshener positioned on a nightstand, ready for use.

3. OPERATING INSTRUCTIONS

3.1 Drying Your Devices

1. **Prepare Devices:** Before placing your hearing aids or other devices into the dryer, remove any batteries. For hearing aids, open the battery doors to allow for better air circulation.
2. **Place Devices:** Open the lid of the Renew Dryer & Freshener. Place your hearing aids, cochlear implants, or wireless earbuds onto the blue drying tray. The spacious interior can accommodate up to four hearing aids or a combination of devices.
3. **Select Drying Cycle:** Use the slide switch labeled "Lo Hum," "Norm Hum," "Hi Hum" to select the appropriate drying time based on your environment's humidity level.
 - **Lo Hum (2 hours):** For low humidity environments or minimal moisture.
 - **Norm Hum (4 hours):** For normal humidity levels.
 - **Hi Hum (8 hours):** For high humidity environments or significant moisture.
4. **Start Drying:** Close the lid. The drying cycle will automatically begin once the lid is securely closed. The unit warms up to 113°F (45°C)

to facilitate moisture evaporation.

5. **UV Refresh Cycle:** To activate the UV light for disinfection, press the "Refresh" button on the control panel. The UV lights will operate for a set period to help eliminate harmful microscopic elements.



Image 3: The spacious interior of the dryer, capable of holding multiple devices like hearing aids and wireless earbuds.

DRYING FOR ANY HUMIDITY ENVIRONMENT

Choose from 3 dry time settings based on your local humidity.



Image 4: Adjusting the humidity setting on the control panel to optimize drying for different environments.

3.2 Battery Management

- **Battery Storage:** The Renew Dryer & Freshener features a dedicated magnetic strip for battery storage. Place your hearing aid batteries on this strip while your devices are drying to prevent them from getting lost.
- **Battery Tester:** A built-in battery tester allows you to check the charge level of your hearing aid batteries. Place a battery on the metal platform next to the battery indicator lights and push it against the metal stick. The indicator lights will show the remaining charge.

EASY BATTERY TRACKING

- Magnetic battery holding strip
- Integrated battery tester
- Rechargeable hearing aid safe



Image 5: The battery management features, including the magnetic strip for storage and the integrated battery tester.

3.3 Official Product Video

Your browser does not support the video tag.

Video 1: An overview of the Serene Innovations Renew Dryer & Freshener, demonstrating its features and ease of use.

4. MAINTENANCE

To ensure optimal performance and hygiene, regular cleaning of your Renew Dryer & Freshener is recommended.

- **Removable Tray:** The blue drying tray is removable for easy cleaning. Simply lift the tray out of the unit.
- **Cleaning the Tray:** Wash the removable tray with mild soap and water. Rinse thoroughly and allow it to air dry completely before placing it back into the unit.
- **Cleaning the Unit:** Wipe the exterior of the dryer with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.



EASY TO CLEAN

- Simply remove the inner tray
- Wipe off and clean with water and soap



Image 6: Demonstrating the removal of the inner tray for easy cleaning.

5. TROUBLESHOOTING

If you encounter any issues with your Serene Innovations Renew Dryer & Freshener, please refer to the following common solutions:

• Device Not Drying:

- Ensure the power adapter is securely plugged into both the unit and the electrical outlet.
- Verify that the lid is fully closed. The drying cycle will not start if the lid is open.
- Check the selected humidity setting. For very humid conditions, ensure the "Hi Hum" setting is selected.

• UV Lights Not Activating:

- Press the "Refresh" button firmly.
- Ensure the unit is powered on and the lid is closed.

• Battery Tester Not Working:

- Ensure the battery is correctly placed on the metal platform and pressed against the metal stick.
- Verify that the battery is making proper contact with both metal points.

If the problem persists after attempting these troubleshooting steps, please contact Serene Innovations customer support for further assistance.

6. SPECIFICATIONS

Model Number	DB-100
Product Dimensions	5.8 x 3.5 x 2 inches
Item Weight	1.04 Pounds
Drying Temperature	Approximately 113°F (45°C)
Drying Cycles	2, 4, or 8 hours (selectable)
UV Light	Yes, for disinfection (Refresh function)
Battery Tester	Integrated
Battery Storage	Magnetic strip
Power Source	AC Adapter (included)

7. WARRANTY AND SUPPORT

Serene Innovations is committed to customer satisfaction. While specific warranty details are not provided in this manual, the product is backed by Serene Innovations' commitment to quality and service.

For warranty information, technical support, or any questions regarding your Renew Dryer & Freshener, please contact Serene Innovations customer service directly. Contact information can typically be found on the product packaging or the official Serene Innovations website.

