

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Veridian Healthcare](#) /
- › [Veridian Healthcare Ear Thermometer 09-340 Instruction Manual](#)

## Veridian Healthcare 09-340

# Veridian Healthcare Ear Thermometer

**MODEL: 09-340**

### Introduction and Product Overview

---

The Veridian Healthcare 09-340 Ear Thermometer is designed for accurate and convenient body temperature measurement from the ear. This device provides instant 1-second readouts, making it suitable for all ages, including children and sleeping patients. It features a backlit display for easy reading in various lighting conditions and a fever alert icon to quickly indicate elevated temperatures. The thermometer also includes a 10-memory recall function to track temperature changes over time.

This tympanic infrared thermometer is manufactured to meet the standards of healthcare professionals, ensuring reliable measurements for home use. No probe covers are required for operation.



Image: The Veridian Healthcare Ear Thermometer shown alongside its retail packaging, highlighting key features like instant readout and fever alert icons.

## Setup and Initial Use

Before first use, ensure the CR2032 battery is correctly installed. The thermometer comes with one CR2032 battery included. Refer to the battery compartment on the back of the device for proper orientation.

### Included Components:

- Veridian Healthcare Ear Thermometer
- Storage Case
- One CR2032 Battery
- Instruction Manual (English/Spanish)



Image: A visual representation of the items included in the product box, such as the thermometer, storage case, and battery.

## Operating Instructions

Follow these simple steps to take a temperature reading:

1. **Activate:** Press the SCAN button to turn on the ear thermometer. The display will illuminate.
2. **Position:** Gently insert the probe into the ear canal. Ensure proper, yet comfortable, positioning for accurate results. The thermometer is designed for comfortable readings.
3. **Scan:** Press the SCAN button again to initiate the temperature reading. The measurement will be displayed within 1 second.

The backlit display provides convenient lighted readings, ideal for use at night or in low-light conditions. The device also features a fever alert icon: a smiley face indicates no fever, while a frowny face indicates a fever is detected. The thermometer stores up to 10 previous readings for easy tracking.

# COMFORTABLE READINGS



Image: A child receiving a temperature reading with the ear thermometer, demonstrating comfortable usage.

# HOW IT WORKS

MEASUREMENTS ARE  
COMPARABLE TO ORAL READINGS



IN ONE SECOND,  
512 INFRARED ENERGY  
MEASUREMENTS ARE TAKEN,  
THEN CONVERTED INTO A  
SINGLE TEMPERATURE READING.

PROPER POSITIONING OF THIS  
EAR THERMOMETER IN THE EAR  
CANAL PROVIDES SAFE AND  
ACCURATE RESULTS EVERY TIME.

Image: An illustration detailing the infrared measurement process, where 512 infrared energy measurements are taken in one second to provide a single temperature reading.



Image: Close-up of the thermometer display showing the fever alert icons: a happy face for no fever and a sad face for fever.



*Your family's health is important to you.  
A reliable thermometer can help identify fevers early,  
track illnesses and assist in making important health decisions.*

*Infrared thermometer technology offers a safe,  
fast alternative to traditional thermometers,  
with clinically accurate, comparable results in an instant.*



Image: The thermometer's backlit LCD display, emphasizing its visibility for silent readings and night mode.

## **Maintenance and Care**

---

To ensure the longevity and accuracy of your Veridian Healthcare Ear Thermometer, it is recommended to wipe it clean after each use. Use a soft cloth dampened with alcohol or a mild disinfectant. Do not immerse the thermometer in water or any other liquid. Store the thermometer in its provided storage case when not in use to protect it from dust and damage.

## **Troubleshooting**

---

If you encounter any issues with your thermometer, first ensure the battery is properly installed and has sufficient charge. If the display shows an error message or provides inconsistent readings, try replacing the battery. For persistent issues or questions not covered in this manual, please contact Veridian Healthcare customer support.

## **Specifications**

---

Feature	Detail
Model Number	09-340
Product Dimensions	4.5 x 1.5 x 1 inches; 1.6 ounces
Display Type	Digital
Connectivity Technology	Infrared
Special Feature	Alarm
Age Range (Description)	All Ages
Outer Material	Plastic
Specification Met	ASTM
Unit Count	1 Count

## Warranty and Support

---

The Veridian Healthcare Ear Thermometer is backed by a **5-year limited warranty**. For warranty claims, technical support, or any product-related inquiries, please refer to the contact information provided in the full instruction manual or visit the official Veridian Healthcare website.

## Important Safety Information

---

This device is intended for measuring human body temperature from the ear. It is for adult use or use under adult supervision only. Always read the included instructions thoroughly before use to ensure best performance and accurate results. Always discuss changes in your medication or health condition with a healthcare professional. Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.