

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [Homecraft](#) /

› [Homecraft Inflatable Ring Cushion 091090984 User Manual](#)

## Homecraft 091090984

# Homecraft Inflatable Ring Cushion User Manual

Model: 091090984 | Brand: Homecraft

## 1. PRODUCT OVERVIEW

---

The Homecraft Inflatable Ring Cushion is a durable rubber cushion designed to provide pressure relief, particularly for the coccyx area. Its ring shape and contoured, textured surface help prevent slipping and ensure comfort. This manual provides essential information for the proper setup, operation, maintenance, and troubleshooting of your cushion.



Image 1.1: Homecraft Inflatable Ring Cushion, demonstrating its design and inflation valve.

## 2. FEATURES

---

- Made from durable rubber for long-lasting use.
- Equipped with a special valve for easy inflation.
- Adjustable firmness: inflate to desired level using a needle valve adapter and a bicycle-style pump (not included).
- Features a slightly contoured and textured surface to prevent slipping.
- Ring-shaped design specifically for coccyx pressure relief.
- Diameter: 45.7 cm (18 inches).

## 3. SETUP INSTRUCTIONS

---

Follow these steps to prepare your inflatable ring cushion for use:

1. Unpack the cushion from its packaging.
2. Locate the inflation valve on the cushion.
3. Insert a needle valve adapter (not included) into the valve opening.
4. Using a bicycle-style pump (not included), inflate the cushion to your desired firmness. **Do not overinflate**, as this can damage the cushion or reduce comfort.
5. Remove the pump and adapter, ensuring the valve is securely closed to prevent air leakage.

## 4. OPERATING INSTRUCTIONS

---

To use the Homecraft Inflatable Ring Cushion:

1. Place the inflated cushion on a flat, stable surface, such as a chair, wheelchair, or bed.
2. Ensure the ring opening is positioned to relieve pressure from the coccyx or other sensitive areas.
3. Sit gently on the cushion, distributing your weight evenly.
4. Adjust your position as needed for maximum comfort and pressure relief.
5. If the cushion feels too firm or too soft, adjust the air pressure by partially deflating or adding more air as described in the Setup Instructions.

## 5. MAINTENANCE

---

Proper maintenance ensures the longevity and hygiene of your cushion:

- **Cleaning:** Wipe the cushion surface with a damp cloth and mild soap. Do not use abrasive cleaners, harsh chemicals, or solvents, as these can damage the rubber. Rinse thoroughly with clean water and allow to air dry completely before storage or reuse.
- **Storage:** Store the cushion in a cool, dry place away from direct sunlight, extreme temperatures, and sharp objects. For long-term storage, the cushion can be partially deflated to save space.
- **Inspection:** Regularly inspect the cushion for any signs of wear, punctures, or damage. Discontinue use if damage is found to prevent injury or further product degradation.

## 6. TROUBLESHOOTING

---

Refer to the table below for common issues and their solutions:

Issue	Solution
Cushion deflates quickly.	Check the inflation valve to ensure it is securely closed. Inspect the cushion surface thoroughly for any punctures or leaks. If a leak is found, the cushion may need to be replaced.
Difficulty inflating the cushion.	Ensure the needle valve adapter is correctly inserted into the valve. Verify that the pump is functioning properly and providing adequate air pressure.
Cushion feels uncomfortable or provides insufficient support.	Adjust the air pressure to achieve the desired firmness. Ensure the cushion is correctly positioned to target the intended pressure relief area.

## 7. SPECIFICATIONS

---

Key technical specifications for the Homecraft Inflatable Ring Cushion:

Specification	Detail
Brand	Homecraft

Specification	Detail
Model Number	091090984
Color	Orange
Diameter	45.7 cm (18 inches)
Material	Rubber
Features	Inflatable
Item Weight	699 g (1.54 lbs)
Product Dimensions (L x W x H)	45.7 x 45.7 x 0.25 cm (18 x 18 x 0.1 inches)
Shape	Round
Recommended Age Range	Adult
Recommended Product Uses	Medical

## 8. WARRANTY AND SUPPORT

---

Warranty and support information for the Homecraft Inflatable Ring Cushion is not available in the provided product data. For details regarding product warranty, returns, or customer support, please refer to the manufacturer's official website or contact their customer service directly.