

Ken-A-Vision T-17541C

Ken-A-Vision T-17541C Digital CoreScope 2 Compound Microscope

USER MANUAL

1. Introduction

This manual provides essential instructions for the safe and effective use of your Ken-A-Vision T-17541C Digital CoreScope 2 Compound Microscope. Designed for biological and educational applications, this digital microscope features a monocular head, LED illumination, and a 1.3-megapixel camera for capturing still images and video. Please read this manual thoroughly before operating the device and retain it for future reference.

2. Safety Information

- **Electrical Safety:** Ensure the power supply is connected to a grounded outlet. Do not operate the microscope with wet hands or in damp environments. Disconnect power before cleaning or servicing.
- **Handling:** Always carry the microscope by its base and arm. Avoid sudden impacts or vibrations.
- **Optical Care:** Do not touch optical surfaces with bare fingers. Use only approved lens cleaning materials.
- **Ventilation:** Ensure adequate ventilation around the microscope to prevent overheating.
- **Children:** This device is not a toy. Adult supervision is recommended when used by children.

3. Package Contents

Verify that all items listed below are present in your package:

- Ken-A-Vision CoreScope 2 Digital Microscope
- Power Supply
- Applied Vision 4 Software (for Windows, Mac, and Linux)
- Microscope Cover

- Allen Wrench
- Instructions (this manual)

4. Product Overview

The Ken-A-Vision T-17541C CoreScope 2 is a robust digital compound microscope designed for clarity and ease of use. Key components include:

- **Monocular Head:** Features a fixed 45-degree incline for comfortable viewing and 360-degree rotation.
- **Eyepiece:** 10x widefield eyepiece with a pointer for easy observation.
- **Objectives:** Reverse-mounted 4x, 10x, and 40xS DIN achromatic lens objectives, protected against humidity and climate variations.
- **Floating Stage:** Round stage with spring-loaded clips to secure slides and a stop mechanism to prevent damage.
- **Focus Knobs:** Separate coaxial coarse and fine focus controls for precise adjustments.
- **LED Illumination:** Built-in bright and cool white LED light with dimmer and on/off switch.
- **Digital Camera:** Integrated 1.3-megapixel (MP), 720P high-definition (HD) CMOS camera with a built-in USB port for direct connection to a computer or monitor.



Figure 1: Front-side view of the Ken-A-Vision T-17541C Digital CoreScope 2 Compound Microscope, showing the monocular head, objective lenses, floating stage, focus knobs, and connected USB cable.



Figure 2: Angled view of the Ken-A-Vision T-17541C Digital CoreScope 2 Compound Microscope, highlighting the compact design and the position of the objective lenses and stage.

5. Setup

1. **Unpacking:** Carefully remove the microscope and all accessories from the packaging.
2. **Placement:** Place the microscope on a stable, level surface away from direct sunlight, excessive heat, or vibrations.
3. **Power Connection:** Connect the provided power supply to the microscope and then plug it into a standard electrical outlet (120-240V). The microscope can also be powered via a low voltage USB computer connection.
4. **Software Installation:** For digital imaging, install the Applied Vision 4 software on your computer. The software is compatible with Windows (XP SP2 or later, Vista, 7), Mac (OS 10.5 or later), and Linux. Follow the on-screen instructions during installation.
5. **USB Connection:** Connect the microscope's built-in USB port to your computer using a USB cable. This connection enables the digital camera functionality and can also provide power.
6. **Remove Cover:** Remove the microscope cover before use.

6. Operating Instructions

1. **Power On:** Turn on the microscope using the on/off switch for the LED illumination.
2. **Adjust Illumination:** Use the dimmer control to adjust the brightness of the LED light to a comfortable level for viewing your specimen.
3. **Place Specimen:** Place your prepared slide on the floating stage, securing it with the spring-loaded clips.
4. **Select Objective:** Rotate the nosepiece to select the desired objective lens (4x, 10x, or 40xS). Start with the lowest magnification (4x) for initial viewing.
5. **Focusing:**

- Use the **coarse focus knob** (larger knob) to bring the specimen into approximate focus.
 - Then, use the **fine focus knob** (smaller knob) for precise focusing and to achieve a sharp image.
6. **Eyepiece Adjustment:** Adjust the variable diopter on the eyepiece for optimal visual clarity.
7. **Rotate Head:** The monocular head can be rotated 360 degrees for shared viewing or comfortable positioning.
8. **Digital Imaging (with software):**
- Ensure the microscope is connected to your computer via USB and the Applied Vision 4 software is running.
 - The software will display the live feed from the 1.3 MP HD CMOS camera.
 - Use the software interface to capture still images or record video. Refer to the Applied Vision 4 software manual for detailed instructions on its features.
9. **Power Off:** When finished, turn off the LED illumination and disconnect the power supply. Cover the microscope with the provided dust cover.

7. Maintenance

- **Cleaning the Body:** Wipe the microscope body with a soft, damp cloth. Avoid using harsh chemicals or solvents.
- **Cleaning Lenses:** Use a soft lens brush to remove dust. For smudges, use a specialized lens cleaning solution and lens paper. Do not use abrasive materials.
- **Storage:** Always cover the microscope with the dust cover when not in use to protect it from dust and debris. Store in a dry environment.
- **Humidity Control:** The objective lenses are designed with humidity and climate control protection, but prolonged exposure to extreme conditions should be avoided.

8. Troubleshooting

- **No Illumination:** Check if the power supply is properly connected and the on/off switch is engaged. Ensure the dimmer is not set to its lowest setting.
- **Image is Blurry:** Adjust both the coarse and fine focus knobs. Ensure the objective lens is properly clicked into place. Clean the eyepiece and objective lenses if smudges are present.
- **No Image on Computer Screen:** Verify the USB cable is securely connected to both the microscope and the computer. Ensure the Applied Vision 4 software is installed and running correctly. Check your computer's device manager to confirm the camera is recognized.
- **Specimen Not Centered:** Gently adjust the position of the slide on the floating stage.

9. Specifications

Feature	Specification
Model Number	T-17541C
Viewing Configuration	Monocular with 45-degree incline, 360-degree rotation
Eyepiece	10x widefield with pointer
Objective Lenses	4x, 10x, 40xS DIN Achromatic

Feature	Specification
Maximum Magnification	400x
Illumination	Bright and cool LED white light with dimmer and on/off switch
Focus	Coaxial coarse and fine focus
Digital Camera	1.3 MP, 720P HD CMOS, USB 2.0
Software Compatibility	Applied Vision 4 (Windows XP SP2+, Vista, 7; Mac OS 10.5+; Linux)
Stage	Round floating stage with 2 spring-loaded clips
Power Source	Low voltage USB computer connection or external USB power plug, 120-240V
Material	Aluminum
Product Dimensions (H x W x D)	40.0 x 30.48 x 22.86 cm (15.75 x 12 x 9 inches)
Weight	2.31 kg (5.1 lb)
Certificates	ISO:9001 – 2000, CE, CSA, RoHS

10. Warranty and Support

The Ken-A-Vision T-17541C Digital CoreScope 2 Compound Microscope is manufactured by Ken-A-Vision Manufacturing, headquartered in Kansas City, MO. For warranty information, technical support, or service inquiries, please refer to the official Ken-A-Vision website or contact their customer service department directly. Please have your model number (T-17541C) and purchase date available when contacting support.