

Ken-A-Vision T-29033

Ken-A-Vision T-29033 Research Microscope User Manual

Model: T-29033



1. INTRODUCTION

This manual provides comprehensive instructions for the setup, operation, and maintenance of your Ken-A-Vision T-29033 Research Microscope. Designed for advanced high school, college, university, and clinical laboratory settings, this microscope offers durable construction and high-quality optical performance. Please read this manual thoroughly before using the instrument to ensure proper function and longevity.

2. SAFETY INFORMATION

Observe the following safety precautions to prevent damage to the microscope and ensure user safety:

- Always handle the microscope with care. Avoid sudden impacts or vibrations.
- Connect the microscope to a grounded power outlet.
- Do not attempt to disassemble the optical or electrical components. Refer servicing to qualified personnel.
- Keep the microscope away from direct sunlight, high temperatures, dust, and corrosive chemicals.
- Unplug the power cord before cleaning or moving the microscope.

3. PACKAGE CONTENTS

Verify that all components are present upon unpacking:

- Ken-A-Vision T-29033 Research Microscope main body
- Binocular head with 10x 20mm eyepieces (2)
- Objective lenses: 4x, 10x, 40xS, 100xS DIN plan objectives
- Power cord
- Dust cover
- User Manual (this document)

4. COMPONENT IDENTIFICATION

Familiarize yourself with the various parts of the Ken-A-Vision T-29033 Research Microscope:



Figure 1: Ken-A-Vision T-29033 Research Microscope. This image displays the complete microscope from a front-right perspective, showing the binocular eyepieces, revolving nosepiece with objective lenses, mechanical stage, and focusing knobs on the side.

1. **Eyepieces:** 10x 20mm, with eyeglass guards and variable adjustment for eye strengths.
2. **Binocular Head:** 30-degree incline, 360-degree rotation.

3. **Revolving Nosepiece:** Holds the objective lenses.
4. **Objective Lenses:** 4x, 10x, 40xS, 100xS DIN plan objectives.
5. **Mechanical Stage:** 180mm x 150mm, with scratch-resistant coating and vernier scales for precise specimen movement.
6. **Condenser:** N.A. 1.25, with rack and pinion adjustment, filter holder.
7. **Coaxial Coarse and Fine Focus Knobs:** For adjusting focus.
8. **Illumination System:** Halogen tungsten bulb (12V 20W) with variable brightness, Kohler illumination.
9. **Base:** Stable aluminum single cast body.

5. SETUP

1. **Unpacking:** Carefully remove the microscope and all accessories from the packaging. Retain the original packaging for future transport or storage.
2. **Placement:** Place the microscope on a stable, level surface, away from direct sunlight, heat sources, and vibrations.
3. **Head Installation:** If the binocular head is separate, gently insert it into the top of the microscope body and secure it with the exposed thumb screw.
4. **Eyepiece Insertion:** Insert the two 10x 20mm eyepieces into the eyepiece tubes of the binocular head.
5. **Power Connection:** Connect the power cord to the microscope's power input and then to a grounded electrical outlet (AC 85V - 240V).

6. OPERATING INSTRUCTIONS

1. **Power On:** Flip the power switch located on the base of the microscope. Adjust the illumination intensity using the variable brightness control knob.
2. **Specimen Placement:** Place a prepared microscope slide onto the mechanical stage. Use the stage clips to hold it securely.
3. **Stage Movement:** Use the X-Y stage control knobs to position the specimen directly under the objective lens.
4. **Objective Selection:** Rotate the revolving nosepiece to select the lowest power objective (e.g., 4x). Ensure the objective clicks firmly into place.
5. **Focusing:**
 - Look through the eyepieces.
 - Use the **coarse focus knob** to bring the specimen into approximate focus.
 - Use the **fine focus knob** for precise focusing and to achieve a sharp image.
6. **Interpupillary Distance Adjustment:** Adjust the distance between the eyepieces by gently pushing or pulling them apart until you see a single, circular field of view.
7. **Diopter Adjustment:** If one eye is clearer than the other, use the diopter adjustment ring on one eyepiece to compensate for differences in vision.
8. **Higher Magnification:** Once focused at low power, rotate the nosepiece to a higher power objective (e.g., 10x, 40xS). The microscope is parfocal, meaning only minor adjustments with the fine focus knob should be needed.
9. **Using the 100xS Oil Immersion Objective:**
 - After focusing with the 40xS objective, rotate the nosepiece halfway between the 40xS and 100xS objectives.

- Place a small drop of immersion oil directly onto the specimen slide where the light passes through.
- Rotate the 100xS objective into the oil drop. The objective should be immersed in the oil.
- Use only the fine focus knob to bring the image into sharp focus. Do not use the coarse focus knob with the 100xS objective.
- After use, clean the 100xS objective and the slide immediately with lens paper and a suitable lens cleaning solution to remove all immersion oil.

10. **Condenser and Iris Diaphragm Adjustment:** Adjust the condenser height and the iris diaphragm to optimize contrast and resolution for each objective. Generally, higher magnification requires more light and a wider diaphragm opening.

7. MAINTENANCE

- **Cleaning Optics:** Use only specialized lens paper and lens cleaning solution for objective lenses and eyepieces. Gently wipe in a circular motion. Never touch optical surfaces with bare fingers.
- **Cleaning Body:** Wipe the microscope body with a soft, damp cloth. Avoid using harsh chemicals or solvents.
- **Dust Protection:** Always cover the microscope with the provided dust cover when not in use.
- **Storage:** Store the microscope in a cool, dry, and dust-free environment. Ensure the lowest power objective is in place and the stage is lowered.
- **Bulb Replacement:** If the halogen bulb burns out, unplug the microscope and allow it to cool. Refer to the diagram in the bulb compartment for replacement instructions. Use only the specified 12V 20W halogen bulb.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
No illumination	Power cord disconnected, bulb burnt out, power switch off, brightness control set to minimum.	Check power connection, replace bulb, turn on power switch, adjust brightness control.
Image unclear/blurry	Incorrect focus, objective lens dirty, eyepiece dirty, immersion oil on dry objective, slide upside down.	Adjust fine focus, clean lenses, ensure correct objective use (oil for 100xS only), check slide orientation.
Dark field of view	Iris diaphragm closed too much, condenser too low, light intensity too low, objective not fully engaged.	Open iris diaphragm, raise condenser, increase light intensity, rotate nosepiece until objective clicks.
Dust or spots in field of view	Dust on eyepiece, objective, or slide.	Clean eyepieces, objectives, and slides with appropriate cleaning materials. Rotate eyepiece to determine if dust is on it.

9. SPECIFICATIONS

Model Number	T-29033
Eyepieces	2 - 10x 20mm with eyeglass guards, eyepiece pointer
Head	Binocular, 30-degree incline, 360-degree rotation

Objective Lenses	4x, 10x, 40xS, 100xS DIN plan objectives
Nosepiece	Reversed 4-hole with ball bearing action stop
Stage	Vernier mechanical stage, 180mm x 150mm, scratch-resistant coating
Condenser	Factory centered N.A. 1.25 with rack and pinion, filter holder
Focus	Coaxial coarse and fine adjustment
Illumination	Halogen tungsten bulb, 12V 20W, AC 85V - 240V, variable brightness, Kohler illumination
Frame Material	Aluminum single cast body
Certifications	ISO:9001-2000, CE, CSA, RoHS
Dimensions	20.5 x 14.4 x 11.1 inches
Weight	18.25 Pounds

10. WARRANTY INFORMATION

The Ken-A-Vision T-29033 Research Microscope is covered by a **10-year limited warranty**. This warranty covers defects in materials and workmanship under normal use. Please refer to the warranty card included with your product for specific terms, conditions, and limitations. Keep your proof of purchase for warranty claims.

11. SUPPORT

For technical assistance, replacement parts, or service inquiries, please contact Ken-A-Vision customer support. When contacting support, please have your microscope model number (T-29033) and serial number (if applicable) ready.

For additional resources and product information, visit the official Ken-A-Vision website.



Related Documents - T-29033



Ken-A-Vision 910-171-210 Furniture Camera: Parts List and Assembly Instructions

Detailed parts list and step-by-step assembly instructions for the Ken-A-Vision 910-171-210 Furniture Camera system. Includes component identification and installation guide.



Excelitas A-Zoom Micro Microscope: Reference Manual

Comprehensive reference manual for the Excelitas A-Zoom Micro Microscope, covering optical, mechanical, and operational instructions, safety guidelines, features, specifications, mounting, power connection, objective integration, head options, polarizers, servicing, and warranty information.



HIKMICRO HABROK Series Quick Start Guide

Get started quickly with the HIKMICRO HABROK Series multi-spectrum binoculars. This guide provides essential information on setup, features, and operation for thermal and digital day & night vision.




ACCU-SCOPE EXM-150 Darkfield Annulus: Installation and Operation Guide

Detailed instructions for installing and operating the ACCU-SCOPE EXM-150 Darkfield Annulus, explaining darkfield illumination principles and usage with EXM-150 series microscopes.



HIKMICRO HABROK 4K Series Multi-Spectrum Binocular Quick Start Guide

Your essential guide to the HIKMICRO HABROK 4K Series Multi-Spectrum Binocular. Learn about setup, features, and operation for this advanced thermal and digital imaging device.



AmScope SM-1 Series Stereo Microscope User Manual

Comprehensive user manual for the AmScope SM-1 Series stereo microscopes, covering setup, operation, specifications, and troubleshooting for models like SM-1B/T, SM-1TS/BS, and SM-1(B/T)-PL.

