

[Manuals.plus](#) /

› [Easy Touch](#) /

› Easy Touch Glucose Test Strips (50 Count) Instruction Manual

Easy Touch 807050

Easy Touch Glucose Test Strips (50 Count) Instruction Manual

INTRODUCTION

This manual provides detailed instructions for the proper use and care of your Easy Touch Glucose Test Strips. These strips are designed for self-testing blood glucose levels at home or on the go, in conjunction with the Easy Touch Glucose Monitoring System. Please read this manual thoroughly before use to ensure accurate results and safe operation.

PRODUCT OVERVIEW

The Easy Touch Glucose Test Strips are an essential component of the Easy Touch Glucose Monitoring System. Each pack contains 50 test strips.



Image showing the packaging for Easy Touch Glucose Test Strips, 50 count. The box features the Easy Touch logo and indicates 'Glucose Test Strips'.

Key Features

- Ideal for self-testing at home or on the go, requiring only a tiny blood sample.
- Delivers results in just 5 seconds.
- Compatible exclusively with the Easy Touch Glucose Monitoring System.
- Each pack includes 50 test strips.

SETUP

The Easy Touch Glucose Test Strips are designed for use with the Easy Touch Glucose Monitoring System. Ensure you have the Easy Touch meter, a lancing device, and lancets ready.

Video demonstrating the components of the Easy Touch Glucose Monitoring System kit, including the meter, lancing device, lancets, batteries, carrying case, and user guides. It also shows how to insert a test strip into the meter.

Preparing the Lancing Device

Refer to your Easy Touch Lancing Device manual for specific instructions on loading a lancet and adjusting the penetration depth. The kit includes 10 twist lancets.

How to use your Test Strip

Apply Blood Sample

Insert Into Meter



A visual guide illustrating the two main steps for using a test strip: inserting it into the meter and then applying a blood sample to the designated area.

OPERATING INSTRUCTIONS

Performing a Blood Glucose Test

1. **Prepare the Meter:** Insert an Easy Touch Test Strip into the test strip port of your Easy Touch Glucose Monitoring System meter. The meter will power on automatically.
2. **Obtain a Blood Sample:** Use the lancing device to obtain a small blood sample (0.5 μ L). For greater comfort, alternate between testing your palms, forearms, and fingertips.
3. **Apply Blood Sample:** Touch the tip of the test strip to the blood sample. The strip will draw the blood in automatically.
4. **Read Results:** The meter will display your blood glucose result in just 5 seconds. No coding is required.

How to use your meter

3 Easy Steps!



**Insert Test Strip,
Meter will Turn On**



**Apply Tiny Blood
Sample (0.5 μ L)**



**Accurate Results in
Just 5 Seconds**

A three-step diagram showing how to use the Easy Touch glucose meter: 1) Insert test strip, meter turns on. 2) Apply tiny blood sample (0.5 μ L). 3) Accurate results in just 5 seconds.

MAINTENANCE

Test Strip Storage

Store Easy Touch Glucose Test Strips in their original vial with the cap tightly closed. Keep them in a cool, dry place between 39°F and 86°F (4°C and 30°C). Do not refrigerate or freeze. Discard any remaining test strips 3 months after opening the vial or on the expiration date, whichever comes first.



A close-up image of an Easy Touch Glucose Test Strips vial, showing the label with storage instructions and warnings.

Meter Cleaning

Clean the exterior of your Easy Touch meter with a soft, damp cloth. Do not immerse the meter in water or use harsh cleaning agents.

Control Solution Testing

Use Easy Touch Control Solution to check your meter's and test strips' performance. Perform a control solution test when you open a new vial of test strips, suspect the meter or strips are not working correctly, or if your test results do not match how you feel. Control solution is sold separately.

TROUBLESHOOTING

If you encounter issues with your Easy Touch Glucose Test Strips or meter, consider the following:

- **Incorrect Results:** Ensure the test strip is not expired and the vial has been properly sealed. Verify that the blood sample size is adequate (0.5 μ L). Perform a control solution test to check meter and strip accuracy.
- **Meter Not Turning On:** Check if the batteries are correctly inserted and have sufficient charge. Replace batteries if necessary.

- **Error Messages:** Refer to your Easy Touch Glucose Monitoring System user manual for specific error code explanations and solutions.

For further assistance, contact MHC Medical Products customer support.

SPECIFICATIONS

Feature	Detail
Product Name	Easy Touch Glucose Test Strips
Model Number	807050
Sample Volume	0.5 Microliters
Test Time	5 seconds
Number of Strips	50 per vial
Compatibility	Easy Touch Glucose Monitoring System only
Storage Temperature	39°F - 86°F (4°C - 30°C)
UPC	084960807053

Note: Specifications are subject to change without notice.

IMPORTANT SAFETY INFORMATION

Always read the entire user manual for your Easy Touch Glucose Monitoring System before performing any tests. Consult with your healthcare professional for guidance on managing your diabetes and interpreting your blood glucose results. Do not make any medical decisions without consulting your doctor.

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

WARRANTY AND SUPPORT

MHC Medical Products stands behind its promise to deliver quality products. For warranty information or customer support, please refer to the contact details provided in your Easy Touch Glucose Monitoring System kit or visit the official Easy Touch website.

If you're not satisfied, let us know, and we'll make it right.

For visual setup assistance, search for 'Easy Touch Glucose Meter Setup' on [YouTube](#).

Visit the official [Easy Touch Store](#) for more products and support.