

Omron HEM-7200

Omron HEM-7200 Blood Pressure Monitor User Manual

Model: HEM-7200

1. INTRODUCTION

This manual provides essential instructions for the safe and effective operation of your Omron HEM-7200 Blood Pressure Monitor. Please read all instructions thoroughly before use and retain this manual for future reference. The Omron HEM-7200 is designed for home use to measure systolic and diastolic blood pressure and pulse rate. It features advanced averaging technology, irregular heartbeat detection, and body movement detection for accurate readings.



Figure 1: The Omron HEM-7200 Blood Pressure Monitor. This image shows the main unit with its digital LCD display showing blood pressure readings (systolic 118 mmHg, diastolic 78 mmHg) and pulse rate (70 bpm). An arm cuff is visible to the left of the main unit.

2. WHAT'S IN THE BOX

Verify that all components are present in the package:

- Main unit (Omron HEM-7200 Blood Pressure Monitor)
- Medium size arm cuff (fits arm circumference 22-32 cm)
- Storage case
- Instruction manual
- 4 AA batteries

An AC adapter (6V....4W) is optional and not included.

3. SETUP

1. **Install Batteries:** Open the battery compartment on the back of the main unit. Insert the 4 AA batteries, ensuring correct polarity (+/-). Close the compartment cover.
2. **Attach Arm Cuff:** Connect the air plug of the arm cuff to the air jack on the left side of the main unit. Ensure a secure connection.
3. **Set Date and Time:** Refer to the detailed instructions in the full manual for setting the date and time, which is important for accurate memory logging.

4. OPERATING INSTRUCTIONS

4.1 Taking a Measurement

1. **Prepare for Measurement:** Sit comfortably with your back supported and feet flat on the floor. Rest your arm on a table so the cuff is at heart level. Avoid talking or moving during measurement.
2. **Apply the Arm Cuff:** Wrap the cuff firmly around your upper arm, approximately 1-2 cm above your elbow. Ensure the air tube is centered on the inside of your arm and points towards your hand. The cuff should fit an arm circumference of 22-32 cm.
3. **Start Measurement:** Press the **START/STOP** button. The cuff will begin to inflate automatically.
4. **Read Results:** Once the measurement is complete (approximately 45 seconds), your systolic pressure, diastolic pressure, and pulse rate will be displayed on the LCD screen.
5. **Power Off:** Press the **START/STOP** button again to turn off the unit, or it will automatically power off after a period of inactivity.

4.2 Advanced Features

- **Advanced Averaging Technology:** The monitor automatically averages the last 3 readings taken within a 10-minute period to provide a more accurate representation of your blood pressure.
- **Irregular Heartbeat Indicator:** If an irregular heartbeat is detected during measurement, a specific symbol will appear on the display. If this symbol appears frequently, consult your doctor.
- **Blood Pressure Level Indicator:** The device includes a hypertension indicator (a digital heart symbol) that blinks continuously if your blood pressure level is outside the standard range (above 135 mmHg systolic or 85 mmHg diastolic).
- **Body Movement Detection:** The monitor detects body movement during measurement, which can affect accuracy. If movement is detected, the device will indicate this, prompting a re-measurement if necessary.
- **Memory Function:** The HEM-7200 stores up to 60 sets of readings with the date and time. To review past readings, use the memory buttons on the unit.

5. MAINTENANCE

- **Cleaning:** Wipe the main unit and cuff with a soft, dry cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Storage:** Store the monitor and cuff in the provided storage case in a clean, dry place, away from direct sunlight and extreme temperatures.
- **Battery Replacement:** Replace all four AA batteries when the low battery indicator appears on the display.
- **Cuff Care:** Ensure the cuff is not excessively creased or damaged. The cuff is designed for repeated use but should be inspected periodically.

6. TROUBLESHOOTING

If you encounter issues with your Omron HEM-7200, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
No display when powered on	Batteries are dead or incorrectly installed.	Replace all batteries with new ones, ensuring correct polarity.
Cuff does not inflate or inflates too little	Air plug is not securely connected; air leakage from cuff.	Ensure the air plug is fully inserted. Check cuff for damage.
Inaccurate or inconsistent readings	Incorrect cuff placement; body movement during measurement; talking during measurement; incorrect posture.	Re-read "4.1 Taking a Measurement" and ensure proper technique. Avoid movement and talking. Ensure cuff is at heart level.
Irregular Heartbeat Symbol appears frequently	Indicates frequent irregular heartbeats.	Consult your doctor for medical advice.

If the problem persists after attempting these solutions, contact Omron customer support.

7. SPECIFICATIONS

Feature	Detail
Model	HEM-7200
Display	Digital display
Measurement Method	Oscillometric method
Measurement Range (Pressure)	0 – 299 mmHg
Measurement Range (Pulse)	40 – 180 beats/min
Accuracy (Pressure)	±3 mmHg
Accuracy (Pulse Rate)	±5% of reading
Inflation	Fuzzy-logic controlled by electric pump

Feature	Detail
Deflation	Automatic pressure release valve
Pressure Detection	Capacitive pressure sensor
Power Supply	4 AA batteries 1.5V or AC/DC adaptor (optional, 6V....4W)
Battery Life	Approx. 1500 measurements with new alkaline batteries
Operating Temperature/Humidity	+10°C to +40°C, 30% to 85% RH
Storage Temperature/Humidity	-20°C to +60°C, 10% to 95% RH
Weight of Main Unit	Approx. 340g without batteries
Main Unit Dimensions	123 (w) x 85 (h) x 141 (l) mm
Cuff Size	Approx. 146mm x 446mm (Medium cuff: fits arm circumference 22-32cm)
Cuff Material	Nylon & polyester

8. IMPORTANT HEALTH INFORMATION

Understanding your blood pressure and monitoring it regularly is crucial for health management.

8.1 Why Measure Blood Pressure at Home?

- **White-coat Hypertension:** Readings may be high in a doctor's office but lower at home due to anxiety.
- **Masked Hypertension:** Readings may be low in a doctor's office but higher at home.
- **Morning Hypertension:** Blood pressure can be higher in the morning.
- Home monitoring allows for measurements in a relaxed state, providing a more accurate picture over time.

8.2 Understanding High Blood Pressure

Normal blood pressure (normotension) is generally defined as less than 130 mmHg systolic and less than 85 mmHg diastolic. Elevated blood pressure places increased strain on blood vessels, potentially leading to damage. This can result in serious conditions such as stroke (internal bleeding in the brain) or ischemic heart disease (e.g., heart attack, angina) due to increased workload on the heart. Maintaining blood pressure within a normal range is vital for preventing these health risks.

8.3 Helpful Blood Pressure Reminders

- Regularly have your blood pressure checked by a healthcare professional.
- Follow your doctor's advice regarding diet and lifestyle changes if you have high blood pressure.
- If medication is prescribed, take it regularly as directed. Report any new symptoms to your doctor.
- Measure and record your blood pressure at home consistently. This data is valuable for your doctor in evaluating your condition and treatment effectiveness.
- Collaborate with your healthcare team to manage your blood pressure effectively.

9. WARRANTY AND SUPPORT


The Omron HEM-7200 Blood Pressure Monitor comes with a **1-year manufacturer's warranty** from the date of purchase. For warranty claims, technical support, or further assistance, please refer to the contact information provided in the full

instruction manual or visit the official Omron website.



Related Documents - HEM-7200


	OMRON BP7360 7 Series Upper Arm Blood Pressure Monitor Instruction Manual Comprehensive instruction manual for the OMRON BP7360 7 Series Upper Arm Blood Pressure Monitor. Learn about setup, usage, safety, troubleshooting, and features like AFib detection for home health monitoring.
	OMRON HEM-6181 Wrist Blood Pressure Monitor: User Guide & Safety Instructions Comprehensive instruction manual for the OMRON HEM-6181 Wrist Blood Pressure Monitor. Learn how to accurately measure blood pressure and pulse rate at home, understand readings, and troubleshoot common issues.
	OMRON JPN1 Automatic Blood Pressure Monitor Instruction Manual Instruction manual for the OMRON JPN1 Automatic Blood Pressure Monitor, covering setup, usage, troubleshooting, and technical specifications.
	OMRON M7 Intelli IT & X7 Smart Blood Pressure Monitors: Features, Specifications, and Usage Guide Comprehensive information on OMRON's M7 Intelli IT and X7 Smart automatic upper arm blood pressure monitors, covering features, technical specifications, usage instructions, safety warnings, and maintenance.
	OMRON JPN610T Automatic Blood Pressure Monitor User Manual Comprehensive user manual for the OMRON JPN610T Automatic Blood Pressure Monitor, covering setup, operation, safety precautions, troubleshooting, and specifications. Learn how to accurately measure your blood pressure at home.
	OMRON X2 Smart Blood Pressure Monitor: Accurate Home Health Monitoring The OMRON X2 Smart Automatic Upper Arm Blood Pressure Monitor offers simple, accurate blood pressure and pulse rate measurement at home, featuring Bluetooth connectivity for easy health tracking.



[\[pdf\]](#) User Manual Quick Start Guide Instructions Guide

5 series Blood Pressure Monitor Omron Healthcare Instruction Manual printed in English and Spanish Quick Start Guide The BP742 IntelliSense® Automatic is intended for home use Please read this instruction manual thoroughly before using the unit keep future reference For specific information about your own blood pressure CONSULT YOUR DOCTOR ENGIM 02032012 omronhealthcare ca INSTRUCTION MANUAL 5 series Blood Pressure Monitor Model BP742 ENGLISH ESPAOL tTaAbBILeE oOfF cCO ... ain Unit Dimensions: Cuff Size: Cuff Circumference: Memory: Contents: Optional Accessories: BP742 **HEM-7200-Z** LCD Digital Display Pressure: 0 to 299 mmHg, Pulse: 40 to 180/min. Pressure: 3mmHg or 2 ...

lang:en score:43 filesize: 7.56 M page_count: 30 document date: 2012-01-31




[\[pdf\]](#) User Manual Quick Start Guide Instructions Guide

110114 HEM 7200 BP742CAN IM EN r1 indd administrator Omron Healthcare 5 series Blood Pressure Monitor Instruction Manual printed in English and French Quick Start Guide The Automatic is intended for home use Please read this instruction manual thoroughly before using the unit keep future reference For specific information about your own blood pressure CONSULT YOUR DOCTOR comes with following components Standard Adult Arm Cuff 4 AA Alkaline Batteries SAFETY INFORMATION This device measuring pulse rate adult population Do not on infants or persons who cannot express their intentions ENG 01272011 omronhealthcare ca

INSTRUCTION MANUAL 5 series Blood Pressure Monitor Model BP742CAN ENGLISH TAABBLE OF CONTENTS Bef ... Unit Dimensions: Cuff Size: Cuff Circumference:: Memory: Contents: Optional Accessories: BP742CAN **HEM-7200-CA** LCD Digital Display Pressure: 0 to 299 mmHg, Pulse: 40 to 180/min. Pressure: 3mmHg or 2...


lang:en score:39 filesize: 7.7 M page_count: 27 document date: 2011-01-27



[\[pdf\]](#) User Manual

untitled Model M3 Instruction Manual 6433000 archive org manualzilla id Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto 617-0002 JAPAN EU-r ... essure Monitor Model M3 Instruction Manual English Franais Deutsch Italiano Espaol Nederlands IM-**HEM-7200-E2** V -01-08/2011 5333267-0A 1 F L c d 2 b J a Contents Thank you for purchasing the...

lang:en score:38 filesize: 848.57 K page_count: 30 document date: 2013-09-20



[\[pdf\]](#) User Manual Quick Start Guide Instructions Guide

321727 2 HEM 7200 Z B M En 5 series CAUTION Blood Pressure Monitor SAFETY INFORMATION BEFORE TAKING A manual It is mainly designed for general household use Please read the Safety Information in this instruction before using unit To assure correct of product basic safety measures should always be followed including warnings and cautions listed OPERATING THE DEVICE INSTRUCTION MANUAL Model BP742 INTRODUCTION Thank you purchasing OMRON® Automatic IM omronhealthcare ca

INSTRUCTION MANUAL 5 series Blood Pressure Monitor Model BP742 INTRODUCTION Thank you for purchasin ... free: FOR CUSTOMER SERVICE www.omronhealthcare.com 1-800-634-4350 SPECIFICATIONS Model: BP742 **HEM-7200-Z** Display: LCD Digital Display Measurement Range: Pressure: 0 to 299 mmHg, Pulse: 40 to...

lang:en score:36 filesize: 4.24 M page_count: 2 document date: 2013-01-15

[\[pdf\]](#) User Manual

Model M3 Instruction Manual Medaval Contents English Français Deutsch Italiano Español Nederlands
IM HEM 7200 E 05 11 2010 5329491 4C Manufacturer OMRON HEALTHCARE Co Ltd 24 Yamanouchi
Yamanoshita cho Ukyo ku Kyoto 615 0084 JAPAN EU representative EUROPE B V Kruisweg 577 2132
NA Hoofddorp THE NETHERLANDS omron healthcare Production Проверьте следующие компоненты!
M01 100917 Digital Automatic Blood Pressure Monitor The is a compact fully automatic blood pressure
monitor operating on the devices uses its advanced IntelliSense technology Intended Use Pages from A
M main Please read this instruction manual thoroughly before using unit keep for future reference For
specific information about your own CONSULT YOUR Omron EN medaval ie docs manuals s cdnmpro
129501539 content Manuale

Manufacturer OMRON HEALTHCARE Co., Ltd. 24, Yamanouchi Yamanoshita-cho,
Ukyo-ku, Kyoto 615-0084 JAP ... ensure Monitor Model M3 Instruction Manual English
Franais Deutsch Italiano Espaol Nederlands IM-**HEM-7200**-E-05-11/2010 5329491-
4C 1 F L c d 2 b J a Contents Thank you for purchasing the OMR...
lang:en **score:36** filesize: 1.39 M page_count: 26 document date: 2011-04-15





White Medical
Meranti Lodge
Hillmorton Lane
Clifton upon Dunsmore
Rugby, CV23 0BA

May 2013
All prices exclude VAT
FREE standard delivery
on ALL orders


t: 01788 553904 f: 01788 560820
enquiries@white-medical.co.uk
www.white-medical.co.uk



White Medical is the main Omron supplier into the NHS


[\[pdf\]](#) User Manual Borchure

A D AD 6121A Instruction manual 7333244 archive org manualzilla id
May 2013 All prices exclude VAT FREE standard delivery on ALL orders White
Medical Meranti Lodge Hi ... enquiries white-medical.co.uk www.white-medical.co.uk
4 Omron Blood Pressure Monitors Omron M3 **HEM-7200**-E8 V 36.00 VAT This fully
automatic blood pressure device provides an easy way to check fo...
lang:en **score:35** filesize: 1.84 M page_count: 48 document date: 2013-05-03



[OMRON JPN1 Automatic Blood Pressure Monitor Instruction Manual](#)

Instruction manual for the OMRON JPN1 Automatic Blood Pressure Monitor, covering
setup, usage, troubleshooting, and technical specifications.
lang:en **score:34** filesize: 2.09 M page_count: 1 document date: 2023-10-18




[\[pdf\]](#) User Manual Quick Start Guide Instructions Guide

301671 HEM 7200 Z BP742 A IM WEB indd Omron Healthcare 5 series Blood Pressure Monitor Instruction Manual printed in English and Spanish Quick Start Guide The IntelliSense® Automatic is intended for home use Your new digital blood pressure monitor uses the oscillometric method of measurement This means detects your s movement through 9 DISPLAY SYMBOLS American Heart Association recommends following guideline upper limit normal readings a general because depends on age morbidity Standard Adult Arm Cuff 4 AA Alkaline Batteries BEFORE TAKING MEASUREMENT To ensure reliable reading follow these recommendations 1 Avoid eating drinking alcohol smoking exercising bathing 30 minutes before taking INSTRUCTION MANUAL Model Read Special Conditions section instruction manual if systolic known to bp742 im eng 09152010 dfb5985b b0e4 4286 a7eb 7517669b72f7 omronhealthcare ca stream manualsonline d 80 media bp

INSTRUCTION MANUAL 5 series Blood Pressure Monitor Model BP742 ENGLISH
ESPAOL TABLE OF CONTENTS Before using the Monitor Introduction

..... 3 Safety Information

lang:en score:32 filesize: 1.93 M page_count: 28 document date: 2010-09-10




[\[pdf\]](#) User Manual

untitled Model M3W Instruction Manual IM HEM 7202 E V 02 01 2012 IT 2739 6640 hospitaltrade it allegatiART |||

Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto, 617-0002 JAPAN EU- ... ore di pressione automatico digitale Espoal Modello M3W Nederlands Manuale di istruzioni i IM-**HEM-7200**-E-01-10/09 5323430-0A Sommario Grazie per aver acquistato il misuratore di pressione aut...

lang:en score:32 filesize: 792.73 K page_count: 30 document date: 2013-09-20



[\[pdf\]](#) User Manual

Full name medidores monitores de presion arterial tensiometros digitales hem7200 omron manual espanol viaindustrial manuales |||

Medidores monitores de presin arterial tensimetros digitales HEM 7200 OMRON Manual Espoal www.viaindustrial.com Medidores monitores de presin arterial tensimetros digitales HEM 7200 OMRON Manual Espoal www.viaindustrial.com Medidores monitores de presin arterial tensimetros digitales HEM 7200 OMRO...

lang:es score:31 filesize: 2.53 M page_count: 14 document date: 2017-03-08



[\[pdf\]](#) User Manual Guide

Automatic Blood Pressure Monitor Brijlal Main unit M size cuff storage case Instruction Manual 4 AA batteries Optional AC adaptor R MEDIUM CUFF CR24 LARGE CL24 2371340 familybrands co nz restricted

DELUXE Automatic Blood Pressure Monitor Blood Pressure Level Indicator Detects irregular heartbeat ... y applies the right amount of pressure for fast, accurate and more comfortable measurements Omron **HEM-7200** Intellisense Blood Pressure Monitor provides easy and comfortable blood pressure measuremen...

lang:en score:31 filesize: 1.18 M page_count: 2 document date: 2009-10-20

[\[pdf\]](#) User Manual Catalog

Manual de instrucciones TENSIMETRO DIGITAL M3 Ortopedia Online tienda en Madrid España IM HEM 7200 E2 V 01 08 2011 ES omron healthcare data catalog 3 4431 1 |||

Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto 617-0002 JAPAN EU-r ... essure Monitor Model M3 Instruction Manual English Franais Deutsch Italiano Espaol Nederlands IM-**HEM-7200**-E2 V -01-08/2011 5333267-0A 1 F L c d 2 b J a Contents Thank you for purchasing the...

lang:en score:30 filesize: 1.78 M page_count: 30 document date: 2013-03-22



[\[pdf\]](#) User Manual

HEM 7200 E V CoverAside main D M02 20120208 SANTÉ FRANCE Digital Automatic Blood Pressure Monitor Model Manual ish Français Deutsch Italiano Español Nederlands 0 30 9B ий Check following components! Vérifier les composants suivants ! Prüfen Sie folgende Teile des Lieferumfangs! Controllare i componenti indicati di seguito! ¡Compruebe los siguientes componentes! Controleer de volgende onderdelen! Проверьте следующие Omron E8 medaval ie docs manuals

Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto, 617-0002 JAPAN EU- ... essure Monitor Model M3 Instruction Manual English Franais Deutsch Italiano Espaol Nederlands IM-**HEM-7200**-E V -05-11/2011 5329530-9D English Franais Digital AultnosmtrautMiccotBidoleonloMMd aP3n... lang:en score:28 filesize: 1.43 M page_count: 26 document date: 2012-07-19





[\[pdf\]](#) User Manual

311034 HEM 7200BR IM r5 indd Mona Manual 7200 CPAP Med MANUAL DE INSTRUÇÕES 5332859 2A BR CONTEÚDO Obrigado por adquirir o Monitor de Pressão Arterial Automático OMRON 111760774 img americanas br produtos 01 02 manual 80 submarino manuais

MANUAL DE INSTRUES Monitor de Presso Arterial Automtico Modelo**HEM-7200** 5332859-2A **HEM-7200**-BR CONTEDO Obrigado por adquirir o Monitor de Presso Arterial Automtico OMRON **HEM-7200**. O OMRON**HEM-7200** um monitor de presso arterial automtico completamente automtico, compacto, que opera conforme o princp...

lang:pt score:28 filesize: 2.07 M page_count: 2 document date: 2011-12-14






[pdf] User Manual Specifications Guide Accessories

Provided by ManualsAndMore Download User Manual Omron M6W Information Sheet manualsandmore

|||

-NEW PRODUCT INFORMATION SHEETS- Specifications subject to change
OMRON HEALTHCARE EUROPE B.V. Sale ... Approx. Weight 340 g 9.5 Kg Approx.
Dimensions W x D x H 123 x 141 x 85 587 x 292 x 216 GP-**HEM-7200**-E2 V -03-
10/2011 Part 5 Instruction Manual Design GP-HEM-7213-E V -01-11/2011 Part 6
Te...

lang:en score:27 filesize: 769.45 K page_count: 13 document date: 2019-07-23




[pdf] User Manual

Untitled OMRON HEM 7200 SPHYGMOMANOMETER MANUAL manuallib 12 |||

...

lang:en score:26 filesize: 986.42 K page_count: 2 document date: 2011-11-18



[pdf] User Manual

untitled Click here to the Device Manual Omron M3W HEM 7202 E V Medaval medaval ie docs manuals

|||

Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko,
Kyoto, 617-0002 JAPAN EU- ... omtre automatique Modle M3W Mode d emploi
English Franais Deutsch Italiano Espao! Nederlands i IM-**HEM-7200**-E-01-10/09
5323430-0A Table des matires Merci d avoir fait l acquisition du tensiomtre a...
lang:en score:25 filesize: 774.66 K page_count: 30 document date: 2013-09-20



[pdf] User Manual

HEM 7200 C m healthcare omron co jp 80 support manual

...

lang:en score:24 filesize: 1.34 M page_count: 20 document date: 2012-11-22



[pdf] Instructions

Medaval rating of Blood Pressure Monitors recommended by BPMs the Japanese Society Hypertension
5th March 2021 Page 1 3 medaval ie docs lists Warning Wrist monitor accuracy requires strict adherence
to manufacturer instructions Omron HBP 1300 E HEM 9700T 5 mar — |||

BPMs recommended by the Japanese Society of Hypertension - 5th March 2021
Page 1 of 3 Medaval rati ... d Validation Dretec BM-200 Nissei DSK-1011 Omron
HEM-7114 Dretec BM-201 Nissei DSK-1031 Omron **HEM-7200** V1 Dretec BM-202
Nissei DSK-1051 Omron HEM-7210 Fukuda Denshi FORM-5 Nissei WSK-1021...
lang:en score:24 filesize: 1.24 M page_count: 3 document date: 2021-03-05

[\[pdf\]](#) Catalog

IM HEM 7200 E2 V 01 08 2011 TR omron healthcare 80 data catalog 3 4450 1

retici AB temsilcisi retim Tesisi Yan Kurulu OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, ... bileenleri kontrol edin Trke Dijital Otomatik Kan Basinci Im Cihazı Model M3 Kullanım Kilavuzu IM-**HEM-7200**-E2 V -01-08/2011 1 F L c d 2 b J a indekiler OMRON M3 Dijital Otomatik Kan Basinci ...

lang:id **score:24** filesize: 1.75 M page_count: 28 document date: 2013-03-22

[\[pdf\]](#) User Manual

HEM 7200 AP3 B M Untitled Omron Healthcare Asia Pacific occurs consult your local OMRON representative Model JPN1 Instruction Manual 6 Technical Data Product Description Automatic Blood Pressure Monitor Title Created Date 8 31 2015 9 35 24 AM d236a97d40f34a55b321112e6e063f97 omronhealthcare ap Content products 80 resources

lang:id **score:24** filesize: 8.5 M page_count: 2 document date: 2016-06-10

[\[pdf\]](#) User Manual

HEM 7200 D qxd G4 120 m healthcare omron co jp support manual

lang:id **score:23** filesize: 2.1 M page_count: 20 document date: 2016-06-09

[\[pdf\]](#) Instructions

Medaval BPMs Validated for Pregnancy Page 1 of 2 html Blood Pressure Monitors Jan 24 2023 —

Warning Wrist monitor accuracy requires strict adherence to manufacturer instructions If the is not positioned correctly HealthSTATS BPro G2 pages bpm medaval ie docs lists |||

BPMs Validated for Pregnancy - 24th January 2023 Page 1 of 2 Blood Pressure Monitors Validated for ... ell Touch 3G Microlife WatchBP Office Target BP 3MD14 Ranor RN-032C Beurer BM 47 Omron BP742 **HEM-7200**-Z Riester RBP-100 Beurer BM 58 Omron BP760 HEM-7220-Z Riester ri-champion N Braun E...


lang:en **score:23** filesize: 252.24 K page_count: 2 document date: 2023-01-24

[\[pdf\]](#) User Manual

untitled Model M3 Manual de instrucțiuni baterii „AA dup cum este indicat apoi pune i la loc capacul acestui Dac data ora nu sunt setate în timpul m sur torii sau efectuarea utilizare tensiometru Omron RO s cdnmpo 129501539 content Manuale Tensiometre |||

Produtor Reprezentant pentru UE Unitate de productie Subsidiar OMRON HEALTHCARE Co., Ltd. 53, Kunot ... ificai urmtoarele componente Romn Tensiometru digital automat Model M3 Manual de instrucțiuni IM-**HEM-7200**-E2 V -01-08/2011 1 F L c d 2 b J a Cuprins V mulumim c ai achiziionat tensiometrul ... lang:it **score:23** filesize: 747.2 K page_count: 28 document date: 2014-01-22

[\[pdf\]](#)



Manuale uso Omron M3 HEM 7131 E MigliorPrezzo it MedisanShop OMRON è un misuratore di pressione compatto completamente automatico il cui funzionamento si basa sul principio oscillometrico L apparecchio misura in sfigmomanometro srsItd

AfmBOoWqVk3ocFs5TjgpkJS9ZDjWQqcXTDOuYr1ac4Vx1BWYXU8uoq2 medisanshop media setlab allegati prodotto |||

Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto 617-0002 JAPAN EU-r ... ensure Monitor Model M3 Instruction Manual English Français Deutsch Italiano Español Nederlands IM-**HEM-7200**-E2 V -01-08/2011 5333267-0A 1 F L c d 2 b J a Contents Thank you for purchasing the... lang:en score:22 filesize: 854.84 K page_count: 30 document date: 2017-09-04

[\[pdf\]](#)

Manuale uso Omron M3 HEM 7131 E MigliorPrezzo it MedisanShop4 1 Icone e messaggi di errore 98 4 2 Risoluzione dei Non vi è una definizione universalmente riconosciuta per l'ipotesioneManuale sfigmomanometromedisanshop media setlab allegati prodotto sfigmomanometro srsItd

AfmBOopuVZlItEpoOXRvg3T HrCmPySOjQe46Bh1dlfykEF XHENvDGg MedisanShopManuale sione4

Installing Replacing the Batteries Setting Buzzer Date and Time 5 3 Using Unit Manuale

sfigmomanometroManuale XHENvDGgmedisanshop AfmBOoqraOMu ZGlf7094XbX

HvRuQ9hTtOeYcm07tSR5gpk t38Y4jh Unit OMRON un misuratore pressione compatto completamente automatico il cui funzionamento si basa sul principio oscillometrico L apparecchio misura in Manuale t38Y4jhmedisanshop AfmBOopASFvbLz8x0SHTJF8KvRQhhtagYETQPZeISkA4txydD0pDw3AJ in La più elevata nel ciclo detta Pressione sistolica mentre quella bassa diastolica Entrambi i valori quello Sistolico e Manuale AfmBOopASFvbLz8x0SHTJF8KvRQhhtagYETQPZeISkA4txydD0pDw3AJmedisanshop

AfmBOoQdJ vySTj8FE1c1CoxGikExAoGULo2Ui9MzQKEvPYQWzLgVXD e La

AfmBOopH1y4DLabOu70IPviFkL4tWyN Ca6IMFAU3yqw1t6lqA62f9medisanshop

vySTj8FE1c1CoxGikExAoGULo2Ui9MzQKEvPYQWzLgVXDmedisanshop

AfmBOoqp71CRLmPul8n6gTmZ0YQED5GpsOCaOXUECagXtJDC0HQXnqw8

AfmBOoYrPHTBjJZZNr3BE4LV0Urqcp3Nus92 Ea6z4MK8jO71V6uroPmedisanshop e 4 Errore visualizzato Causa Soluzione Le pulsazioni rilevate sono irregolari o deboli Rimuovere bracciale Attendere 3 Manuale

AfmBOoqp71CRLmPul8n6gTmZ0YQED5GpsOCaOXUECagXtJDC0HQXnqw8medisanshop

AfmBOoq0W8bDRX0TQdelt9TnEB 1npi7Y2pZXEmkJwoPtld4ZgltlWc8 3 OMRON

1npi7Y2pZXEmkJwoPtld4ZgltlWc8medisanshop

AfmBOookozzLiMfo3U7yBCAJEe9RsU5HJg3h2cxT5RkBs dqK j5Wi8W2 in Indicatore sgonfiaggio

Indicatore della memoria battito Premere pulsante Impostare avviso acustico ora j5Wi8W2medisanshop

AfmBOorixMvjFLSj0URr5tuO6qcqevuRQB4GfmP7yQ zboHEuE7QxFF5 arteriosa può risultare differente se misurata al braccio destro sinistro pertanto anche misurazione possono differenti Manuale


sfigmomanometrov1medisanshop zboHEuE7QxFF5medisanshop AfmBOorBAVKNIvxaWauo4YX


QrQjIASrtvNOZu7gBOU f2ql zBJTkg differenti Se display viene simbolo batteria in esaurimento spegnere quindi sostituirle contemporaneamente tutte le batterieManuale sfigmomanometrov1Manuale

zBJTkgmedisanshop AfmBOoqCptZWOjgRNNF7cDFRjyMGmgC43TtwZxXZDBGf1aLzUpr8Q 4T |||

Manufacturer OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto 617-0002 JAPAN EU-r ... ensure Monitor Model M3 Instruction Manual English Français Deutsch Italiano Español Nederlands IM-**HEM-7200**-E2 V -01-08/2011 5333267-0A 1 F L c d 2 b J a Contents Thank you for purchasing the...

lang:en score:22 filesize: 854.84 K page_count: 30 document date: 2017-09-04





[\[pdf\]](#) Instructions

Blood Pressure Monitors Validated for Pregnancy Medaval Th ese monitors have been clinically validated in 5 mar 2021 — Warning Wrist monitor accuracy requires strict adherence to manufacturer instructions Microlife WatchBP Office Central Twin200 CBP BPMs 5th March Page 1 of 2 These pregnancy and medaval ie docs lists



[\[pdf\]](#) User Manual Catalog

Omron blood pressure monitors wrist BLOOD PRESSURE MONITORING OMRON o ers you true quality accuracy and reliability Whether are looking for an upper arm or type device we have the omronkatalog 2013 ajmedical se PRODUCT CATALOGUE HEALTHCARE Hem 637 Manual Pdf Wrist Blood Pressure Monitor with Advanced Positioning |||

OMRON PRODUCT CATALOGUE 2013 HEALTHCARE AJ Medical Uddvgen 32A 181 30 Liding Tel: 08-767 29 70 Fax: ... whether applying it on your own or with assistance. 6 HEM-7223-E HEM-7213-E V HEM-7202-E V **HEM-7200**-E2 V M3W, M3 Our fully automatic upper arm monitors are designed with both comfort and ac... lang:en score:20 filesize: 4.35 M page count: 64 document date: 2013-08-26



[\[pdf\]](#) Instructions

Medaval rating of Blood Pressure Monitors recommended by the BPMs BIHS UK and IRL 5th March 2021
Page 1 of 4 British Irish Hypertension medaval ie docs lists Warning Wrist monitor accuracy requires strict
adherence to manufacturer instructions HoMedics BPA 540 Omron M2 HEM 7117 E Society 5 mar —
SocietyBPMs 17 February 2025Kinetik TMB 1117 S H L HL168KB Kinetik 1491 HL168ZA 1970 HL858BA
2080 HL858DK 2083 N HL858GABPMs HypertensionBPMs Hypertension5 —BPMs IRLBPMs
IRLmedaval |||
BPMs recommended by the BIHS UK and IRL - 5th March 2021 Page 1 of 4 Medaval
rating of Blood Pre ... CT40 Meditech ABPM-06 Transtek TMB-1491 Omron HBP-
1300 HBP-1300-E Withings BP-800 Omron M3 **HEM-7200**-E Yuwell YE690A 2021
Medaval Ltd. 48 Serpentine Ave, Dublin 4, IRELAND. 353 1 664 376...
lang:en score:20 filesize: 5.32 M page count: 4 document date: 2021-03-05

Clinical validations Overview OMRON blood pressure monitors June 2013					
Upper Arm					
Name	Model	ESB	Validation Protocol	BHS	Special cases
M1 Classic	HEM-401-C-E	Passed ¹⁾			
M1 Plus	HEM-401-C-E	Passed ¹⁾			
M1 Compact	HEM-401-E	Passed ¹⁾			
M1	HEM-401-E	Passed ¹⁾			
HEM-ACULAR	HEM-400-ES-E	Passed ¹⁾			
M10 Plus	HEM-742-E	Passed ¹⁾			
M10 Basic	HEM-742-E	Passed ¹⁾			
M10 Compact	HEM-742-E	Passed ¹⁾			
M3 Professional	HEM-7101-E	Passed ¹⁾			
M2 Basic	HEM-7116-E/E8	Passed ¹⁾			
M2	HEM-7117-E/E8	Passed ¹⁾			
M2	HEM-7117-E/V	Passed ¹⁾			
M3	HEM-7200-E/E8	Passed ¹⁾			
M3	HEM-7202-E	Passed ¹⁾			
M3V	HEM-7202-E/V	Passed ¹⁾			
750T	HEM-759-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
750CP-1	HEM-759P-E2	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
M5	HEM-7001-E	Passed ¹⁾	Passed ¹⁾	Obese ¹⁾ Children ¹⁾ Elderly ¹⁾	
M5	HEM-7211-E/E8	Passed ¹⁾	Passed ¹⁾	Obese ¹⁾ Children ¹⁾ Elderly ¹⁾	
M5V	HEM-7211-E/V	Passed ¹⁾	Passed ¹⁾	Obese ¹⁾ Children ¹⁾ Elderly ¹⁾	
M41	HEM-750-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
M51	HEM-750-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
M7	HEM-780-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Pregnancy ¹⁾	


[pdf] Declaration of Conformity Catalog

Microsoft Word Clinical validation overview of OMRON devices June 03062013 doc meijan Validation list

2013 2 omron healthcare 80 data catalog 3 3393 1

Clinical validations Overview OMRON blood pressure monitors June 2013 Upper Arm Name M1 Classic M1 ... 102-E HEM-7101-E V HEM-7051-E HEM-7116-E/E8/ HEM-7116-E2 HEM-7117-E/E8 HEM-7117-E V HEM-7119-E V **HEM-7200-E/E8/ HEM-7200-E2 V HEM-7202-E V HEM-759-E HEM-759P-E2 HEM-7001-E HEM-7211-E/E8 HEM-7211-...**

lang:en **score:18** filesize: 44.75 K page_count: 5 document date: 2013-06-04



[pdf] Instructions Guide

110114 HEM 7200 BP742CAN IM FRCA r1 indd administrator Série 5 Tensiomètre FRN 012720111

01272011 omronhealthcare.ca 80

GUIDE DE L UTILISATEUR Srie 5 Tensiomt্রে Modle BP742CAN CANADIAN

FRENCH TTAABBLEEDOEFSCMOANTTIEN ... onez sans frais au :

www.omronhealthcare.ca 1-800-634-4350 26 SPCIFICATIONS Modle : BP742CAN

HEM-7200-CA Affichage : Affichage numrique cristaux liquides Plage de mesure :

Tension : 0 29...

lang:fr **score:18** filesize: 7.35 M page_count: 28 document date: 2011-01-27

Clinical validations Overview OMRON blood pressure monitors March 2013					
Upper Arm					
Name	Model	ESB	Validation Protocol	BHS	Special cases
M1 Classic	HEM-401-C-E	Passed ¹⁾			
M1 Plus	HEM-401-C-E	Passed ¹⁾			
M1 Compact	HEM-401-E	Passed ¹⁾			
M1	HEM-401-E	Passed ¹⁾			
HEM-ACULAR	HEM-400-ES-E	Passed ¹⁾			
M10 Plus	HEM-742-E	Passed ¹⁾			
M10 Basic	HEM-742-E2	Passed ¹⁾			
M10 Compact	HEM-742-E	Passed ¹⁾			
M3 Professional	HEM-7101-E	Passed ¹⁾			
M2 Basic	HEM-7116-E/E8	Passed ¹⁾			
M2	HEM-7117-E/E8	Passed ¹⁾			
M2	HEM-7117-E/V	Passed ¹⁾			
M3	HEM-7200-E/E8	Passed ¹⁾			
M3	HEM-7202-E	Passed ¹⁾			
M3V	HEM-7202-E/V	Passed ¹⁾			
750T	HEM-759-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
750CP-1	HEM-759P-E2	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
M5	HEM-7001-E	Passed ¹⁾	Passed ¹⁾	Obese ¹⁾ Children ¹⁾ Elderly ¹⁾	
M5	HEM-7211-E/E8	Passed ¹⁾	Passed ¹⁾	Obese ¹⁾ Children ¹⁾ Elderly ¹⁾	
M5V	HEM-7211-E/V	Passed ¹⁾	Passed ¹⁾	Obese ¹⁾ Children ¹⁾ Elderly ¹⁾	
M41	HEM-750-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
M51	HEM-750-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Elderly ¹⁾	
M7	HEM-780-E	Passed ¹⁾	Passed ¹⁾	Children ¹⁾ Obese ¹⁾ Pregnancy ¹⁾	

[pdf]

Clinical validation overview of OMRON devices haamar *10 Validation the Omron M5 I R5 and HEM 907

automated blood pressure monitors in elderly individuals according to International Protocol the Clinical

March 25032013s cdnmpo 129501539 content Tensiometre 25032013 |||

Clinical validations Overview OMRON blood pressure monitors March 2013 Upper

Arm Name M1 Classic M1 ... 102-E HEM-7101-E V HEM-7051-E HEM-7116-E/E8/

HEM-7116-E2 HEM-7117-E/E8 HEM-7117-E V HEM-7119-E V **HEM-7200-E/E8/**

HEM-7200-E2 V HEM-7202-E V HEM-759-E HEM-759P-E2 HEM-7001-E HEM-7211-

E/E8 HEM-7211-...

lang:en **score:18** filesize: 157.39 K page_count: 4 document date: 2013-05-26


K121932		Omron HEM-7200-Z (BP742)
Premarket Notification 510(k) Section 5-510(k) Summary	Device Name Page 1 of 7 05/2012	
Omron Healthcare, Inc. 1122 N. 1st Ave PMB 500 Lake Forest, IL 60045, USA		AUG 2 9 2012
Official Contact:	Tel: 847-247-9626 Fax: 847-247-9626	
Proprietary or Trade Name:	Mirna DiPano - Director, Quality & Regulatory	
Common/Usual Name:	Model HEM-7200-Z (BP742)	
Classification Name/Code:	UNN - Noninvasive blood pressure measurement 21078, 870, 110 Class II	
Device:	Model HEM-7200-Z (BP742)	
Predecessor Device:	Omron HEM-741C (G1, N - K12151)	
Device Description: The device is an automatic non-invasive blood pressure system. The device is battery powered and can also be powered from an AC adapter. The device inflates a cuff with an air pump, then deflates the cuff to an electronically controlled valve. During deflation the cuff pressure is monitored and pulse waveforms data is extracted. The extracted pulse waveforms data is then analyzed by software which determines pulse rate as well as systolic and diastolic pressure. The device is intended to be used with specified Omron cuffs in order sizes to encompass some ranging from 9 to 17 inches in circumference. The device also simulates the appearance of irregular heartbeats during measurement.		
Intended Use: Over-the-counter		
Intended Population: This device is intended for use on adults		
		Page 1 of 7 PDF Page 15 of 542

[pdf]

K121932 accessdata fda gov cdrh docs 12 |||

Premarket Notification 510 k Section 5-510 k Summary Omron Healthcare, Inc. 1925 West Field Court Lake Forest, IL 60045 USA K1 21 932 51 0 k Summary Page 1 ofs5 6/27/2012 Tel - 847-247-5626 Fax - 847-247-5626 Omron **HEM-7200-Z** BP742 AUG 2 9 2012 Official Contact: Mirna DiPano- Director, Qual...

lang:en score:15 filesize: 193.84 K page_count: 8 document date: 0000-00-00



[pdf]

Simedo Model M3 Gebruiksaanwijzing Omron handleiding nl bloeddrukmeterswebshop media ||| Digitale automatische bloeddrukmeter Model M3 Gebruiksaanwijzing Nederlands i IM-**HEM-7200-E**-01-10/09 5323430-0A Inhoud Dank u voor de aankoop van de OMRON M3 digitale automatische bloeddrukmeter. De OMRON M3 is een compacte, volledig automatische bloeddrukmeter die werkt volgens het oscillometr... lang:nl score:15 filesize: 1.35 M page_count: 29 document date: 2013-04-19

BPMs recommended by the Japanese Society of Hypertension 17 February 2025		Page 1 of 4
Medaval™ Certified Accuracy		
Medical rating of Blood Pressure Monitors recommended by the Japanese Society of Hypertension Monitors are listed according to clinical validation.		
Medical certification is awarded only to devices with the most accurate measurement technology. For accreditation, proof of reliability and adherence to international standards is continuously validated with independent peer review. Devices without evidence of validation include those where the cited paper refers to a different device and those without evidence of possible health impairment via a validated device. Medical self-regulation device where security proof is validated or experience is provided. Warning: Your monitor accuracy requires strict adherence to manufacturer instructions. If the monitor is not positioned correctly, the results can be very inaccurate. Symbols: () Listed as Unavailable () Listed Potentially () Discontinued		
***** Certified 5-Star Blood Pressure Monitors *****		
ARD UA 401 ()	***** Certified 4-Star Blood Pressure Monitors *****	
Omron HEM-7120 (BP) ()	Omron HEM-7120F () ()	
Omron HEM-7200 () ()	Omron HEM-7200F () ()	
Omron HEM-7200F () ()		
*** 3-Star Blood Pressure Monitors ***		
Omron Evolv (HEM-7600T) ()		
** 2-Star Blood Pressure Monitors **		
ARD TM 2441	Omron HEM-1200 (HEM-1200-F) () ()	
ARD UA 1100SC	Omron HEM-901C	
ARD UA 1200SL	Omron HEM-901C (HEM-1200)	
ARD UA 401 Plus (UA 400A) () ()	Omron HEM-900 ()	
ARD UA 401SL Plus ()	Omron HEM-7200	
ARD UA 21 ()	Omron HEM-7600T ()	
ARD UA 1100SL ()	Omron HEM-7600T ()	
Nissei DM-3000 ()	Omron HEM-4700T ()	
* 1-Star Blood Pressure Monitors *		
ARD TM 2430 () ()	Panola Denshi DM-3000 ()	
ARD UA 1024 ()	Omron HEM-901 (HEM-901-F) ()	
ARD UA 101 () ()		
Other Validated Blood Pressure Monitors		
ARD TM 3015 () ()	ARD UA 401 () ()	
ARD TM 3017 () ()	ARD UA 401C	
ARD TM 2000 () ()	ARD UA 401C	
ARD UA 401 () ()	ARD UA 401C	
ARD UA 700 ()	ARD UA 401C	
ARD UA 701 ()	Omron C1-401C ()	
ARD UA 701 Plus ()	Omron C1-401C ()	
ARD UA 774 () ()	Panola Denshi DM-3000	
ARD UA 774C	Panola Denshi DM-3000	
ARD UA 707 () ()	Omron HEM-900 ()	
ARD UA 401 ()	Omron HEM-901 (HEM-901-F) ()	
ARD UA 401C	Toshiba Elextra ES-115	
© 2025 Medaval Ltd. All Rights Reserved. Data as of September 10, 2025. 10/1/2025 FHM info@medaval.co.uk info@medaval.co.uk		

[pdf] Instructions

BPMs recommended by the Japanese Society of Hypertension 17 February 2025 Feb · Omron HEM 7114 Dretec BM 202 7200 V1 Fukuda Denshi FBX Solar 4500 Sole DSL 8001 HEM medaval ie docs lists ||| BPMs recommended by the Japanese Society of Hypertension - 17 February 2025 Page 1 of 4 Medaval ra ... 021 Omron HEM-6200 Omron HEM-7114 Dretec BM-202 Fukuda Denshi FBX-1000 Fukuda Denshi FORM-5 Omron **HEM-7200** V1 Omron HEM-7210 Panasonic EW-BU17 Fukuda Denshi SC-1800 COL-SC-1800 Nissei DS-R10J N...

lang:en score:15 filesize: 491.89 K page_count: 4 document date: 2025-05-08

