

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [Coloplast](#) /

> [Biatain Adhesive Foam Dressing User Manual](#)

## Coloplast B00833MY5K

### Biatain Adhesive Foam Dressing User Manual

Model: B00833MY5K

#### PRODUCT OVERVIEW

The Biatain Adhesive Foam Dressing is a super absorbent foam dressing designed for wound care. It features a hydrocolloid-based adhesive border, making it suitable for a range of exuding wounds.

This dressing provides padding and protection and can be used throughout the healing process. Its unique 3D foam structure conforms closely to the wound bed, even under compression.



Image: Biatain Adhesive Foam Dressing. This image shows the top view of the square-shaped foam dressing with its adhesive border.

#### KEY FEATURES

- Soft and conformable polyurethane foam dressing with hydrocolloid adhesive.
- Suitable for a wide range of exuding wounds, including leg ulcers, pressure ulcers, second-degree burns, donor sites, postoperative wounds, and skin abrasions.
- Unique 3D foam structure conforms to the wound bed, even under compression.
- Well-suited for use on wounds with healthy peri-ulcer skin and where extra adhesion is needed.
- Clinical studies indicate a wear time of up to 7 days.

#### SETUP AND PREPARATION

Before applying the Biatain Adhesive Foam Dressing, ensure the wound area is properly prepared.

1. Clean the wound thoroughly according to standard wound care protocols.
2. Ensure the skin around the wound is clean and dry to facilitate proper adhesion.
3. Select the appropriate size of Biatain dressing. If necessary, the dressing can be cut to size to fit the wound area, ensuring the foam covers the wound and the adhesive border adheres to healthy skin.

## OPERATING INSTRUCTIONS (APPLICATION)

Follow these steps for effective application of the dressing:

1. Remove the dressing from its sterile packaging.
2. Peel off the protective film from the adhesive side of the dressing.
3. Gently apply the foam side of the dressing directly onto the wound.
4. Press the adhesive border firmly onto the surrounding healthy skin to ensure secure attachment.
5. The dressing can be used under compression devices if indicated for the wound type.
6. The dressing can remain in place for up to 7 days, depending on the wound exudate level and clinical assessment.

## MAINTENANCE AND REMOVAL

Proper maintenance and removal are essential for wound healing and patient comfort.

- **Monitoring:** Regularly inspect the dressing and the surrounding skin for signs of leakage, saturation, or adverse reactions.
- **Change Frequency:** Change the dressing when clinically indicated, such as when exudate reaches the edge of the foam, or after a maximum of 7 days.
- **Removal:** To remove, gently lift one corner of the adhesive border and slowly peel the dressing away from the skin. Dispose of used dressings according to local medical waste regulations.

## TROUBLESHOOTING

Common issues and their resolutions:

Problem	Possible Cause	Solution
Dressing not adhering properly	Skin not clean or dry; excessive hair; incorrect application.	Ensure skin is clean, dry, and free of excessive hair. Press firmly on the adhesive border during application.
Leakage around the dressing	Dressing saturated with exudate; incorrect size for wound.	Change dressing more frequently. Consider a larger size or a dressing with higher absorbency if exudate levels are very high.
Skin irritation around the dressing	Allergic reaction; prolonged wear; improper removal.	Discontinue use if irritation persists. Ensure gentle removal. Consult a healthcare professional.

## SPECIFICATIONS

Attribute	Detail
Brand	Coloplast

Attribute	Detail
Model Number	B00833MY5K
Dimensions	4 x 4 inches (dressing size)
Material	Polyurethane (PU) foam with hydrocolloid adhesive
Sterility Rating	Sterile
Material Feature	Absorbent, Adherent
Wear Time	Up to 7 days

## IMPORTANT SAFETY INFORMATION

This product is intended for wound care as described. Consult a healthcare professional for specific medical advice or if conditions worsen.

*Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This disclaimer is provided for informational purposes and does not apply to the medical device itself but to general health claims.*

## WARRANTY AND SUPPORT

For product inquiries, technical support, or warranty information, please contact Coloplast directly. Refer to the product packaging for specific contact details or visit the official Coloplast website.

Manufacturer: Coloplast

For further assistance, please visit the [Coloplast official website](#).

© 2024 Coloplast. All rights reserved.

This manual is for informational purposes only and should not replace professional medical advice.