

Madoec LC-200

Madoec Mega Motion LC-200 Easy Comfort 3-Position Electric Lift Recliner User Manual

Model: LC-200

1. INTRODUCTION

This manual provides essential instructions for the safe and effective use of your Madoec Mega Motion LC-200 Easy Comfort 3-Position Electric Lift Recliner. Please read this manual thoroughly before assembly and operation to ensure proper function and to prevent injury or damage. Keep this manual for future reference.

The LC-200 is designed to provide comfort and assistance with sitting, reclining, and standing. It features a 3-position electric lift and recline mechanism, allowing for lifted, neutral, and reclined positions. The chair is constructed with furniture-grade hardwood and plywood, and includes an integrated emergency battery backup system.



Figure 1: Madoec Mega Motion LC-200 Lift Recliner

This image displays the Madoec Mega Motion LC-200 Lift Recliner in a slightly elevated position, showcasing its chocolate brown fabric upholstery and the distinctive tufted backrest. The chair's robust base mechanism is partially visible, indicating its lift functionality.

2. SAFETY INFORMATION

WARNING: Failure to follow these safety instructions may result in injury or damage to the product.

- Do not exceed the maximum weight capacity of **325 pounds**.
- Assemble the chair according to the instructions provided to ensure stability and safety.
- Be cautious when using the electric lift and recline functions. Ensure the area around the chair is

clear of obstructions, children, and pets before operating.

- Keep the power cord away from heat sources, sharp objects, and areas where it could be tripped over or damaged.
- In case of a power outage, use the integrated emergency battery backup system. Refer to the "Operating Instructions" section for details.
- Always turn off the chair by unplugging it from the wall outlet when not in use or before cleaning and maintenance.
- This chair is intended for indoor use only.
- Do not allow children to play on or around the chair.
- Avoid placing liquids or heavy objects on the chair's arms or backrest.

3. PACKAGE CONTENTS

Before beginning assembly, ensure all components are present and undamaged. If any parts are missing or damaged, do not proceed with assembly and contact customer support.

- Chair Base Assembly
- Backrest
- Power Cord
- Hand Control Remote
- User Manual (this document)

4. SETUP AND ASSEMBLY

Assembly typically requires two adults. Follow these steps carefully:

1. **Unpack Components:** Carefully remove all components from the packaging. Place them on a soft, clean surface to prevent scratching.
2. **Attach Backrest:** Align the metal brackets on the backrest with the corresponding posts on the chair base. Slide the backrest down firmly until it locks into place. You should hear a click.
3. **Connect Power:** Locate the power cord and connect it to the power transformer, then connect the transformer to the chair's power input. Plug the power cord into a grounded electrical outlet.
4. **Install Battery Backup (Optional/Emergency):** If your model includes an emergency battery backup, locate the battery compartment (usually on the transformer or under the chair). Insert two 9V batteries (not included) as indicated. This system allows limited use during a power outage.
5. **Test Operation:** Before sitting, use the hand control remote to test the lift and recline functions to ensure proper assembly and connection.



Figure 2: Chair Dimensions

This image illustrates the key dimensions of the Madoec Mega Motion LC-200 Lift Recliner: 41.5 inches in height, 34 inches in width, and 38.5 inches in depth. These measurements are crucial for planning placement and ensuring the chair fits comfortably in its intended space.

5. OPERATING INSTRUCTIONS

The Madoec Mega Motion LC-200 is operated using the wired hand control remote. Ensure the chair is plugged into a live electrical outlet before use.

5.1 Hand Control Remote

The remote typically features two buttons:

- **UP Button:** Press and hold to raise the chair to a standing assist position or to bring it to an upright seated position from a reclined state.
- **DOWN Button:** Press and hold to recline the chair or to lower it from a standing assist position to a seated position.

5.2 Chair Positions

The LC-200 offers three primary positions:

1. **Lifted Position:** Press and hold the UP button until the chair reaches a comfortable height for easy entry or exit.

2. **Neutral/Seated Position:** Press the DOWN button to lower the chair from the lifted position, or the UP button to raise it from a reclined position, until it is fully upright.
3. **Reclined Position:** Press and hold the DOWN button to smoothly recline the backrest and elevate the footrest. Release the button when the desired recline angle is achieved.



Figure 3: Lifted Position

This image shows the Madoec Mega Motion LC-200 Lift Recliner in its fully lifted position, designed to assist users in standing up or sitting down with ease. The chair is tilted forward and elevated.



Figure 4: Reclined Position

This image depicts the Madoec Mega Motion LC-200 Lift Recliner in a comfortable reclined position, with the footrest fully extended. This position is ideal for relaxation, reading, or napping.

5.3 Emergency Battery Backup

The integrated emergency battery backup system allows for limited operation of the chair during a power failure. This system requires two 9V batteries (not included). When power is lost, the chair can be moved to an upright or reclined position using the remote. The battery backup is for emergency use only and is not intended for continuous operation.

6. MAINTENANCE

6.1 Cleaning

- **Fabric:** For general cleaning, use a soft, damp cloth with mild soap and water. Test any cleaning solution on an inconspicuous area first. Do not saturate the fabric.
- **Frame:** Wipe the metal frame components with a dry cloth. Avoid abrasive cleaners.
- **Power Cord:** Regularly inspect the power cord for any signs of damage. If damaged, discontinue use and contact customer support for a replacement.

6.2 General Care

- Ensure the chair is placed on a level surface.
- Avoid exposing the chair to direct sunlight or extreme temperatures, which can damage the fabric

and internal components.

- Periodically check all connections and fasteners to ensure they are secure.

7. TROUBLESHOOTING

If you encounter issues with your LC-200 lift recliner, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Chair does not move.	<ul style="list-style-type: none">• Power cord unplugged.• Outlet not functioning.• Loose connection at transformer or chair.• Hand control remote malfunction.	<ul style="list-style-type: none">• Ensure power cord is securely plugged into a live outlet.• Test the outlet with another device.• Check all cable connections.• Contact customer support if remote is suspected faulty.
Chair moves slowly or intermittently.	<ul style="list-style-type: none">• Low battery in emergency backup (if in use).• Overload (exceeding weight capacity).• Obstruction in mechanism.	<ul style="list-style-type: none">• Replace 9V batteries in the emergency backup system.• Ensure user weight is within 325 lbs capacity.• Inspect for any objects obstructing the mechanism.
Noise during operation.	<ul style="list-style-type: none">• Normal operational sounds.• Loose components.	<ul style="list-style-type: none">• Some motor noise is normal.• Check for loose bolts or connections.

If the problem persists after attempting these solutions, please contact Madoec customer support for further assistance.

8. SPECIFICATIONS

Feature	Detail
Model Name	LC-200
Brand	Madoec
Color	Brown (Chocolate Fabric)
Product Dimensions (D x W x H)	38.5"D x 34"W x 41.5"H
Item Weight	130 pounds
Maximum Weight Capacity	325 pounds
Reclining Positions	3 (Lifted, Neutral, Reclined)
Mechanism	Electric Power Lift and Recline
Frame Material	Furniture grade hardwood / plywood construction, Heavy-duty steel lift mechanism

Feature	Detail
Special Features	Integrated Emergency Battery Backup
UPC	619084064800

9. WARRANTY AND SUPPORT

The Madoec Mega Motion LC-200 Lift Recliner comes with a heavy-duty steel lift mechanism that has a lifetime warranty. For specific warranty details on other components and for customer support, please refer to the warranty card included with your purchase or contact Madoec customer service directly.

For assistance with assembly, operation, troubleshooting, or warranty claims, please contact the manufacturer's customer support. Contact information can typically be found on the product packaging or the manufacturer's official website.