

[Manuals.plus](#) /

› [Fisher & Paykel](#) /

› Fisher & Paykel FlexiFit 432 Cushion/Seal Instruction Manual

Fisher & Paykel 4332435628

Fisher & Paykel FlexiFit 432 Cushion/Seal Instruction Manual

Model: 4332435628

1. PRODUCT OVERVIEW

This manual provides essential instructions for the Fisher & Paykel FlexiFit 432 Full Face Mask Seal and Foam Cushion Kit. This kit is designed as a replacement component for the FlexiFit 432 Full Face Mask, ensuring proper fit and function for respiratory therapy.



Image 1.1: The Fisher & Paykel FlexiFit 432 Cushion/Seal Kit. This image displays the replacement cushion and seal components, typically made of silicone and foam, designed to fit the FlexiFit 432 Full Face Mask.

Key Features:

- Replacement Full Face Mask Seal and Foam Cushion Kit.
- Available in multiple sizes: Small, Medium, Large, and Extra Large to ensure a proper fit.
- Designed for use with the Fisher & Paykel FlexiFit 432 Full Face Mask.

2. SETUP AND INSTALLATION

Proper installation of the cushion and seal is crucial for effective therapy and comfort. Follow these steps to replace your FlexiFit 432 Full Face Mask cushion and seal:

1. **Remove the Old Cushion:** Carefully detach the worn or damaged cushion and seal from the mask frame. Refer to your original mask manual for specific removal instructions if needed.
2. **Inspect the Mask Frame:** Ensure the mask frame is clean and free from debris or damage before installing the new cushion.
3. **Align the New Cushion:** Position the new FlexiFit 432 cushion and seal onto the mask frame. Ensure that the alignment tabs or grooves on the cushion match those on the mask frame.
4. **Secure the Cushion:** Gently press the cushion into place around the entire perimeter of the mask frame until it is securely seated. Verify there are no gaps or loose areas.
5. **Perform a Seal Check:** Once installed, connect the mask to your CPAP device and perform a brief seal check as described in your CPAP device's user manual. Adjust the headgear if necessary to achieve a comfortable and leak-free seal.

3. OPERATING INSTRUCTIONS

The FlexiFit 432 Cushion/Seal is an integral part of your full face mask system. Once installed, its operation is directly tied to the overall function of the mask and your CPAP therapy.

- **Proper Fit:** Ensure the mask, with the new cushion, fits snugly but comfortably on your face. The cushion should create a seal without excessive pressure or discomfort.
- **Headgear Adjustment:** Adjust the headgear straps to achieve an optimal seal. Avoid overtightening, which can lead to skin irritation and discomfort.
- **Leak Management:** If leaks occur, gently reposition the mask or adjust the headgear. Persistent leaks may indicate an incorrect size or a need for further adjustment.
- **Breathing Comfort:** The cushion is designed to provide a soft interface for comfortable breathing during therapy. If you experience discomfort, re-evaluate the fit and consider if a different cushion size might be more appropriate.

4. MAINTENANCE AND CLEANING

Regular cleaning and maintenance of your FlexiFit 432 Cushion/Seal are essential for hygiene, comfort, and extending the life of the product. Daily cleaning is recommended.

1. **Daily Cleaning:** Gently wash the cushion and seal daily with warm water and a mild, non-abrasive soap. Avoid soaps with moisturizers, perfumes, or antibacterial agents, as these can damage the material.
2. **Rinsing:** Thoroughly rinse all soap residue from the cushion and seal under running water.
3. **Drying:** Allow the cushion and seal to air dry completely out of direct sunlight. Do not use a clothes dryer or heat to dry.
4. **Weekly Cleaning:** For a more thorough cleaning, you may soak the cushion and seal in a solution of warm water and mild soap for approximately 10-15 minutes before rinsing and air drying.
5. **Inspection:** Regularly inspect the cushion and seal for signs of wear, tear, discoloration, or hardening.

Replace the cushion if any damage is observed or if the seal becomes compromised.

5. TROUBLESHOOTING

If you encounter issues with your FlexiFit 432 Cushion/Seal, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Air Leaks Around Cushion	Improper fit, worn cushion, incorrect size, headgear too loose/tight.	Re-adjust headgear, reposition mask, check cushion for damage, ensure correct cushion size, replace cushion if worn.
Skin Irritation/Red Marks	Mask overtightened, cushion not clean, incorrect cushion size.	Loosen headgear, clean cushion thoroughly, ensure correct cushion size, consider cushion liners if irritation persists.
Cushion Feels Stiff/Hard	Accumulation of oils/debris, material degradation over time.	Clean cushion more frequently and thoroughly, replace cushion as it may be at the end of its lifespan.

6. SPECIFICATIONS

Detailed specifications for the Fisher & Paykel FlexiFit 432 Cushion/Seal Kit:

Item Model Number: 4332435628

Manufacturer: Fisher & Paykel

Package Dimensions: Approximately 4.4 x 4 x 1.5 inches

Item Weight: Approximately 2.24 ounces

ASIN: B00BCOTRSE

Material: Silicone seal with foam cushion (specific composition may vary)

Available Sizes: Small, Medium, Large, Extra Large

7. WARRANTY AND SUPPORT

For information regarding product warranty, technical support, or to purchase replacement parts, please contact Fisher & Paykel Healthcare directly or consult your authorized medical equipment provider. Keep your purchase receipt for warranty claims.

Contact Information: Refer to the official Fisher & Paykel Healthcare website or product packaging for the most current contact details.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question