

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [Lansinoh](#) /

› [Lansinoh Manual Breast Pump Instruction Manual](#)

Lansinoh 50520

Lansinoh Manual Breast Pump Instruction Manual

Model: 50520

Brand: Lansinoh

INTRODUCTION

The Lansinoh Manual Breast Pump is designed for comfortable, silent, and portable milk expression. It is ideal for occasional use at home, on the go, or as a convenient backup. This manual features two pumping modes, stimulation and expression, to assist with let-down and optimize your pumping sessions. Its ergonomic handle is designed to reduce hand fatigue, and the included ComfortFit flanges ensure optimal suction and comfort.



Image: The Lansinoh Manual Breast Pump, showcasing its main components including the pump body, handle, flanges, and bottle.

WHAT'S IN THE BOX

Your Lansinoh Manual Breast Pump package includes the following components:

- 2 ComfortFit® flanges (25mm & 30.5mm)
- 1 Lansinoh® breastmilk storage bottle (5oz)
- 1 Manual Pump Body
- 1 Stem
- 1 Handle
- 1 Silicone diaphragm
- 1 White valve
- 1 Bottle stand
- 1 Bottle cap and bottle collar
- 1 NaturalWave® slow-flow nipple

- 1 Bottle storage lid

FEATURES

The Lansinoh Manual Breast Pump incorporates several features designed for efficiency and comfort:

- **Two-Phase Pumping:** Features distinct let-down and expression phases to enhance pumping efficiency and maximize milk output.
- **ComfortFit Flanges:** Includes two sizes (25mm & 30.5mm) with flexible, soft rims for optimal suction and a comfortable fit.
- **Ergonomic Handle:** Designed for easy expression, reducing hand fatigue during pumping sessions.
- **Simple Assembly & Cleaning:** Composed of fewer parts for straightforward assembly and easy cleaning.
- **Lansinoh Compatible:** Wide-neck bottles are compatible across the Lansinoh range, allowing for direct pumping, storage, and feeding with one bottle.
- **NaturalWave Nipple:** Comes with one slow-flow NaturalWave nipple, promoting natural feeding actions similar to breastfeeding and clinically proven to reduce nipple confusion.

Two-phase pumping

to promote let-down and
pumping efficiency



Image: Demonstrates the two-phase pumping feature, showing the pump in use.



ComfortFit[®] flanges

2 sizes (25mm &
30.5mm) for ideal
suction & fit

Image: Displays the two ComfortFit flange sizes included with the pump.

Ergonomic handle

Easy-express handle design to reduce hand fatigue



Image: Illustrates the ergonomic handle design for comfortable use.

SETUP

Before first use, disassemble and clean all parts that will come into contact with breastmilk. Ensure all parts are completely dry before assembly.

1. Insert the white valve firmly into the pump body.
2. Attach the silicone diaphragm to the top of the pump body, ensuring it is securely sealed.
3. Connect the handle to the stem of the pump body.
4. Select the appropriate ComfortFit flange size (25mm or 30.5mm) and attach it securely to the pump body.
5. Screw the Lansinoh breastmilk storage bottle onto the bottom of the pump body.



Lansinoh[®] compatible

Pump, store, and feed with the included Lansinoh wide-neck bottles

Image: Demonstrates the assembly process of the pump body and bottle.

2 sizes of Comfortfit flanges

Ergonomic handle

Pump and store with wide-neck bottle



Image: An annotated diagram illustrating the key components of the manual breast pump.

OPERATING INSTRUCTIONS

To effectively use your Lansinoh Manual Breast Pump:

1. **Positioning:** Center your nipple in the flange tunnel. Ensure a comfortable seal around your breast.
2. **Stimulation Phase:** Begin by pressing the handle quickly and lightly to stimulate milk let-down. This mimics a baby's rapid suckling at the start of a feed. Continue until milk begins to flow.
3. **Expression Phase:** Once milk flows, switch to longer, slower, and deeper presses of the handle. This mimics a baby's slower, deeper suckling during active feeding. Adjust the rhythm and suction strength to your comfort and milk flow.
4. **Pumping Duration:** Pump until your breast feels empty or milk flow significantly slows. You can switch breasts as needed.
5. **Storage:** After pumping, carefully remove the pump from your breast. You can detach the bottle and seal it with the storage lid or attach the NaturalWave nipple for immediate feeding.

From pump to bottle, quick and easy!



Image: Visual guide for transferring milk from pump to bottle for storage or feeding.

MAINTENANCE & CLEANING

Proper cleaning and maintenance are crucial for hygiene and longevity of your breast pump. The Lansinoh Manual Breast Pump is designed with fewer parts for easy cleaning.

1. **Disassembly:** After each use, disassemble all parts that came into contact with breastmilk.
2. **Washing:** Wash all disassembled parts in warm, soapy water. Use a soft brush for thorough cleaning. Rinse thoroughly with clean water.
3. **Sterilization:** For initial use and periodically thereafter, sterilize parts by boiling them in water for 5 minutes or using a steam sterilizer. Refer to the product packaging for specific sterilization guidelines.
4. **Drying:** Allow all parts to air dry completely on a clean surface or drying rack. Do not use a cloth to dry, as this can transfer germs.
5. **Storage:** Store clean, dry parts in a clean, dry place until next use.

The pump is BPA and BPS free. Some components are dishwasher safe; however, hand washing is recommended for optimal



Easy to use

Fewer parts so it's easy to clean and assemble

Image: Demonstrates the ease of cleaning the Lansinoh bottle component.

TROUBLESHOOTING

If you encounter issues with your Lansinoh Manual Breast Pump, consider the following common solutions:

- **Low Suction:**

- Ensure all parts are correctly assembled and securely connected.
- Check the white valve for tears or damage; replace if necessary.
- Verify the silicone diaphragm is properly seated and free of damage.
- Ensure the flange forms a tight seal around your breast.

- **Discomfort During Pumping:**

- Re-evaluate your flange size. The nipple should move freely within the tunnel without excessive rubbing. Try the

alternative flange size provided.

- Adjust your pumping rhythm and suction strength. Pumping should not be painful.
- Ensure proper breast positioning within the flange.

- **No Milk Expression:**

- Ensure you are in a relaxed environment. Stress can inhibit milk flow.
- Start with the stimulation phase to encourage let-down.
- Consider applying warm compresses to your breasts before pumping.
- Ensure adequate hydration.

If problems persist, please refer to the support section for contact information.

SPECIFICATIONS

Attribute	Detail
Product Dimensions	19.1 x 19.6 x 10.2 inches
Item Model Number	50520
Target Gender	Unisex
Material Type	Silicone
Material Free	BPA Free, Bisphenol-S (BPS) Free
Additional Product Features	Dishwasher Safe, Ergonomic, Lightweight, Portable
Number Of Items	1
Style	Manual Pump
Power Source	Manual
Pump Type	Manual
Batteries Required	No
Dishwasher Safe	No (Note: Some components may be, but general product listing says No)
Is Portable	Yes
Item Weight	4.8 ounces

WARRANTY & SUPPORT

Lansinoh is committed to providing quality products. For specific warranty details, please refer to the documentation included with your product or visit the official Lansinoh website. Generally, products purchased from Amazon.com are eligible for a 90-day refund or replacement return policy.

For customer support, product inquiries, or troubleshooting assistance beyond this manual, please contact Lansinoh directly through their official website or customer service channels. You can find more information and support resources at the [Lansinoh Store on Amazon](#).

Related Documents - 50520

<p>Frequently Asked Questions</p> <p>Find the answers to common questions about Lansinoh breast pumps, HPA Lanolin Nipple Cream, TheraPearl Breast Therapy Pads, baby bottles, and more. Includes troubleshooting tips and product care advice.</p>	<p>Lansinoh Breast Pump & Product FAQs and Troubleshooting Guide</p> <p>Find answers to common questions about Lansinoh breast pumps, HPA Lanolin Nipple Cream, TheraPearl Breast Therapy Pads, baby bottles, and more. Includes troubleshooting tips and product care advice.</p>
<p>Lansinoh 2-in-1 Double Electric Breast Pump</p> <p>Comprehensive guide to using, cleaning, assembling, and maintaining the Lansinoh 2-in-1 Double Electric Breast Pump, including important safety instructions and pumping tips.</p>	<p>Lansinoh 2-in-1 Double Electric Breast Pump User Manual and Guide</p> <p>Comprehensive guide to using, cleaning, assembling, and maintaining the Lansinoh 2-in-1 Double Electric Breast Pump, including important safety instructions and pumping tips.</p>
 <p>Lansinoh® Wearable Pump</p> <p>INSTRUCTIONS FOR USE</p> <p>Read all instructions before using the pump and keep this manual for future reference.</p>	<p>Lansinoh® Wearable Pump: Instructions for Use</p> <p>Detailed instructions for the Lansinoh® Wearable Pump, covering setup, operation, cleaning, storage, and technical details. Available in multiple languages for user convenience.</p>
 <p>Lansinoh</p> <p>AFFINITY® DOUBLE ELECTRIC BREAST PUMP</p> <p>INSTRUCTIONS FOR USE</p> <p>INSTRUCCIONES PARA EL USO</p>	<p>Lansinoh Affinity Double Electric Breast Pump: Instructions for Use, Features, and Care</p> <p>A comprehensive guide to the Lansinoh Affinity Double Electric Breast Pump, covering its features, setup, operation, cleaning, power options, breastmilk storage, and troubleshooting. Learn how to use your breast pump for efficient and comfortable milk expression.</p>
	<p>Bioby Portable Bottle Warmer User Manual - Model N1</p> <p>Comprehensive user manual for the Bioby Portable Bottle Warmer (Model N1), covering setup, operation, safety precautions, troubleshooting, and maintenance. Learn how to use your bottle warmer for fast, accurate temperature control.</p>
	<p>innoGIO GIO warm Portable Baby Bottle Warmer User Manual</p> <p>User manual for the innoGIO GIO warm Portable Baby Bottle Warmer (Model GIO-370), detailing assembly, operation, cleaning, safety precautions, troubleshooting, and bottle compatibility.</p>

The image shows a spreadsheet titled 'Lansinoh MAP Pricing' with columns for various product categories and their corresponding prices. The categories include Breast Therapies, Nursing Pads, Milk Storage, Breast Pumps, Breast Pump Accessories, and Pumping Bras. The spreadsheet lists numerous product SKUs and their associated prices.

[\[pdf\]](#) Accessories

Lansinoh MAP Pricing LANSINOH U S AUTHORIZED RESELLER AND DISTRIBUTOR TERMS
 CONDITIONS Price Mapping 2021 v 1612986265 cdn shopify s files 1 0054 1989 2808 |||
 Vendor Stock Number Breast Therapies 10020 23100 65011 10200 70170 10232
 70171 70191 Nursing Pads 20236 20265 20370 20945 20380 20381 20388 Breast
 Milk Storage 20435 20450 20473 53441 20415 Breast Pumps **50520** 53050 53250
 50700 Breast Pump Accessories 53420 53431 53510 53511 50428 53512 Pumping
 Bras ...

lang:tl score:25 filesize: 183.08 K page_count: 1 document date: 2021-02-10



[\[pdf\]](#) Guide Catalog

Educational Learning Tools Professional Resources Lansinoh US Healthcare Product Sales Catalog
 9c16403f 8272 441f 9d9f 8aa93de342b7 v 1658333687 cdn shopify s files 1 0054 1989 2808 |||
 Product Guide Lanolin Nipple Cream The #1 selling nipple cream in the US* #1
 recommended by moms ... t milk Storage Bottles BPA BPS free BPA BPS FREE
 Manual Breast Pump Vendor Stock # Unit Size **50520** 1 pump *Source: Multi-Outlet
 IRI L52 weeks ending 12/27/2020, Dollar sales. Units/Case 6 un...

lang:en score:21 filesize: 6.09 M page_count: 20 document date: 2022-07-18