

[Manuals.plus](#) /

> [Mobiclinic](#) /

> Mobiclinic Delta Foldable Shower Chair Instruction Manual

Mobiclinic QA-00357/01

Mobiclinic Delta Foldable Shower Chair

Model: QA-00357/01

Brand: Mobiclinic

INTRODUCTION

This instruction manual provides important information for the safe and effective use of your Mobiclinic Delta Foldable Shower Chair. Please read this manual thoroughly before assembly and use, and keep it for future reference.



Image: The Mobiclinic Delta Foldable Shower Chair, highlighting its adjustable height, product weight (2.55 kg), supported weight (100 kg), and materials (steel and PVC). Dimensions are also shown.

SAFETY INFORMATION

Your safety is paramount. Adhere to the following guidelines:

- **Maximum Weight Capacity:** Do not exceed the maximum supported weight of **100 kg (220 lbs)**.
- **Stable Surface:** Always place the shower chair on a flat, stable, and non-slip surface.
- **Leg Adjustment:** Ensure all four legs are adjusted to the same height and are securely locked in place before each use. Uneven legs can cause instability.
- **Non-Slip Feet:** Verify that the rubber non-slip caps on the legs are intact and clean to prevent slipping in wet conditions.
- **Regular Inspection:** Periodically inspect the chair for any signs of wear, damage, or loose components. Discontinue use if any damage is found.
- **Proper Use:** Use the chair only for its intended purpose as a shower or bath seat. Do not stand on the chair.

SETUP

Follow these steps to set up your Mobiclinic Delta Foldable Shower Chair:

1. **Unfold the Chair:** Carefully unfold the shower chair until it is fully open and the frame locks into

position.

2. **Adjust Leg Height:** Press the push-pins on each leg and slide the inner leg tube to the desired height. Ensure the push-pin snaps securely into one of the adjustment holes. Repeat for all four legs, making sure they are all set to the same height. The adjustable height range is between 34 cm and 47 cm.
3. **Verify Stability:** Before sitting, gently press down on the seat to ensure the chair is stable and does not wobble. All legs should be firmly on the ground.



Image: Visual guide for quick assembly: 1. Unfold the chair. 2. Place the legs. 3. Ready to use!

OPERATING INSTRUCTIONS

The Mobiclinic Delta Shower Chair is designed for comfort and safety during showering or bathing:

- **Seating:** Sit gently on the chair, ensuring your weight is evenly distributed. The ergonomic, concave seat provides comfort and support.
- **Drainage:** The seat features drainage holes to prevent water accumulation, promoting hygiene and quick drying.
- **Portability:** When not in use, the chair can be easily folded for compact storage and transport, making it ideal for small spaces or travel.



Image: Illustration of the drainage holes on the shower chair seat, designed to prevent water accumulation and facilitate quick drying.

MAINTENANCE

To ensure the longevity and hygiene of your shower chair:

- **Cleaning:** Clean the chair regularly with mild soap and water. Avoid abrasive cleaners or harsh chemicals that could damage the materials. Rinse thoroughly and wipe dry.
- **Drying:** Allow the chair to air dry completely after each use to prevent mold or mildew buildup.
- **Storage:** Store the folded chair in a dry place away from direct sunlight when not in use.
- **Component Check:** Periodically check all screws, bolts, and leg adjustments to ensure they are tight and secure. Replace any worn or damaged non-slip feet immediately.

TROUBLESHOOTING

If you encounter any issues with your shower chair, refer to the following common solutions:

- **Chair feels unstable or wobbly:**
 - Ensure all four legs are adjusted to the exact same height.
 - Check that the push-pins on each leg are fully engaged in the adjustment holes.
 - Verify that the non-slip rubber feet are clean and making full contact with the floor.
 - Ensure the floor surface is even and free of debris.
- **Difficulty folding/unfolding:**
 - Check for any obstructions in the folding mechanism.
 - Ensure the chair is clean and free of soap residue that might hinder movement.

If problems persist, contact Mobiclinic customer support.

SPECIFICATIONS

Feature	Detail
Brand	Mobiclinic
Model Number	QA-00357/01
Material	Aluminum, Plastic, Stainless Steel
Product Dimensions (D x W x H)	45D x 51.5W x 52.5H cm (unfolded)
Adjustable Height Range	34 cm to 47 cm
Product Weight	2.55 kg
Recommended Maximum Weight	100 kg

DELTA

Technische kenmerken





Kleur
Wit



netto gewicht
2,55kg



Afmetingen
31x54x34-47 cm



Ondersteund gewicht
100kg



Materiaal
Roestvrij staal en PVC



MOBICLINIC
Spanish Leading Company Since 1985

Image: Technical specifications and dimensions of the Mobiclinic Delta Shower Chair, including a 100 kg weight capacity and adjustable height.

WARRANTY AND SUPPORT

For warranty information or technical support, please contact Mobiclinic directly through their official website or the retailer where the product was purchased. Keep your proof of purchase for any warranty claims. Mobiclinic is a leading manufacturer of clinical and hospital furniture, technical aids, and orthopedic equipment, committed to providing quality and trust since 1985.