

McKesson Stretch Vinyl Exam Gloves

McKesson Stretch Vinyl Exam Gloves Instruction Manual

Model: Stretch Vinyl Exam Gloves

1. INTRODUCTION

This manual provides essential information for the proper use, care, and storage of McKesson Stretch Vinyl Exam Gloves. These gloves are designed for general examination purposes and offer protection in various non-sterile environments. Please read these instructions carefully before use.

2. PRODUCT FEATURES

McKesson Stretch Vinyl Exam Gloves are engineered with user comfort and protection in mind. Key features include:

- **Material:** Made from polyvinyl chloride for flexibility and durability.
- **Fit and Dexterity:** Designed to provide a superb fit, enhancing dexterity during tasks.
- **Powder-Free:** Eliminates powder-related irritation and contamination.
- **Latex-Free:** Suitable for individuals with latex sensitivities or allergies.
- **Beaded Cuff:** Features a beaded cuff to prevent the glove from rolling down the wrist during use.
- **Ambidextrous Design:** Can be worn on either hand for convenience.
- **Non-Sterile:** Intended for non-surgical and general examination applications.
- **Disposable:** Single-use only to maintain hygiene and prevent cross-contamination.

Flexible and stretchy to fit comfortably



Image 1: McKesson Stretch Vinyl Exam Gloves box, emphasizing flexible and stretchy comfort.



- Non-sterile
- Beaded cuff
- Ambidextrous

Image 2: Hands wearing the gloves, illustrating key features like non-sterile, beaded cuff, and ambidextrous design.

3. USAGE INSTRUCTIONS

3.1. Preparation

1. **Inspect Packaging:** Before opening, check the glove box for any signs of damage or tampering. Do not use if the packaging is compromised.
2. **Hand Hygiene:** Wash and thoroughly dry your hands before donning gloves. This helps prevent contamination and ensures a better fit.
3. **Select Correct Size:** Ensure you select the appropriate glove size for a comfortable and secure fit. An ill-fitting glove can compromise dexterity and protection.

3.2. Application (Donning)

1. **Open Box:** Open the dispenser box carefully to access the gloves.
2. **Remove Glove:** Grasp a single glove by the cuff. Avoid touching the main body of the glove to maintain cleanliness.
3. **Slide Hand In:** Gently slide your hand into the glove, ensuring your fingers are properly seated in the finger stalls.
4. **Adjust Cuff:** Pull the cuff up to cover your wrist, ensuring a snug fit.
5. **Repeat:** Repeat the process for the other hand.

Powder-free

Disposable



Size medium

Latex-free

Image 3: Box and gloves illustrating key attributes like powder-free, disposable, medium size, and latex-free.

3.3. Removal (Doffing)

1. **Grasp Cuff:** Pinch the outside of one glove near the wrist.
2. **Peel Off:** Peel the glove away from your body, turning it inside out as you remove it. Hold the removed glove in your gloved hand.
3. **Remove Second Glove:** Slide an ungloved finger under the cuff of the remaining glove.
4. **Encase First Glove:** Peel the second glove off, turning it inside out and encasing the first glove within it.
5. **Dispose:** Dispose of the bundled gloves immediately in an appropriate waste receptacle.
6. **Hand Hygiene:** Always wash your hands thoroughly after removing gloves.

4. CARE AND STORAGE

- **Storage Conditions:** Store gloves in their original packaging in a cool, dry place, away from direct sunlight, heat, and fluorescent lighting.
- **Temperature:** Avoid extreme temperatures, as this can affect the integrity and shelf life of the gloves.
- **Chemical Exposure:** Keep gloves away from ozone, oils, greases, petroleum-based products, and solvents, which can degrade the vinyl material.
- **Handling:** Handle boxes carefully to prevent damage to the gloves inside.

5. SAFETY INFORMATION

- **Single Use Only:** These gloves are for single use only. Reusing gloves can lead to cross-contamination and reduced protective barrier integrity.
- **Not for Sterile Procedures:** These are non-sterile gloves and should not be used for sterile surgical procedures or where a sterile barrier is required.
- **Allergies:** While latex-free, individuals with extreme sensitivities should exercise caution. Discontinue use if irritation occurs.
- **Chemical Resistance:** Vinyl gloves offer limited chemical resistance. Consult a chemical resistance chart if working with hazardous chemicals.
- **Puncture Resistance:** Vinyl gloves are not puncture-proof. Exercise caution when handling sharp objects.
- **Proposition 65 Warning:** This product may expose you to chemicals known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

6. SPECIFICATIONS

Attribute	Detail
Product Name	McKesson Stretch Vinyl Exam Gloves
Material	Polyvinyl Chloride (Vinyl)
Sterility	Non-Sterile
Powder Content	Powder-Free
Latex Content	Latex-Free
Color	Ivory
Cuff Type	Beaded Cuff
Hand Compatibility	Ambidextrous
Length	9.84 inches (standard cuff length)
Care Instructions	Hand Wash (referring to hand hygiene before/after use)
Origin	Imported

7. TROUBLESHOOTING

- **Gloves Tearing During Donning:**
 - Ensure hands are completely dry before putting on gloves.
 - Select the correct glove size. Gloves that are too small are more prone to tearing.
 - Avoid sharp fingernails or jewelry that could snag the glove material.
- **Gloves Slipping Off:**
 - Ensure the beaded cuff is properly pulled up over the wrist.
 - Verify that the glove size is appropriate; gloves that are too large may slip.

- **Skin Irritation:**

- Although latex-free and powder-free, some individuals may experience irritation. Discontinue use if irritation persists and consult a healthcare professional.
- Ensure proper hand hygiene before and after glove use.

8. DISPOSAL

Used gloves should be disposed of immediately after single use in accordance with local regulations for waste disposal, especially if they have come into contact with biohazardous materials or chemicals. Do not flush gloves down the toilet.

9. OFFICIAL PRODUCT VIDEO

Your browser does not support the video tag.

Video 1: An official product video from SimplyMedical showcasing the features and benefits of McKesson Stretch Vinyl Exam Gloves. This video provides a visual overview of the gloves' fit, durability, and comfort.

10. MANUFACTURER INFORMATION

Brand: McKesson

Model: Stretch Vinyl Exam Gloves

Manufacturer: McKesson Corporation

For further assistance or inquiries, please refer to the contact information provided on the product packaging or the official McKesson website.