

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Equate](#) /
- › [Equate Extra Comfort Mint Dental Floss User Manual](#)

Equate 701812382

Equate Extra Comfort Mint Dental Floss User Manual

Model: 701812382

INTRODUCTION

The Equate Extra Comfort Mint Dental Floss is designed for effective plaque removal and comfortable use. This manual provides instructions for proper usage and care to ensure optimal oral hygiene.





Image: Packaging of Equate Extra Comfort Mint Dental Floss. The image displays the product in its retail packaging, featuring two green floss dispensers. Each dispenser contains 43.7 yards (40 meters) of mint-flavored dental floss, with the total package containing 87.4 yards (80 meters). The packaging highlights features such as "Gentle on gums," "Effectively removes plaque," and "Slides easily between tight teeth."

SETUP

The Equate Extra Comfort Mint Dental Floss comes ready for use. No assembly or special setup is required. Simply remove the floss dispenser from its packaging.

OPERATING INSTRUCTIONS

Follow these steps for effective dental flossing:

1. Break off about 18 inches (45 cm) of floss.
2. Wind most of the floss around your middle fingers, leaving about one or two inches of floss to work with.
3. Hold the floss tightly between your thumbs and forefingers.
4. Gently guide the floss between your teeth using a rubbing motion.

5. Once the floss reaches the gum line, curve it into a C-shape against one tooth.
6. Gently slide the floss into the space between the gum and the tooth.
7. Hold the floss firmly against the tooth and move it up and down, scraping the side of the tooth.
8. Repeat this process for the other side of the tooth and for all teeth, using a new section of floss for each tooth.
9. Avoid snapping or forcing the floss to prevent injury to your gums.

MAINTENANCE

To maintain the quality of your dental floss:

- Store the floss dispenser in a cool, dry place away from direct sunlight and extreme temperatures.
- Keep the dispenser cap closed when not in use to protect the floss from dust and contaminants.
- Dispose of used floss properly in a waste bin; do not flush down the toilet.

TROUBLESHOOTING

| Problem | Possible Cause | Solution |
|---|--|--|
| Floss shreds or breaks easily. | <ul style="list-style-type: none"> • Excessive force applied. • Very tight contact points between teeth. • Rough edges on teeth or dental work. | <ul style="list-style-type: none"> • Use a gentler, rubbing motion to guide the floss. • Consult your dentist if you consistently experience shredding, as it may indicate very tight contacts or rough surfaces that need professional attention. • Ensure you are not using too long a section of floss, which can weaken it. |
| Difficulty sliding floss between teeth. | <ul style="list-style-type: none"> • Tight contact points. • Incorrect technique. | <ul style="list-style-type: none"> • Use a gentle sawing motion to work the floss through tight contacts. • Ensure the floss is held taut. • This floss is designed for "Extra Comfort" and "Slides easily between tight teeth," so if issues persist, consult a dental professional. |

SPECIFICATIONS

| Feature | Detail |
|--------------------|--------------------------|
| Brand | Equate |
| Model Number | 701812382 |
| Item Form | String |
| Flavor | Mint |
| Material | Nylon |
| Package Dimensions | 6.5 x 3.25 x 0.05 inches |
| Item Weight | 0.81 ounces |

| Feature | Detail |
|----------------------|----------------------|
| Manufacturer | Wal-Mart Store, Inc. |
| Date First Available | October 12, 2014 |

SUPPORT AND LEGAL INFORMATION

For questions or concerns regarding your Equate Extra Comfort Mint Dental Floss, please refer to the contact information provided on the product packaging or visit the official Equate website.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is for oral hygiene use only.