

Evenflo 849409000514

Evenflo Advanced Disposable Nursing Pads User Manual

Model: 849409000514

INTRODUCTION

Thank you for choosing Evenflo Advanced Disposable Nursing Pads. These pads are designed to provide superior absorbency and comfort for breastfeeding mothers, ensuring dryness and discretion. This manual provides instructions for proper use and care.

KEY FEATURES

Evenflo Advanced Disposable Nursing Pads incorporate advanced technology for optimal performance.

- **Best-in-class absorbency technology:** Designed for exceptional absorption and performance.
- **SuperiorShield™ protection:** Guards against leaks, with soft material for discreet and silent movement.
- **FreeFlow™ breathability:** Helps keep skin dry and comfortable.
- **Embossed channels:** Aid in directing milk flow onto the pad and away from the nipple.
- **Ultra-soft material:** Gentle and comfortable against sensitive skin.
- **Secure positioning:** Features two adhesive strips for reliable placement.
- **Individually wrapped:** For hygiene and convenience on the go.



Image: Packaging of Evenflo Advanced Disposable Nursing Pads, highlighting key features like absorbency, protection, and breathability, with a count of 40 individually wrapped pads.

USAGE INSTRUCTIONS

Setup and Application

1. **Preparation:** Ensure your hands are clean and dry before handling the nursing pad.
2. **Unwrap:** Open the individual wrapper to remove one nursing pad.
3. **Positioning:** Place the pad inside your bra, centering it over your nipple. The contoured shape is designed to fit discreetly.
4. **Secure:** Remove the backing from the two adhesive strips on the back of the pad. Press the pad firmly against your bra to secure it in place.

Daily Use

Wear the nursing pads as needed throughout the day and night to absorb breast milk leaks and maintain dryness.

Maintenance and Disposal

- **Replacement:** Change nursing pads frequently to maintain hygiene and prevent skin irritation. Replace immediately if they become wet or soiled.
- **Disposal:** Dispose of used nursing pads in a waste bin. Do not flush.

TROUBLESHOOTING

If you experience any discomfort or irritation while using the nursing pads, discontinue use and consult a healthcare professional. Ensure pads are changed regularly to prevent moisture buildup, which can lead to skin issues.

SPECIFICATIONS

Product Name	Evenflo Advanced Disposable Nursing Pads
Model Number	849409000514
Quantity	40 Count
Dimensions	12.8 x 14.5 x 18.3 inches (packaging)
Weight	8.48 ounces (packaging)
Manufacturer	Evenflo

WARRANTY AND SUPPORT

For specific warranty information or customer support regarding Evenflo products, please refer to the official Evenflo website or contact their customer service directly. Contact details are typically available on the product packaging or the manufacturer's website.

Important Note: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is for personal hygiene and comfort.

