



[Manuals.plus](#) /

> [Medline](#) /

> Medline NON256025H Iodoform Gauze Packing Strip User Manual

Medline NON256025H

Medline NON256025H Iodoform Gauze Packing Strip User Manual

Brand: Medline | Model: NON256025H

1. INTRODUCTION

This manual provides essential information and instructions for the safe and effective use of the Medline NON256025H Iodoform Gauze Packing Strip. This product is designed for professional medical use in wound care. Please read all instructions carefully before use.

2. PRODUCT DESCRIPTION

The Medline NON256025H is a sterile iodoform gauze packing strip. It is made from gauze material and is impregnated with iodoform. The strip measures 2 inches by 15 feet (approximately 5.1 cm x 4.6 m). This product is intended for packing wounds to facilitate drainage and promote healing.

FOR YOUR PROTECTION
SEALING FOR YOUR PROTECTION
NOTION

CURAD[®]

WE HELP HEAL/NOUS AIDONS VRAIMENT
AYUDAMOS A CURAR

STERILE

IODOFORM PACKING STRIP

MÉCHE IODOFORMÉE
TIRA DE VENDAJE IODOFORM

1
IODOFORM PACKING STRIP
MÉCHE IODOFORMÉE
TIRA DE VENDAJE IODOFORM

REF NON256025

2 in x 5 yd
(5.08 cm x 4.57 m)

latexfree
PRODUCT
Sans latex
Sin látex

Do not autoclave
Ne pas autoclaver
No autoclave
Do not use
Ne pas utiliser
No use
LOT
MEDLINE
Made in USA
Curad[®] is a
Manufactured
Medline
Made in USA
1-800-433-1111

Figure 1: Medline NON256025H Iodoform Gauze Packing Strip. The image shows a roll of yellow-tinted gauze, indicating the iodoform impregnation.

3. INTENDED USE

The Medline NON256025H Iodoform Gauze Packing Strip is recommended for:

- Wound care
 - Packing open wounds
 - Managing draining wounds
 - Treating boils and abscesses
 - Addressing fistulas
 - Packing tunneling wounds
-

4. WARNINGS AND PRECAUTIONS

- **Sterility:** This product is sterile. Do not use if the packaging is opened or damaged.
 - **Single Use:** This product is for single use only. Do not re-sterilize or reuse.
 - **Allergies:** Discontinue use if signs of irritation, redness, or allergic reaction appear. Consult a healthcare professional.
 - **Professional Use:** This product should be used by or under the direction of a healthcare professional.
 - **Storage:** Store in a cool, dry place away from direct sunlight.
 - **Disposal:** Dispose of used product according to local medical waste regulations.
-

5. SETUP

1. **Hand Hygiene:** Thoroughly wash and sanitize hands. Don sterile gloves if required by protocol.
 2. **Wound Assessment:** Assess the wound to determine the appropriate length and packing technique.
 3. **Prepare Area:** Ensure the wound area is clean and dry according to standard wound care procedures.
 4. **Open Packaging:** Carefully open the sterile packaging without contaminating the gauze strip.
-

6. OPERATING INSTRUCTIONS (APPLICATION)

Follow these steps for applying the iodoform gauze packing strip:

1. Gently grasp one end of the sterile iodoform gauze strip.
2. Carefully pack the strip into the wound cavity, ensuring it fills the space without excessive pressure.
3. Avoid over-packing, which can impede circulation and healing.
4. Leave a small portion of the strip exposed outside the wound for easy removal.
5. Secure the packing with an appropriate secondary dressing, such as a sterile gauze pad and adhesive tape, to hold it in place and absorb exudate.

6. Document the procedure, including the length of gauze used and the date of application.

Removal: Remove the packing strip gently, following healthcare professional instructions, typically during dressing changes. If the strip adheres, moisten it with sterile saline solution before removal.

7. MAINTENANCE AND STORAGE

This product requires no maintenance. For storage:

- Store in its original, unopened, and undamaged packaging.
 - Keep in a cool, dry environment.
 - Protect from direct sunlight and extreme temperatures.
-

8. TROUBLESHOOTING

If any of the following issues occur, discontinue use and consult a healthcare professional:

- Signs of infection (increased redness, swelling, pain, pus).
 - Allergic reaction or skin irritation around the wound.
 - Unexpected bleeding or increased wound exudate.
 - The packing strip becomes dislodged or difficult to remove.
-

9. SPECIFICATIONS

Brand	Medline
Model Number	NON256025H
Material	Gauze
Size	2 inches x 15 feet (5.1 cm x 4.6 m)
Sterility Rating	Sterile
Item Format	Strips
Recommended Uses	Wound care, open wounds, draining wounds, boils, abscesses, fistulas, tunneling wounds
Manufacturer	Medline Industries Healthcare
UPC	080196733067
GTIN	10080196733064

10. DISPOSAL

Used iodoform gauze packing strips are considered medical waste. Dispose of all used materials in accordance with local, state, and federal regulations for biohazardous waste. Do not dispose of in general household waste.

11. WARRANTY AND SUPPORT

For information regarding product warranty or technical support, please contact Medline customer service directly. Contact details can typically be found on the product packaging or the official Medline website.
