

Toorx AHF-075

Toorx AHF-075 Neoprene Wrist Protection User Manual

Model: AHF-075

1. PRODUCT OVERVIEW

The Toorx AHF-075 Neoprene Wrist Protection is designed to provide comfortable and ergonomic support for your wrist. It offers a degree of stability and gentle compression, which can help alleviate discomfort. The wrist protection is made from durable neoprene material and features an adjustable strap for easy application and a secure fit.

Key Features:

- Comfortable and ergonomic design for wrist protection.
- Provides support and stability through light compression.
- Easy to apply and adjust with a secure strap closure.
- Constructed from high-quality neoprene material.

2. SAFETY INFORMATION

Please read all instructions carefully before using the Toorx AHF-075 Neoprene Wrist Protection. Improper use may lead to discomfort or injury.

- **Consult a Physician:** If you experience persistent pain, swelling, or discomfort, consult a medical professional before using this product.
- **Proper Fit:** Ensure the wrist protection is snug but not overly tight. Excessive compression can restrict blood circulation.
- **Discontinue Use:** If you experience numbness, tingling, increased pain, or skin irritation while wearing the product, remove it immediately and consult a physician.
- **Allergies:** If you have known allergies to neoprene or similar materials, exercise caution or consult a doctor before use.
- **Not for Open Wounds:** Do not apply the wrist protection over open wounds or broken skin.
- **Children:** Keep out of reach of small children.

3. APPLICATION INSTRUCTIONS

Follow these steps to correctly apply your Toorx AHF-075 Neoprene Wrist Protection:

1. **Position the Thumb Loop:** Slide your thumb through the smaller loop of the wrist protection. Ensure the main body of the support is positioned over the back of your hand and wrist.



Image 1: Initial placement of the wrist protection with the thumb through the loop.

2. **Wrap Around the Wrist:** Wrap the main section of the wrist protection around your wrist. The material should cover the wrist joint comfortably.
3. **Secure the Strap:** Use the adjustable strap closure to secure the wrist protection. Pull the strap firmly but not too tightly, ensuring a comfortable and supportive fit. The strap allows for easy

adjustment to achieve the desired level of compression and support.



Image 2: The wrist protection fully applied and secured around the wrist.

4. **Check for Comfort:** Move your hand and wrist to ensure the protection feels secure and comfortable, without restricting movement or causing discomfort. Adjust the strap as needed.

4. USAGE GUIDELINES

The Toorx AHF-075 Neoprene Wrist Protection is intended for providing support and gentle compression to the wrist. It can be used during physical activities, for general wrist support, or to help alleviate minor wrist discomfort.

- Wear the wrist protection for the duration recommended by your healthcare professional or as needed for comfort and support.
- Avoid wearing the product for extended periods if it causes discomfort or skin irritation.
- Ensure the product does not interfere with your normal activities or cause any pain.

5. CARE AND MAINTENANCE

Proper care will extend the life of your Toorx AHF-075 Neoprene Wrist Protection.

- **Hand Wash Only:** Hand wash in cold water with a mild soap.
- **Do Not Bleach:** Avoid using bleach or harsh detergents.
- **Air Dry:** Lay flat to air dry. Do not tumble dry, iron, or dry clean.
- **Storage:** Store in a cool, dry place away from direct sunlight.

6. TROUBLESHOOTING

If you encounter any issues while using the Toorx AHF-075 Neoprene Wrist Protection, refer to the following common solutions:

- **Discomfort or Too Tight:** Loosen the adjustable strap. Ensure the product is not applied too tightly, which can restrict circulation.
- **Not Enough Support:** Tighten the adjustable strap slightly until the desired level of compression and support is achieved, without causing discomfort.
- **Skin Irritation:** Ensure the product is clean. If irritation persists, discontinue use and consult a medical professional. This could indicate an allergy or improper fit.
- **Product Slipping:** Re-adjust the strap to ensure a snug fit. The neoprene material should provide good grip, but an improper fit can lead to slipping.

If problems persist, please contact customer support or consult a medical professional.

7. SPECIFICATIONS

Model Number:	AHF-075
Brand:	Toorx
Material:	Neoprene
Product Dimensions:	44 x 36 x 36 cm (approximate packaging dimensions); 100 g (weight)
ASIN:	B00VAOZXF2
Global Trade Item Number (GTIN):	08029975991856

8. WARRANTY AND SUPPORT

Warranty Information:

The Toorx AHF-075 Neoprene Wrist Protection comes with a **2-year warranty** from the date of purchase. This warranty covers manufacturing defects under normal use. It does not cover damage resulting from misuse, accidents, unauthorized modifications, or normal wear and tear.

Customer Support:

For any questions, concerns, or warranty claims regarding your Toorx AHF-075 Neoprene Wrist Protection, please contact your retailer or the official Toorx customer support. Please have your proof of purchase and model number (AHF-075) ready when contacting support.

For more information, visit the official Toorx website: www.toorx.it (Note: This is a placeholder link, please verify the actual Toorx website for support).