

ResMed SlimLine

ResMed SlimLine Replacement Tubing Instruction Manual

Model: SlimLine

1. PRODUCT OVERVIEW

The ResMed SlimLine Replacement Tubing is designed to provide a comfortable and efficient continuous positive airway pressure (CPAP) therapy experience. This lightweight and flexible tubing is compatible with ResMed Air10 and S9 series machines, ensuring a reliable connection for your respiratory care needs.

- **Lightweight Design:** The SlimLine Non-Heated Tubing is sleek and light for a more comfortable CPAP experience.
- **Compatibility Assured:** This CPAP accessory seamlessly fits ResMed Air10 and S9 machines, ensuring a reliable connection.
- **Slim & Flexible:** Engineered for space-saving, the slim profile of this CPAP tubing offers greater flexibility and convenience.
- **Durable CPAP Supplies:** Constructed to be long-lasting, the authentic ResMed brand SlimLine Air Tubing is a dependable component of your CPAP equipment.



Image 1: ResMed Official Logo.



Image 2: The ResMed SlimLine Replacement Tubing, coiled for storage.

2. SETUP AND INSTALLATION

Follow these steps to correctly install your ResMed SlimLine Replacement Tubing:

1. **Disconnect Old Tubing:** Gently pull the old tubing straight out from the air outlet of your CPAP machine.
2. **Inspect New Tubing:** Ensure the new SlimLine tubing is free from any damage or obstructions.
3. **Connect to Machine:** Locate the end of the new tubing with the small chip or slot. Align this chip with the corresponding slot on the air outlet of your ResMed Air10 or S9 machine. Push the tubing in firmly until it is flush and secure.
4. **Connect to Mask:** Attach the other end of the tubing to your CPAP mask. Different mask types may have varying connection mechanisms; ensure a secure fit according to your mask's instructions.



Image 3: ResMed SlimLine Tubing connected to a CPAP machine.

3. OPERATING INSTRUCTIONS

Once the SlimLine tubing is securely connected to your CPAP machine and mask, you can begin your therapy. Refer to your specific ResMed CPAP machine's user manual for detailed operating instructions, including how to start and stop therapy, adjust settings, and monitor your progress.

- Ensure all connections are tight before starting therapy to prevent air leaks.
- Position the tubing to avoid kinks or obstructions that could impede airflow.
- If using a humidifier, ensure the water chamber is filled to the appropriate level.

4. MAINTENANCE AND CLEANING

Regular cleaning of your SlimLine tubing is essential for hygiene and optimal performance. Follow these general guidelines:

- **Daily Cleaning:** Disconnect the tubing from your machine and mask. Wash the tubing in warm water with a mild soap. Rinse thoroughly with clean water.
- **Drying:** Hang the tubing in a clean, dry place away from direct sunlight to air dry completely before reattaching. Ensure no water remains inside the tubing.
- **Inspection:** Regularly inspect the tubing for any signs of wear, tears, or damage. Replace the tubing if any damage is observed to maintain effective therapy.
- **Replacement:** It is recommended to replace CPAP tubing every 3-6 months, or as advised by your healthcare provider or equipment manufacturer, to ensure optimal hygiene and performance.

5. TROUBLESHOOTING

If you encounter issues with your SlimLine tubing, consider the following:

- **Air Leaks:** Check all connections between the tubing, machine, and mask. Ensure they are securely fitted. Inspect the tubing for any visible holes or cracks.
- **Reduced Airflow:** Ensure the tubing is not kinked or obstructed. Check if the tubing is clean and free from debris.
- **Condensation (Rainout):** If you experience water buildup in the tubing, consider using a heated humidifier with a heated tube (if your machine supports it and you have the appropriate heated tubing), or a tubing wrap/cover to insulate the air.
- **Unusual Noises:** Check for loose connections or damage to the tubing.

If problems persist, consult your healthcare provider or the ResMed customer support.

6. SPECIFICATIONS

| Feature | Detail |
|----------------------|--|
| Product Dimensions | 12.28 x 10.08 x 2.36 inches; 3.53 ounces |
| Manufacturer | ResMed |
| ASIN | B00YFVG9EG |
| Date First Available | May 28, 2015 |

7. SAFETY INFORMATION

Important Safety Warning Regarding Magnets:

Keep the mask magnets at a safe distance of at least 6 inches (150 mm) away from implants or medical devices that may be adversely affected by magnetic interference. This warning applies to you or anyone in close physical contact with your mask. The magnets are in the frame and lower headgear clips, with a magnetic field strength of up to 400mT. When worn, they connect to secure the mask but may inadvertently detach while asleep.

Legal Disclaimer:

Masks with magnetic components are contraindicated for use by patients where they, or anyone in close physical contact while using the mask, have the following:

- Active medical implants that interact with magnets (i.e., pacemakers, implantable cardioverter defibrillators (ICD), neurostimulators, cerebrospinal fluid (CSF) shunts, insulin/infusion pumps)
- Metallic implants/objects containing ferromagnetic material (i.e., aneurysm clips/flow disruption devices, embolic coils, stents, valves, electrodes, implants to restore hearing or balance with implanted magnets, ocular implants, metallic splinters in the eye)

8. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the documentation provided with your original ResMed CPAP machine or visit the official ResMed website. Keep your purchase receipt for any warranty claims.