

Chicco 1166000500

Chicco Mamma Donna Postpartum Modelling Sheath - Instruction Manual

Model: 1166000500

1. PRODUCT OVERVIEW

The Chicco Mamma Donna Postpartum Modelling Sheath is designed to provide gentle support and shaping for the abdominal area after childbirth. This adjustable garment is crafted for comfort and effectiveness, aiding in the recovery process.

Key Features:

- Adjustable fit for personalized comfort and support.
- Designed for postpartum use.

2. WEARING INSTRUCTIONS

To ensure optimal comfort and effectiveness, follow these steps when wearing your Chicco Mamma Donna Postpartum Modelling Sheath:

1. **Preparation:** Ensure the sheath is clean and dry before wearing.
2. **Step In:** Step into the sheath and pull it up over your hips and abdomen.
3. **Positioning:** Adjust the sheath so it sits comfortably around your waist and covers the desired abdominal area. The high-waisted design is intended to provide comprehensive support.
4. **Adjustability:** Utilize the adjustable features to achieve a snug yet comfortable fit. Avoid overtightening, which can cause discomfort or restrict circulation. The garment should provide gentle compression without feeling restrictive.
5. **Smooth Out:** Smooth out any wrinkles or bunched fabric to ensure an even distribution of support and a seamless look under clothing.

If you experience any discomfort, readjust the sheath or consider consulting the sizing guide to ensure you have the correct size.



Image 1: Side view of the Chicco Mamma Donna Postpartum Modelling Sheath. This image shows the high-waisted design and the gentle contouring it provides to the abdominal area.



Image 2: Front view of the Chicco Mamma Donna Postpartum Modelling Sheath. This image highlights the front panel design and the overall coverage.



Image 3: Back view of the Chicco Mamma Donna Postpartum Modelling Sheath. This image illustrates the rear coverage and the seamless appearance.

3. CARE AND MAINTENANCE

Proper care will extend the life and maintain the effectiveness of your Chicco Mamma Donna Postpartum Modelling Sheath.

- **Washing:** Hand Wash Only. Use mild detergent and cold water.
- **Drying:** Do not tumble dry. Lay flat to dry in a shaded area to preserve elasticity and shape.
- **Bleaching:** Do not bleach.
- **Ironing:** Do not iron.
- **Dry Cleaning:** Do not dry clean.

4. SPECIFICATIONS AND SIZING

Refer to the following specifications and sizing chart to ensure the correct fit for your Chicco Mamma Donna Postpartum Modelling Sheath.

Product Specifications:

- **Brand:** Chicco
- **Model Number:** 1166000500
- **Product Dimensions:** 11.6 x 5 x 21.3 cm; 99.79 g
- **Age Range:** Adult
- **Number of Items:** 1
- **Special Feature:** Adjustable



Sizing Chart:

Circonferenza Fianchi	Circonferenza Vita	Taglia Internazionale	Taglia
-	-	-	1
da 90 a 94 cm	da 60 a 64 cm	S	2
da 95 a 99 cm	da 65 a 69 cm	M	3
da 100 a 103 cm	da 70 a 74 cm	L	4
da 104 a 108 cm	da 75 a 79 cm	XL	5
da 109 a 112 cm	da 80 a 84 cm	-	6
da 113 a 117 cm	da 85 a 89 cm	-	7

Image 4: Sizing chart for the Chicco Mamma Donna Postpartum Modelling Sheath. This chart provides measurements for hip and waist circumference to help determine the appropriate size.

Chicco Mamma Donna Sizing Guide

Hip Circumference	Waist Circumference	International Size	Size
-	-	-	1

 <p>The image shows the Chicco Hoopla baby bouncer, a blue and grey portable seat with a canopy and wheels. The Chicco logo is visible in the bottom right corner of the image.</p>	<p>Chicco Hoopla Baby Bouncer: User Manual & Safety Guide</p> <p>Comprehensive user manual for the Chicco Hoopla baby bouncer. Learn about safe assembly, usage modes (reclined cradle and chair), adjustable features, and maintenance for this versatile baby product.</p>
 <p>The image shows a row of various Chicco NaturalFeeling baby bottle teats, each with a different shape and flow rate.</p>	<p>Chicco NaturalFeeling Teats - Instructions for Use, Cleaning, and Sterilization</p> <p>Comprehensive guide for Chicco NaturalFeeling baby bottle teats, covering usage instructions, cleaning procedures, sterilization methods, and flow rate descriptions for safe and optimal feeding.</p>