

Manuals.plus /

- › Omron /
- › Omron 7 Series Wireless Upper Arm Blood Pressure Monitor User Manual

Omron BP7350

Omron 7 Series Wireless Upper Arm Blood Pressure Monitor

Model: **BP7350** | Brand: **Omron**

1. INTRODUCTION AND OVERVIEW

The Omron 7 Series wireless upper arm home blood pressure monitor (BP7350) is designed for accurate and convenient blood pressure measurement at home. It features a horizontally designed dual-display monitor and can store up to 120 readings (60 readings per user for two users). The device offers Bluetooth connectivity to the free Omron Connect app, allowing for unlimited memory, graphical trend tracking, and support for unlimited users. Omron is a leading brand in home blood pressure monitoring, recommended by doctors and pharmacists for its clinically accurate readings.



Figure 1: Omron 7 Series Wireless Upper Arm Blood Pressure Monitor and ComFit Cuff.

2. SAFETY INFORMATION

Always consult with your physician or healthcare provider before making any decisions about your health. This device is intended for home use to measure blood pressure and pulse rate. It is not intended to be a diagnostic device. Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

- Do not use this device for any purpose other than measuring blood pressure.
- Keep out of reach of children.
- Avoid strong electromagnetic fields near the device.
- Do not disassemble or attempt to repair the device.

3. PACKAGE CONTENTS

- Omron 7 Series Wireless Upper Arm Blood Pressure Monitor (BP7350)
- Easy Wrap ComFit Cuff (fits arms 9" to 17" in circumference)
- 4 AA Batteries
- Instruction Manual (this document)
- Quick Start Guide

4. SETUP

4.1 Installing Batteries

1. Open the battery cover on the back of the monitor.
2. Insert the 4 AA batteries, ensuring correct polarity (+/-).
3. Close the battery cover securely.

4.2 Connecting the Arm Cuff

- Insert the air plug of the arm cuff firmly into the air jack on the left side of the monitor.
- Ensure the connection is secure to prevent air leakage during measurement.

4.3 Pairing with Omron Connect App (Optional)

For enhanced tracking and unlimited storage, you can sync your monitor with the Omron Connect app.

1. Download the free Omron Connect app from your device's app store (iOS and Android compatible).
2. Enable Bluetooth on your smartphone or tablet.
3. Follow the in-app instructions to pair your Omron BP7350 monitor.

5. OPERATING INSTRUCTIONS

5.1 Proper Cuff Placement

Proper cuff placement is crucial for accurate readings. The Easy Wrap ComFit Cuff is designed for simple application.

1. Sit upright with your back supported and feet flat on the floor.
2. Place the cuff on your left upper arm, approximately 1-2 cm (0.5-1 inch) above your elbow. The air tube should run down the inside of your arm.
3. Ensure the cuff is snug but not too tight. The indicator on the cuff helps confirm correct placement.
4. Rest your arm on a table so the cuff is at heart level.

5.2 Taking a Measurement

1. Select User 1 or User 2 using the user selection switch if multiple users are tracking readings.
2. Press the **START/STOP** button. The cuff will begin to inflate automatically.
3. Remain still and do not talk during the measurement.
4. Once the measurement is complete, the cuff will deflate, and your systolic, diastolic, and pulse readings will be displayed on the dual LCD screen.

5.3 Understanding Your Readings

- **Dual LCD Screen:** Displays your last reading next to your current reading for convenient comparison.
- **Irregular Heartbeat Symbol:** Alerts you to the number of irregular heartbeats (up to 3) detected during measurement.
- **Hypertension Indicator:** Displays if your systolic measurement is 130 mmHg or above, and/or diastolic measurement is 80 mmHg or above.

6. PRODUCT FEATURES

- **Dual LCD Screen:** Displays your last reading next to your current reading for convenient comparison; stores 120 total blood pressure readings for 2 users (60 per user); includes our easy wrap ComFit cuff that fits arms 9" to 17" in circumference.
- **Enhanced Irregular Heartbeat Symbol:** Alerts you to the number of irregular heartbeats (up to 3) during measurements, while the Hypertension Indicator displays if systolic measurement is 130 mmHg or above, and/or diastolic measurement is 80 mmHg or above.
- **Bluetooth Connectivity:** With the free Omron Connect app (on select iOS & Android devices; visit [omronhealthcare.com/connect](https://www.omronhealthcare.com/connect) for compatibility) stores unlimited readings and users; works with Amazon Alexa enabled devices using the Omron health skill.
- **Trusted Brand:** Omron is the #1 recommended home blood pressure monitor brand by doctors and pharmacists for clinically accurate home monitoring, and the #1 selling manufacturer of home blood pressure monitors for over 40 years.
- **Important Note:** Blood pressure varies constantly due to many factors including stress, time of day, and how you wrap the cuff, may affect your blood pressure.

7. MAINTENANCE AND CARE

- Clean the monitor and cuff with a soft, dry cloth. Do not use abrasive cleaners or immerse in water.
- Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.
- Replace all 4 AA batteries when the low battery indicator appears on the display.

8. TROUBLESHOOTING

If you encounter issues with your Omron BP7350 monitor, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
"E1" Error on Display	Cuff not properly attached or air leakage.	Ensure the cuff's air plug is securely inserted into the monitor. Re-wrap the cuff snugly on your arm.
Inaccurate Readings	Improper cuff placement, movement/talking during measurement, incorrect body position, or recent physical activity.	Review "Proper Cuff Placement" and "Taking a Measurement" sections. Rest for 5 minutes before re-measuring.
Monitor Does Not Turn On	Dead or incorrectly installed batteries.	Check battery polarity or replace all 4 AA batteries.

For further assistance, please refer to the full user manual (PDF) available on the Omron website or contact Omron customer support.

9. SPECIFICATIONS

- **Model Number:** BP7350
- **Product Dimensions:** 7.51"L x 4.13"W x 3.11"Th
- **Item Weight:** 1.9 Pounds
- **Cuff Size:** 22-42 cm (9-17 inches)
- **Power Source:** 4 AA batteries (included)
- **Memory:** 120 readings (60 per user for 2 users)
- **Connectivity:** Bluetooth
- **UPC:** 073796267353

10. WARRANTY AND SUPPORT

The Omron 7 Series Wireless Upper Arm Blood Pressure Monitor unit comes with a **Five-year limited warranty**. The ComFit cuff has a **One-year limited warranty**.

For product support, FAQs, and additional resources, please visit the official Omron website or contact their customer service. You can also utilize the Omron Connect app for further assistance and to manage your health data.

Visit the Omron Store