

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Hartmann](#) /

› [Hartmann Econo-Wrap Reinforced Elastic Bandage \(Model 33400000\) Instruction Manual](#)

Hartmann 33400000

Hartmann Econo-Wrap Reinforced Elastic Bandage (Model 33400000) Instruction Manual

1. PRODUCT OVERVIEW

The Hartmann Econo-Wrap Reinforced Elastic Bandage, Model 33400000, is a professional healthcare product designed for various bandaging needs. This product is latex-free and sterile, ensuring suitability for sensitive skin and medical environments. Each bandage measures 4 inches in width and 4.5 yards in length, providing ample material for secure application. This product is typically supplied in a pack containing 60 individual bandages.



This image displays the Hartmann Econo-Wrap Reinforced Elastic Bandages. In the background, two white cardboard boxes with blue labeling are visible, indicating packaging for multiple bandages. In the foreground, several individual rolled bandages, varying slightly in size, are arranged. These bandages are off-white or beige in color, characteristic of elastic medical wraps.

2. PRODUCT FEATURES

- **Latex-Free:** Suitable for individuals with latex sensitivities, minimizing allergic reactions.
- **Sterile:** Ensures hygienic application in medical settings, reducing the risk of infection.
- **Reinforced Elasticity:** Provides balanced stretch and compression for effective support and comfort.
- **Optimal Dimensions:** Each bandage measures 4 inches (approximately 10.16 cm) in width and 4.5 yards (approximately 4.11 meters) in length.
- **Convenient Packaging:** Available in a pack of 60 individual bandages for extended use.

3. INTENDED USE

The Hartmann Econo-Wrap Reinforced Elastic Bandage is intended for use in various medical and first aid applications, including but not limited to:

- Providing support and compression for sprains and strains.

- Assisting in the reduction of swelling.
- Securing dressings, splints, or cold/hot packs.
- Offering post-operative support to surgical sites.

4. SETUP AND PREPARATION

1. **Inspect Packaging:** Before use, carefully inspect the individual bandage packaging for any signs of damage or compromise to sterility. Do not use if the package is opened or damaged.
2. **Clean and Dry Skin:** Ensure the area of application is clean, dry, and free from oils or lotions.
3. **Position Comfortably:** Position the affected body part in a comfortable and stable position to facilitate proper bandage application.

5. OPERATING INSTRUCTIONS (APPLICATION)

Follow these steps for proper application of the elastic bandage:

1. **Unroll:** Unroll a small section of the bandage to begin.
2. **Start Wrapping:** Begin wrapping from the point furthest from the body (e.g., fingers for a hand injury, toes for a foot injury). This helps promote proper circulation.
3. **Overlap Layers:** Overlap each layer of the bandage by approximately half its width. This ensures even compression and coverage.
4. **Apply Even Pressure:** Apply even, firm pressure as you wrap, but avoid wrapping too tightly. The bandage should feel supportive, not constricting. Ensure two fingers can easily slide under the bandage.
5. **Secure End:** Once the desired area is covered, secure the end of the bandage using appropriate clips or medical tape (not included).
6. **Monitor:** After application, regularly monitor the wrapped area for signs of discomfort, numbness, tingling, discoloration (pale or blue skin), or increased swelling. If any of these symptoms occur, remove the bandage immediately and reapply with less tension, or seek medical advice.

6. CARE AND MAINTENANCE

- **Single-Use Product:** The Hartmann Econo-Wrap Reinforced Elastic Bandage is designed for single use only. Dispose of the used bandage properly according to local regulations.
- **Storage:** Store unused bandages in their original, unopened packaging in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Do Not Reuse or Re-sterilize:** Do not attempt to wash, re-sterilize, or reuse the bandage, as this can compromise its sterility, elasticity, and effectiveness.

7. TROUBLESHOOTING

- **Bandage feels too tight:** If the wrapped area feels numb, tingles, or appears discolored (pale or blue), the bandage is too tight. Immediately remove the bandage and reapply with less tension.
- **Skin irritation or allergic reaction:** Although latex-free, some individuals may experience skin irritation. If irritation, rash, or itching occurs, discontinue use and consult a healthcare professional.
- **Bandage unravels or loses support:** Ensure the end of the bandage is securely fastened with clips or tape. If the bandage loses elasticity prematurely, it may be due to improper storage or extended use beyond its intended single-use purpose.

8. SPECIFICATIONS

Model Number	33400000
Product Type	Reinforced Elastic Bandage
Material	Latex-Free
Sterility	Sterile
Dimensions (per bandage)	4 inches (W) x 4.5 yards (L)
Package Quantity	60 bandages per pack
Manufacturer	Hartmann USA, Inc.
Country of Origin	China

9. WARNINGS AND PRECAUTIONS

- **External Use Only:** This product is for external application only.
- **Open Wounds:** Do not apply the elastic bandage directly to open wounds without a sterile dressing underneath.
- **Medical Consultation:** Always consult a healthcare professional for proper diagnosis and treatment of injuries or medical conditions.
- **Children:** Keep out of reach of children.
- **Adverse Reactions:** If signs of infection (e.g., increased pain, redness, swelling, pus), severe discomfort, or allergic reactions occur, remove the bandage and seek immediate medical attention.

10. MANUFACTURER INFORMATION

Manufacturer: Hartmann USA, Inc.

ASIN: B01DME3VOC

For further inquiries or support, please refer to the manufacturer's official contact channels.