



[Manuals.plus](#) /

> [A+Elite Medical](#) /

> Welch Allyn Thermoscan PRO 6000 Ear Thermometer Probe Covers 200/BX Instruction Manual

**A+Elite Medical AE06000-800-200**

# Welch Allyn Thermoscan PRO 6000 Ear Thermometer Probe Covers Instruction Manual

## INTRODUCTION

---

This instruction manual provides guidance for the proper use and handling of Welch Allyn Thermoscan PRO 6000 Ear Thermometer Probe Covers. These covers are designed to ensure hygienic and accurate temperature readings when used with compatible Welch Allyn Thermoscan ear thermometers.

These probe covers are compatible with Welch Allyn Thermoscan PRO 6000 series thermometers, and also work with older models such as the PRO 4000 and PRO 4250.



Image: A box containing 200 Welch Allyn ThermoScan PRO 6000 Ear Thermometer Probe Covers. The box is white and blue, displaying the product name and quantity.

## SETUP AND INSTALLATION

Follow these steps to properly install a probe cover onto your Welch Allyn ThermoScan ear thermometer:

1. **Unpack the Dispenser Box:** Open the main dispenser box. It contains 10 smaller boxes, each holding 20 probe covers.
2. **Access a Single Box:** Remove one of the smaller 20-count boxes from the dispenser.

3. **Prepare the Thermometer:** Ensure your Welch Allyn Thermoscan ear thermometer is clean and ready for use.
4. **Attach a Probe Cover:** The individual probe cover boxes are designed for easy, hygienic attachment. Without touching the probe cover directly, gently insert the thermometer's probe tip into the opening of the probe cover box. The cover will snap securely onto the probe tip.
5. **Verify Attachment:** Confirm that the probe cover is fully and securely seated on the thermometer's probe tip. It should fit snugly without any gaps.

**Important:** Always use a new, undamaged probe cover for each measurement to ensure accuracy and prevent cross-contamination.

## OPERATING INSTRUCTIONS

---

Once a probe cover is correctly installed, the thermometer is ready for use. Refer to your specific Welch Allyn Thermoscan ear thermometer's instruction manual for detailed operating procedures. The probe cover facilitates hygienic temperature measurement by providing a barrier between the patient and the thermometer probe.

- Ensure the probe cover is clean and free from damage before use.
- After taking a temperature reading, remove and discard the used probe cover immediately.
- Do not reuse probe covers.

## MAINTENANCE AND STORAGE

---

### Disposal of Used Probe Covers:

After each use, remove the probe cover from the thermometer and dispose of it in a waste receptacle. Probe covers are single-use items and should not be cleaned or sterilized for reuse.

### Storage of New Probe Covers:

Store new, unused probe covers in their original dispenser box in a cool, dry place, away from direct sunlight and extreme temperatures. Keep them out of reach of children.

## TROUBLESHOOTING

---

- **Probe Cover Does Not Fit:** Ensure you are using genuine Welch Allyn Thermoscan PRO series probe covers. These covers are specifically designed for Welch Allyn Thermoscan PRO 6000, 4000, and 4250 series thermometers. If the cover still does not fit, check for any damage to the cover or the thermometer probe.
- **Inaccurate Readings:** If you suspect inaccurate readings, first ensure a new, undamaged probe cover is correctly installed. Also, verify that the thermometer itself is functioning correctly according to its own manual.
- **Damaged Covers:** Do not use probe covers that appear torn, punctured, or otherwise damaged. Discard them and use a new one.

## SPECIFICATIONS

---

<b>Product Name</b>	Welch Allyn Thermoscan PRO 6000 Ear Thermometer Probe Covers
---------------------	--

<b>Model Number</b>	AE06000-800-200
<b>UPC</b>	682318776206
<b>Quantity per Box</b>	200 probe covers (10 x 20-count smaller boxes)
<b>Compatibility</b>	Welch Allyn Thermoscan PRO 6000, PRO 4000, PRO 4250 series ear thermometers
<b>Material</b>	Plastic
<b>Color</b>	Clear
<b>Usage</b>	Single-use, disposable
<b>Manufacturer</b>	A+Elite Medical

## WARRANTY AND SUPPORT

Specific warranty information for these probe covers is not available in the provided product data. For details regarding warranty, returns, or technical support, please contact the manufacturer, A+Elite Medical, or the authorized distributor from whom the product was purchased.

For support related to your Welch Allyn Thermoscan ear thermometer, please refer to the thermometer's specific user manual or contact Welch Allyn customer service.