

## MEGA MOTION NM101

# MEGA MOTION Power Recliner Lift Chair NM101 User Manual

Model: NM101

## IMPORTANT SAFETY INFORMATION

---

Please read and understand all instructions and safety warnings before operating your MEGA MOTION Power Recliner Lift Chair. Failure to follow these instructions may result in injury or damage to the product.

- Always read and follow the safety instructions provided with the chair.
- Do not exceed the maximum weight capacity of 325 lbs (147 kg).
- Be cautious when using the lift mechanism to avoid injury. Ensure the area around the chair is clear before operating.
- Keep the chair away from heat sources and liquids to prevent damage to the electrical components.
- Store the chair in a dry place when not in use.
- In case of any malfunction, immediately discontinue use and contact the manufacturer.
- This chair is intended for indoor use only.

## SETUP INSTRUCTIONS

---

Your MEGA MOTION Power Recliner Lift Chair requires some assembly before first use. Follow these steps carefully:

1. **Unpack Components:** Carefully remove all components from the packaging. Ensure all parts are present: chair base, backrest, power supply, and hand control.
2. **Attach Backrest:** Align the metal brackets on the backrest with the corresponding receivers on the chair base. Slide the backrest down firmly until it locks into place. You should hear a click.
3. **Connect Power:** Locate the power cord from the chair and connect it to the power supply unit. Then, plug the power supply unit into a standard wall outlet.
4. **Connect Hand Control:** Plug the hand control cable into the designated port on the chair's power system.
5. **Test Functionality:** Before sitting, briefly test the recline and lift functions using the hand control to ensure proper operation.



Image: Side view of the recliner, illustrating the internal mechanism and potential assembly points for the backrest.

## OPERATING INSTRUCTIONS

---

The MEGA MOTION Power Recliner Lift Chair offers easy operation with its intuitive hand control.

### Hand Control Functions:

- **Recline Button:** Press and hold to smoothly recline the chair back and raise the footrest. Release the button at your desired position. The chair supports infinite reclining positions, including a lay-flat option.
- **Lift Button:** Press and hold to gently raise the chair, assisting you in standing up from a seated position. Release the button when you reach the desired height or standing position.
- **Dual Motor Control:** This chair features dual motors, allowing independent control of the footrest and backrest for personalized comfort. Consult your specific hand control for dedicated buttons for back and footrest adjustment.
- **USB Port:** A USB port is integrated into the hand wand, allowing you to charge your tablet or phone conveniently while relaxing in the chair.



Image: The MEGA MOTION Power Recliner Lift Chair shown in its upright, seated position, ready for use.



Image: The MEGA MOTION Power Recliner Lift Chair fully reclined, demonstrating the comfortable resting position.

### Emergency Battery Backup:

The chair includes an integrated emergency battery backup system. In the event of a power outage, this system allows for a limited number of operations to return the chair to a seated position. Ensure batteries are installed and checked periodically.

## MAINTENANCE

---

Proper maintenance ensures the longevity and optimal performance of your recliner.

- **Cleaning:** Wipe clean with a damp cloth. Avoid harsh chemicals or abrasive cleaners that may damage the textile material.
- **Power Cord:** Regularly inspect the power cord and connections for any signs of damage. Do not use the chair if the cord is frayed or damaged.
- **Moving Parts:** The steel lift mechanism is designed for durability. Avoid obstructing the moving parts during operation.
- **Spills:** Immediately clean any spills to prevent staining or damage to electrical components.

## TROUBLESHOOTING

---

If you encounter issues with your recliner, refer to the following common solutions:

Problem	Possible Cause	Solution
Chair does not operate	No power supply	Ensure the power cord is securely plugged into a working outlet and the power supply unit. Check circuit breaker.
Hand control unresponsive	Loose connection or faulty control	Check if the hand control cable is firmly connected to the chair. If issues persist, contact support.
Chair moves slowly or erratically	Overload or obstruction	Ensure the chair is not overloaded (max 325 lbs). Check for any obstructions under or around the chair mechanism.
Emergency backup not working	Dead or missing batteries	Replace batteries in the emergency backup system. Ensure correct polarity.

If these steps do not resolve the issue, please contact MEGA MOTION customer support.

## SPECIFICATIONS

Feature	Detail
Brand	MEGA MOTION
Model Name	NM101
Color	Navy
Material	Textile (Foam fill)
Product Dimensions (D x W x H)	36"D x 33.5"W x 43.5"H (91.4 cm D x 85.1 cm W x 110.5 cm H)
Seat Depth	22 inches (55.9 cm)
Item Weight	140 Pounds (63.5 kg)
Maximum Weight Recommendation	325 Pounds (147 kg)
Reclining Position Count	3 (Infinite positions possible)
Frame Material	Furniture grade hardwood/plywood construction
Assembly Required	Yes

## WARRANTY INFORMATION

Your MEGA MOTION Power Recliner Lift Chair is covered by the following warranty:

- **Parts and Service:** One (1) year from the date of purchase.
- **Steel Lift Frame:** Limited lifetime warranty.

Please retain your proof of purchase for warranty claims. The warranty covers manufacturing defects under normal use. It does not cover damage due to misuse, accidents, or unauthorized modifications.

## CUSTOMER SUPPORT

---

For technical assistance, warranty claims, or any questions regarding your MEGA MOTION Power Recliner Lift Chair, please contact the manufacturer directly.

Refer to the contact information provided with your purchase documentation or visit the official MEGA MOTION website for support details.