

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

Manuals.plus /

› [Dr. Scholl's](#) /

› Dr. Scholl's Work All-Day Superior Comfort Insoles Instruction Manual

Dr. Scholl's Work All-Day Superior Comfort Insoles

Dr. Scholl's Work All-Day Superior Comfort Insoles Instruction Manual

Model: Work All-Day Superior Comfort Insoles (Men's 8-14)

1. PRODUCT OVERVIEW

The Dr. Scholl's Work All-Day Superior Comfort Insoles are designed to provide lasting comfort and support for individuals with demanding jobs and active lifestyles that involve prolonged standing or walking on hard surfaces.

Key features include:

- **Clinically Proven All-Day Comfort:** Engineered to deliver continuous comfort throughout the workday.
- **Fatigue Reduction:** Helps to evenly distribute foot pressure, reducing fatigue, joint strain, and muscle soreness.
- **Massaging Gel Technology:** Features extended gel cushioning to reduce impact and shock on muscles and joints, particularly on hard surfaces like asphalt, concrete, and wood floors.
- **Responsive Cushioning:** Adapts to movement, providing energy return with each step for enhanced comfort and sustained productivity.
- **Advanced Odor Control:** Incorporates Polygiene StayFresh Technology to manage foot odor and minimize moisture buildup.



Image 1.1: Dr. Scholl's Work All-Day Superior Comfort Insoles, showing the product and its packaging.

2. INSTALLATION AND TRIMMING

Proper installation ensures optimal comfort and fit. Follow these steps to customize your insoles:

1. **Remove Existing Insoles:** Most shoes come with removable insoles. For best results, remove your shoes' current insoles first.
2. **Test for Size:** Place the Dr. Scholl's insole into your shoe to check the fit.
3. **Trim if Necessary:** If the insole is too large, use sharp scissors to trim along the outline located on the bottom of the insole near the toes. Match this outline to your shoe size. Avoid tracing your foot directly, as this may result in an insole that is too small.
4. **Insert Insole:** Place the trimmed insole into your shoe with the gel side facing down.

TRIM TO FIT MEN'S SHOE SIZES 8-14



Image 2.1: Visual guide for trimming the insoles to achieve a custom fit for men's shoe sizes 8-14.

Specially designed for men with **HARD- WORKING FEET**



Image 2.2: An individual demonstrating the insertion of the insole into a work boot.

3. USAGE AND BENEFITS

Once installed, the insoles work to provide continuous comfort and support throughout your day, especially when on your feet for extended periods or on hard surfaces.

- **Enhanced Comfort:** The massaging gel and responsive cushioning adapt to your foot's movement, providing comfort and reducing impact.
- **Reduced Fatigue:** Designed to alleviate tired, aching feet and legs by distributing pressure evenly.
- **Odor Control:** Polygiene StayFresh Technology helps maintain freshness and reduce moisture.

KEEPS YOU FEELING MORE ENERGIZED

even during
the longest
work day



Image 3.1: Illustration of how the insoles contribute to sustained energy during long workdays.

SUPPORTIVE ARCH

Helps to evenly distribute weight



EXTENDED MASSAGING GEL®

Reduces the impact to muscles & joints on hard surfaces

RESPONSIVE-WAVE CUSHIONING™

Provides added comfort & energy return to keep you productive

Image 3.2: Detailed view of the insole's construction, showing the supportive arch, extended massaging gel, and responsive wave cushioning.

MADE FOR HARD SURFACES LIKE

Asphalt, Concrete, Cement, Hard Dirt and Wood Floors



Image 3.3: The insole's effectiveness on hard surfaces such as asphalt, concrete, cement, hard dirt, and wood floors.

4. MAINTENANCE AND CARE

To maintain the effectiveness and hygiene of your Dr. Scholl's Work All-Day Superior Comfort Insoles:

- **Regular Replacement:** It is recommended to replace the insoles every six months or at the first signs of wear. This ensures continued cushioning and support.
- **Cleaning:** Wipe the insoles with a damp cloth and mild soap if needed. Allow them to air dry completely before reinserting into shoes. Do not machine wash or tumble dry.



Image 4.1: The Polygiene StayFresh Technology logo, indicating the insole's odor control feature.

5. TROUBLESHOOTING

This product is designed for straightforward use. If you experience discomfort, ensure the insoles are correctly sized and positioned within your footwear. If issues persist, consult a foot care professional.

6. PRODUCT SPECIFICATIONS

Feature	Detail
Product Name	Dr. Scholl's Work All-Day Superior Comfort Insoles with Massaging Gel
Model	Work All-Day Superior Comfort Insoles (Men's 8-14)
Brand	Dr. Scholl's
Package Dimensions	13.58 x 4.09 x 2.68 inches
Item Weight	13.76 ounces

ASIN	B01M7VGCTI
Manufacturer	Dr. Scholl's
Origin	Imported

7. WARRANTY AND SUPPORT

Dr. Scholl's offers a 90-day money-back guarantee on this product. For further assistance or to redeem the guarantee, please visit [drscholls.com/guarantee](https://www.drscholls.com/guarantee) or contact customer support at 1-800-986-8022.