



[Manuals.plus](#) /

› [MOVINGLIFE](#) /

› ATTO Folding Travel Powered Mobility Scooter User Manual

MOVINGLIFE AT01-100-B2-0

ATTO Folding Travel Powered Mobility Scooter User Manual

Model: AT01-100-B2-0 | Brand: MOVINGLIFE

INTRODUCTION

The MOVINGLIFE ATTO Folding Travel Powered Mobility Scooter is designed for individuals seeking enhanced mobility and independence. This full-size electric scooter combines advanced engineering with user-friendly features, allowing for seamless travel and daily activities. Its innovative folding mechanism transforms the scooter into a compact, portable unit, making it ideal for transport and storage.



Figure 1: The ATTO Folding Travel Powered Mobility Scooter in its fully assembled, ready-to-use configuration. This image highlights its sleek design and compact footprint.

Key Features:

- **Airline Approved:** Equipped with a lightweight lithium battery that meets FAA regulations for safe air travel. The scooter can be separated into two smaller pieces for easy storage in overhead bins.
- **Ultra-Portable:** Folds into a compact, luggage-like trolley for convenient transport and storage.
- **Enhanced Mobility:** Designed to navigate various environments, from shopping centers to public transportation, offering users greater freedom.
- **User-Centric Design:** Features an adjustable seat, tiller, and multiple speed adjustments for a comfortable and personalized experience.

SETUP AND ASSEMBLY

The ATTO scooter is designed for quick and easy setup. Follow these steps to unfold and prepare your scooter for use.

Unfolding the Scooter:

1. Place the folded scooter on a flat, stable surface.
2. Release the locking mechanisms that secure the scooter in its folded position.
3. Gently extend the front and rear sections until they click into place, forming the full scooter frame.
4. Raise the tiller (handlebar assembly) to your desired height and secure it with the locking pin.
5. Adjust the seat height for optimal comfort and lock it in position.



ATTO Drive

Will take you anywhere
you want to go



ATTO Trolley

Folds easily and
wheels like a trolley



ATTO Compact

Opens to full size
in 10 seconds



ATTO Split

Splits into 2 parts,
for easy lifting and storing

Figure 2: Visual representation of the ATTO's versatile modes: Drive, Trolley, Compact, and Split. The scooter can be easily split into two lighter parts for transport or storage.

Splitting the Scooter for Transport:

For easier lifting and storage, the ATTO can be separated into two parts:

1. Ensure the scooter is in its compact folded state.
2. Locate the release lever or button (usually red) near the central hinge.
3. Press the lever and gently pull the two sections apart.
4. Each section can then be lifted and stored independently, such as in an overhead bin or car trunk.



AIRLINE COMPATIBLE
TSA approved for carry-on
lithium Ion battery



TRAVEL EASY
The ATTO folds and unfolds with
ease and easily fits in the back
seat or trunk of your car



GREAT FOR CRUISES
Wheel around the ship and on
the deck, then store easily when
not in use

Figure 3: The ATTO scooter fitting compactly into a car trunk, demonstrating its travel-friendly design.

OPERATING INSTRUCTIONS

Familiarize yourself with the controls and features of your ATTO scooter for safe and efficient operation.

Controls Overview:



Figure 4: Detailed view of the ATTO scooter's control panel, including the horn, forward/reverse controls, battery power indicator, and speed adjustments.

- **Power Button:** Located on the control panel, press to turn the scooter on/off.
- **Forward/Reverse Controls:** Use the levers on the handlebars to move forward or backward.
- **Speed Adjustments:** The scooter features 6 speed settings, allowing you to adjust your pace up to 4 mph.
- **Battery Power Indicator:** An LED display shows the current battery level.
- **Horn:** Press the horn button to alert others.

Battery and Charging:

The ATTO is powered by a lightweight lithium-ion battery. Ensure the battery is fully charged before each use.

- Connect the charger to the charging port on the scooter.
- Plug the charger into a standard electrical outlet.
- The battery indicator will show charging status. A full charge provides approximately 12 miles or 4.5 hours of travel.

Driving Tips:

- Always maintain a comfortable and stable posture while driving.
- Be aware of your surroundings and potential obstacles.
- Use appropriate speed settings for different terrains and environments.
- The scooter features anti-flat airless wheels for a smooth ride.

Your browser does not support the video tag.

Video 1: Official product video demonstrating the ATTO Mobility Scooter's features, including its folding mechanism and ease of use in various settings. This video is provided by the seller, MOVINGLIFE.

MAINTENANCE

Regular maintenance ensures the longevity and optimal performance of your ATTO scooter.

- **Cleaning:** Wipe down the scooter with a damp cloth. Avoid using abrasive cleaners or excessive water.
- **Battery Care:** Store the battery in a cool, dry place. Charge it regularly, even if not in frequent use, to maintain battery health.
- **Inspections:** Periodically check all nuts, bolts, and moving parts for tightness and wear. Ensure the anti-flat airless wheels are free from debris.
- **Storage:** When storing for extended periods, ensure the scooter is clean and the battery is partially charged (around 50-70%).

TROUBLESHOOTING

If you encounter any issues with your ATTO scooter, refer to the common troubleshooting tips below:

Problem	Possible Cause	Solution
Scooter does not power on	Low battery; Power button not pressed correctly; Loose battery connection.	Charge the battery fully; Press and hold the power button; Ensure battery is securely connected.
Scooter moves slowly or inconsistently	Low battery; Incorrect speed setting; Overloaded.	Charge the battery; Adjust speed setting; Reduce load on scooter.
Folding mechanism is stiff	Debris in hinge; Lack of lubrication.	Clear any debris; Apply a small amount of silicone lubricant to hinges.

For issues not listed here or if problems persist, please contact MOVINGLIFE customer support.

SPECIFICATIONS

Detailed technical specifications for the ATTO Folding Travel Powered Mobility Scooter:

Specification	Value
Model Number	AT01-100-B2-0
Product Dimensions (Unfolded)	15.3"D x 16.5"W x 28"H

Specification	Value
Item Weight	62.2 Pounds
Battery Type	Lithium Ion (FAA Approved)
Max Speed	Up to 4 mph
Range per Charge	Approx. 12 miles / 4.5 hours
Material	Aluminum
Color	White/Chrome
Manufacturer	MOVINGLIFE

WARRANTY AND SUPPORT

MOVINGLIFE is committed to providing high-quality products and excellent customer service.

- **Warranty Information:** For specific warranty details, please refer to the documentation included with your product or visit the official MOVINGLIFE website.
- **Customer Support:** If you have any questions, require technical assistance, or need to report an issue, please contact MOVINGLIFE customer support. Contact information can typically be found on the manufacturer's website or on the product packaging.