

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

- › [Accu-Chek](#) /
- › [Accu-Chek FastClix Lancets \(Model 351-2795\) User Manual](#)

## Accu-Chek 351-2795

# Accu-Chek FastClix Lancets User Manual

Model: 351-2795

## INTRODUCTION

---

This manual provides essential instructions for the safe and effective use of Accu-Chek FastClix Lancets. These lancets are designed exclusively for use with the Accu-Chek FastClix lancing device to facilitate blood glucose testing. Please read this manual thoroughly before using the product.

## PRODUCT OVERVIEW

---



This image displays the retail packaging for Accu-Chek FastClix Lancets. The box clearly indicates 'ACCU-CHEK FastClix 204 Lancets' and highlights features such as 'CLIX MOTION', '6 Lancets in a drum', 'Virtually pain free', 'No lancet handling', and 'For single use only'. It also states compatibility for 'all models of Accu-Chek FastClix' and includes the Roche logo.

Accu-Chek FastClix Lancets are sterile, single-use lancets featuring a unique drum system. Each drum contains 6 preloaded lancets, eliminating the need for individual lancet handling. The 30-gauge needles, combined with Clixmotion technology, are designed for gentle and virtually pain-free blood sampling.

- **Quantity:** 204 lancets (34 drums of 6 lancets each)
- **Compatibility:** Exclusively for Accu-Chek FastClix lancing device
- **Gauge:** 30-gauge for fine, gentle sampling
- **Technology:** Clixmotion for smooth lancing
- **Handling:** No individual lancet handling required

## SETUP AND PREPARATION

Follow these steps to prepare your Accu-Chek FastClix lancing device with a new lancet drum:

1. **Remove the Cap:** Hold the lancing device firmly and pull off the cap.
2. **Insert Lancet Drum:** Insert a new lancet drum, white end first, into the lancing device. Push it in until it clicks into place. Do not force it.
3. **Replace the Cap:** Place the cap back onto the lancing device until it clicks securely.
4. **Set Penetration Depth:** Rotate the cap to select your desired penetration depth. Start with a shallower depth (e.g., 2 or 3) and adjust as needed for sufficient blood sample.

***Note:** The lancet drum is designed to be used only once. Once inserted, it cannot be removed and reused. The device automatically advances to a new lancet after each use.*

## OPERATING INSTRUCTIONS

---

To obtain a blood sample using the Accu-Chek FastClix lancing device with FastClix Lancets:

1. **Prime the Device:** Press the priming button (usually located on the side or end of the device) until it clicks. This action primes the device and advances a new lancet.
2. **Position on Skin:** Place the lancing device firmly against the side of your fingertip.
3. **Lance:** Press the release button to lance your finger.
4. **Obtain Blood Sample:** Gently squeeze your finger to form a drop of blood. Avoid excessive squeezing, which can dilute the sample.
5. **Test Blood Glucose:** Apply the blood drop to your blood glucose meter's test strip as instructed by your meter's manual.

## CARE AND MAINTENANCE

---

- **Storage:** Store lancet drums in their original packaging in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Disposal:** Used lancet drums contain sharp, contaminated lancets. Dispose of them safely in a puncture-resistant sharps container immediately after all 6 lancets have been used. Do not dispose of used lancets or drums in household trash. Follow local regulations for sharps disposal.
- **Cleaning:** The lancets themselves do not require cleaning as they are single-use. Refer to your Accu-Chek FastClix lancing device manual for cleaning instructions for the device.

## TROUBLESHOOTING

---

- **No Blood Sample:** If you do not get enough blood, try increasing the penetration depth setting on your lancing device. Ensure the device is held firmly against the skin.
- **Lancet Drum Not Inserting:** Ensure the drum is oriented correctly (white end first) and pushed in until it clicks. Do not force it. If it still doesn't insert, check for any obstructions or damage to the drum or device.
- **Device Not Priming/Firing:** Ensure the cap is securely attached. If the device has been used 6 times, the drum is empty and needs to be replaced. If a new drum is inserted and the device still doesn't prime, refer to your Accu-Chek FastClix lancing device manual.

## PRODUCT SPECIFICATIONS

---

- **Product Name:** Accu-Chek FastClix Lancets

- **Model Number:** 351-2795
- **Lancet Gauge:** 30-gauge
- **Lancets per Drum:** 6
- **Total Lancets:** 204 (34 drums)
- **Sterility:** Sterile, for single use only
- **Dimensions:** 1.38 x 3.15 x 4.72 inches (packaging)
- **Weight:** 2.12 ounces (packaging)
- **Manufacturer:** Abbott Diabetes Care



## IMPORTANT SAFETY INFORMATION





- **Single Use Only:** Lancets are sterile and intended for single use only. Reusing lancets can lead to infection and dulling of the needle, causing more pain.
- **Personal Use:** This product is for individual use only. Do not share lancets or lancing devices with others, as this poses a risk of transmitting infections.
- **Keep Out of Reach of Children:** Store lancets and lancing devices safely away from children.
- **Proper Disposal:** Always dispose of used lancet drums in an approved sharps container.
- **Consult Healthcare Professional:** Always consult your healthcare professional for guidance on blood glucose testing and proper use of medical devices.

**Legal Disclaimer:** Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is a medical device and should be used as directed.

© 2023 Accu-Chek. All rights reserved.

## Related Documents - 351-2795

	<p><a href="#">Accu-Chek Instant Blood Glucose Meter User Guide and Support</a></p> <p>Comprehensive guide to using the Accu-Chek Instant blood glucose meter, including setup, testing, connecting with the mySugr app, and customer support information.</p>
	<p><a href="#">Accu-Chek Guide Me Quick Start Guide: Setup and Blood Glucose Testing</a></p> <p>Concise guide to setting up your Accu-Chek Guide Me blood glucose meter and performing blood glucose tests, including lancing device setup and meter operation.</p>

	<p><a href="#">Accu-Chek Instant Quick Reference Guide</a></p> <p>A quick reference guide for the Accu-Chek Instant blood glucose monitoring system, detailing how to perform a blood glucose test and set up the lancing device.</p>
	<p><a href="#">Accu-Chek Instant Quick Reference Guide: Blood Glucose Monitoring</a></p> <p>A quick reference guide for the Accu-Chek Instant blood glucose meter. Learn how to perform blood glucose tests, set up the device, change batteries, and use the Accu-Chek FastClix lancing device for diabetes management.</p>
	<p><a href="#">Accu-Chek Guide User's Manual: Blood Glucose Monitoring</a></p> <p>Explore the Accu-Chek Guide Blood Glucose Monitoring System with its user-friendly manual. Learn about meter setup, testing procedures, data management, and maintenance for effective diabetes care.</p>
	<p><a href="#">Accu-Chek Guide Blood Glucose Monitoring System: User Guide and Test Strip Information</a></p> <p>Learn how to use the Accu-Chek Guide, Guide Link, and Guide Me blood glucose monitoring systems and test strips for effective diabetes management. Includes instructions, safety information, and performance data.</p>