

MOPEDIA RS826

MOPEDIA Foldable Shower Chair with Backrest RS826 Instruction Manual

1. INTRODUCTION

Thank you for choosing the MOPEDIA RS826 Foldable Shower Chair with Backrest. This manual provides essential information for the safe and effective use of your new shower chair. Please read these instructions carefully before assembly and use, and retain them for future reference.

The MOPEDIA RS826 is designed to provide support and stability in the shower, enhancing safety and comfort for individuals who require assistance. It features a robust yet lightweight painted steel frame, a durable polyethylene seat and backrest, and anti-slip rubber tips on the legs. Its foldable design allows for convenient storage, and the height is adjustable to suit individual user needs.

2. SAFETY INFORMATION AND WARNINGS

Your safety is paramount. Please observe the following warnings and safety guidelines:

- **Read All Instructions:** Always read the label, warnings, and instructions provided with the product before use. The product packaging may contain information different from that shown on our website.
- **Weight Capacity:** Do not exceed the maximum weight capacity of **120 kg (264 lbs)**.
- **Stability Check:** Before each use, ensure the chair is fully unfolded, all legs are securely locked in place, and the anti-slip rubber tips are firmly on the shower floor.
- **Wet Surfaces:** Exercise extreme caution on wet or soapy surfaces. Ensure the shower floor is as dry as possible before entering and exiting the shower.
- **Regular Inspection:** Periodically inspect the chair for any signs of wear, damage, or loose components. Do not use the chair if it is damaged. Contact MOPEDIA support for replacement parts or service.
- **Proper Use:** Use the chair only for its intended purpose as a shower aid. Do not stand on the chair.
- **Children and Pets:** Keep children and pets away from the chair during assembly and use.

3. SETUP AND ASSEMBLY

The MOPEDIA RS826 shower chair is designed for easy setup. It arrives largely pre-assembled and requires unfolding and height adjustment.

1. **Unfolding the Chair:** Carefully remove the chair from its packaging. Gently unfold the chair until the seat and backrest are in their upright position and the frame is fully extended. Ensure all hinges and locking mechanisms (if present) are engaged.
2. **Adjusting Height:** Refer to Section 4, "Height Adjustment," to set the chair to a comfortable and safe height for the user. Ensure all four legs are adjusted to the same height.
3. **Final Check:** Before first use, firmly press down on the seat to ensure the chair is stable and all legs are securely locked at the desired height.



Figure 1: MOPEDIA RS826 Foldable Shower Chair with Backrest in its fully assembled and upright position. The chair features a white polyethylene seat and backrest, with a painted steel frame and adjustable legs equipped with anti-slip rubber tips.

4. HEIGHT ADJUSTMENT

The MOPEDIA RS826 shower chair offers adjustable height settings to accommodate various user needs. The height can be adjusted in 2.5 cm (approximately 1 inch) increments.

1. **Locate Adjustment Pins:** On each leg, locate the spring-loaded adjustment pin.
2. **Press and Slide:** Press the adjustment pin inward and slide the lower leg section up or down to the desired height.
3. **Lock in Place:** Release the pin, ensuring it fully protrudes through one of the adjustment holes, locking the leg securely in place.
4. **Even Adjustment:** Repeat this process for all four legs, ensuring they are all set to the exact same height to maintain stability.
5. **Verify Stability:** After adjustment, always check the chair's stability by gently rocking it and applying downward pressure before use.

5. OPERATING INSTRUCTIONS (USAGE)

To use the MOPEDIA RS826 shower chair safely:

- **Placement:** Position the chair securely in the shower area, ensuring all four anti-slip rubber tips are firmly on the shower floor. Avoid placing it on uneven surfaces or over drains.
- **Sitting:** Carefully approach the chair. Use any available grab bars for support if needed. Turn around and slowly lower yourself onto the seat, ensuring your weight is evenly distributed.
- **During Shower:** Maintain a stable posture. The backrest and armrests provide additional support.
- **Standing Up:** When finished, use the armrests for support to push yourself up slowly. If possible, use grab bars for added stability.
- **After Use:** If desired, the chair can be folded for storage (refer to Section 7).

6. MAINTENANCE AND CLEANING

Regular maintenance and cleaning will prolong the life of your shower chair and ensure hygienic use.

- **Cleaning:** The polyethylene seat and backrest are designed for easy cleaning. Use a mild soap or non-abrasive household cleaner and warm water. Wipe down all surfaces with a soft cloth.
- **Rinsing:** Rinse thoroughly with clean water to remove all soap residue.
- **Drying:** Dry the chair completely with a clean cloth after each cleaning to prevent water spots and maintain hygiene.
- **Rust Resistance:** The painted steel frame and polyethylene components are rust-resistant, but proper drying helps maintain their condition.
- **Inspection:** Periodically check the anti-slip rubber tips for wear. Replace them if they show signs of cracking or significant wear to ensure continued safety.

7. STORAGE

The MOPEDIA RS826 shower chair is foldable for convenient storage when not in use.

1. **Fold Down:** Ensure the chair is clean and dry. Gently fold the chair by reversing the unfolding process.
2. **Secure:** Ensure any folding mechanisms are securely engaged to keep the chair in its folded position.
3. **Location:** Store the folded chair in a dry, cool place away from direct sunlight and extreme temperatures.

8. TROUBLESHOOTING

If you encounter any issues with your MOPEDIA RS826 shower chair, please refer to the following common solutions:

- **Chair Feels Unstable:**
 - Ensure all four legs are adjusted to the same height.
 - Verify that the adjustment pins on each leg are fully engaged and locked into a hole.
 - Check that the anti-slip rubber tips are clean and free of soap residue or debris.
 - Ensure the chair is placed on a flat, even surface in the shower.
- **Difficulty Folding/Unfolding:**
 - Check for any obstructions in the folding mechanism.
 - Ensure no parts are bent or damaged. Do not force the chair if it resists folding or unfolding.

If these steps do not resolve the issue, please contact MOPEDIA customer support.

9. SPECIFICATIONS

- **Model:** RS826
- **Overall Dimensions (L x W x H):** 52 cm x 52 cm x 90 - 100 cm (Adjustable in 2.5 cm increments)
- **Seat Dimensions (L x W):** 40 cm x 40 cm
- **Armrest Height from Seat:** 19 cm
- **Seat Height Adjustment Range:** 49 cm - 59 cm (Adjustable in 2.5 cm increments)
- **Maximum Weight Capacity:** 120 kg (264 lbs)

- **Product Weight:** 4 kg (8.8 lbs)
- **Frame Material:** Painted Steel Alloy
- **Seat and Backrest Material:** Polyethylene
- **Leg Tips Material:** Rubber
- **Features:** Foldable, Lightweight, Backrest, Armrests, Adjustable Height, Anti-slip Feet

10. WARRANTY AND SUPPORT

MOPEDIA is committed to providing high-quality products. For information regarding warranty coverage, technical support, or to purchase replacement parts, please contact MOPEDIA customer service through your retailer or the official MOPEDIA website.

Please have your model number (RS826) and purchase date available when contacting support.

