

MOPEDIA RS838

MOPEDIA Onda Series RS838 Backless Shower Seat with Padded Armrests User Manual

1. INTRODUCTION

Thank you for choosing the MOPEDIA Onda Series RS838 Backless Shower Seat. This manual provides essential information for the safe and effective use, assembly, maintenance, and troubleshooting of your new shower seat. Please read these instructions carefully before initial use and retain them for future reference.

The MOPEDIA Onda Series RS838 is designed to provide a stable and comfortable seating solution in the shower, enhancing safety and independence for users. It features a durable aluminum frame, a comfortable polyethylene seat, padded armrests, and height-adjustable legs with suction-cup rubber tips for secure placement.

2. SAFETY INFORMATION

Important: The product packaging may contain information that differs from what is presented on our website. Always read the label, warnings, and instructions provided with the product before use.

- **Weight Limit:** Do not exceed the maximum weight capacity specified for this shower seat. Refer to the specifications section for details.
- **Stability Check:** Before each use, ensure all legs are firmly on the shower floor and the rubber tips provide adequate grip. The seat must be stable and not wobble.
- **Assembly:** Ensure all components are correctly assembled and securely fastened before use.
- **Surface:** Use the shower seat only on a flat, non-slip surface.
- **Maintenance:** Regularly inspect the seat for any signs of wear, damage, or loose components. Discontinue use immediately if any damage is observed.
- **Cleaning:** Clean the seat regularly to prevent soap scum buildup, which can reduce the effectiveness of the rubber tips.

- **Children:** Keep children away from the shower seat during assembly and ensure they do not play with it unsupervised.

3. PRODUCT COMPONENTS

The MOPEDIA Onda Series RS838 shower seat consists of the following main components:

- Polyethylene Seat (White)
- Aluminum Frame (Legs and Crossbars)
- Padded Armrests (Grey)
- Height Adjustment Pins
- Suction-Cup Rubber Tips (for legs)



This image displays the MOPEDIA Onda Series RS838 backless shower seat. It features a white, perforated seat, an aluminum frame, and grey padded armrests. The legs are adjustable in height and have grey rubber tips for stability.

4. SETUP AND ASSEMBLY

The MOPEDIA Onda Series RS838 shower seat is designed for tool-free assembly. Follow these steps to set up your shower seat:

1. **Unpack Components:** Carefully remove all parts from the packaging and ensure all listed components are present.
2. **Attach Legs:** Insert each leg into the designated slots on the underside of the seat frame. Ensure the spring-loaded pins snap securely into place.
3. **Attach Armrests:** Slide the padded armrests onto the side supports of the seat frame until they click into position. Ensure they are firmly attached.
4. **Adjust Height:** Press the spring-loaded pin on each leg and slide the inner leg tube up or down to the desired height. Ensure all four legs are set to the same height for stability. Release the pin and ensure it locks into one of the adjustment holes.
5. **Final Check:** Before use, gently push down on the seat to ensure it is stable and all connections are secure.

5. OPERATING INSTRUCTIONS

Using your MOPEDIA Onda Series RS838 shower seat correctly ensures safety and comfort.

- **Placement:** Place the shower seat in the desired position within your shower or bath area. Ensure all four suction-cup rubber tips are in full contact with a clean, flat, and non-slip surface.
- **Stability:** Before sitting, apply gentle pressure to the seat to confirm its stability and that it does not slide or wobble.
- **Sitting and Standing:** Use the padded armrests for support when sitting down or standing up. Distribute your weight evenly on the seat.
- **Height Adjustment:** If the seat height needs to be changed, follow the 'Adjust Height' steps in the Setup and Assembly section. Always adjust all four legs to the same height.

6. MAINTENANCE

Regular maintenance will prolong the life of your shower seat and ensure its safe operation.

- **Cleaning:** Clean the shower seat regularly with mild soap and water. Avoid abrasive cleaners, solvents, or harsh chemicals, as these can damage the materials. Rinse thoroughly and dry completely after cleaning.
- **Inspection:** Periodically inspect all components for signs of wear, cracks, corrosion, or loose connections. Pay close attention to the rubber tips, armrest padding, and leg adjustment pins.
- **Rubber Tips:** Ensure the suction-cup rubber tips are clean and free of debris or soap residue to maintain maximum grip. Replace worn or damaged tips immediately.
- **Storage:** If storing the seat for an extended period, ensure it is clean and dry. Store in a cool, dry place away from direct sunlight.

7. TROUBLESHOOTING

If you encounter issues with your MOPEDIA Onda Series RS838 shower seat, refer to the following common problems and solutions:

- **Seat is Unstable/Wobbles:**
 - Ensure all four legs are adjusted to the same height.
 - Verify that the spring-loaded pins on each leg are fully engaged in the adjustment holes.
 - Check that the rubber tips are clean and making full contact with the shower floor.
 - Ensure the shower floor is flat and even.
- **Difficulty Adjusting Height:**
 - Ensure the spring-loaded pin is fully depressed before attempting to slide the leg.
 - Check for any debris or obstruction in the adjustment holes or inside the leg tubes.
- **Armrests Feel Loose:**
 - Ensure the armrests are fully inserted and clicked into their locking positions on the frame.

If the problem persists or you experience any other issues, please contact MOPEDIA customer support or your retailer.

8. SPECIFICATIONS

Feature	Detail
Brand	MOPEDIA
Model	Onda Series RS838
Material	Aluminum (frame), Polyethylene (seat)
Color	White (seat), Silver (frame), Grey (armrests/tips)
Assembly	Tool-free
Features	Height adjustable legs, Padded armrests, Suction-cup rubber tips
Backrest	No backrest

9. WARRANTY AND SUPPORT

For information regarding the warranty of your MOPEDIA Onda Series RS838 shower seat, please refer to the warranty card included with your product packaging or contact your retailer or the manufacturer directly.

If you require technical assistance, have questions about product use, or need to report a defect, please contact MOPEDIA customer support through the contact information provided on the product packaging or the official MOPEDIA website.