

[Manuals.plus](#) /

> [MOPEDIA](#) /

> MOPEDIA Shower Seat with Armrests and Backrest - Model RS839 User Manual

MOPEDIA RS839

MOPEDIA Shower Seat with Armrests and Backrest - Model RS839 User Manual

INTRODUCTION

This manual provides essential information for the safe and effective use of your MOPEDIA Shower Seat, Model RS839. Please read these instructions carefully before assembly and use, and retain them for future reference. The MOPEDIA Shower Seat is designed to provide support and stability during showering. It features a durable aluminum frame, a comfortable polyethylene seat and backrest, and is height-adjustable to suit individual needs. The non-slip rubber tips ensure secure placement in the shower.



Image: MOPEDIA Shower Seat Model RS839. This image shows the complete shower seat with its aluminum frame, white polyethylene seat, backrest, and armrests, designed for stability and comfort in the shower.

1. SETUP AND ASSEMBLY

The MOPEDIA Shower Seat is designed for tool-free assembly. Follow these steps to set up your shower seat:

1. **Unpack Components:** Carefully remove all parts from the packaging. Verify that all components are present: seat, backrest, armrests, frame legs, and rubber tips.
2. **Attach Legs:** Insert the four frame legs into the designated slots on the underside of the seat. Ensure they click securely into place.
3. **Install Backrest:** Slide the backrest into the slots at the rear of the seat until it locks.
4. **Attach Armrests:** Insert the armrests into the side slots of the seat. Ensure they are firmly attached.
5. **Adjust Height:** Press the push-pins on each leg and slide the leg up or down to achieve the desired height. Ensure all four legs are adjusted to the same height for stability. Release the push-pins to lock the legs in position.
6. **Check Stability:** Before first use, place the shower seat on a flat, stable surface and apply gentle pressure to ensure it is stable and does not wobble.

2. OPERATING INSTRUCTIONS

To use the MOPEDIA Shower Seat safely and effectively:

- **Placement:** Position the shower seat on a flat, non-slip surface inside your shower or bath area. Ensure all four rubber tips are firmly in contact with the floor.
- **Sitting:** Approach the seat carefully. Use the armrests for support when sitting down or standing up. Distribute your weight evenly on the seat.
- **Stability Check:** Always ensure the seat is stable before transferring your full weight onto it.
- **Weight Limit:** Do not exceed the maximum weight capacity specified in the 'Specifications' section.
- **Supervision:** If the user has significant mobility limitations, ensure supervision is available during use.

3. CARE AND MAINTENANCE

Proper care will extend the life of your MOPEDIA Shower Seat:

- **Cleaning:** Clean the seat regularly with a mild soap and water solution. Rinse thoroughly and dry with a soft cloth. Avoid abrasive cleaners or harsh chemicals, as these can damage the materials.
- **Inspection:** Periodically inspect all components for wear, damage, or loose connections. Pay close attention to the rubber tips for signs of cracking or excessive wear.
- **Replacement Parts:** If any part is damaged or worn, contact customer support for replacement parts. Do not use the shower seat if any component is compromised.
- **Storage:** Store the shower seat in a dry place when not in use for extended periods.

4. TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Shower seat wobbles or is unstable.	Legs are not adjusted to the same height or not fully locked. Rubber tips are worn or missing. Uneven shower floor.	Ensure all legs are adjusted to the same height and push-pins are fully engaged. Check rubber tips for wear and replace if necessary. Reposition the seat on a flatter surface.
Difficulty adjusting leg height.	Push-pins are stuck or dirty.	Press the push-pin firmly. Clean any debris around the pin. If still difficult, apply a small amount of silicone lubricant.
Seat or backrest feels loose.	Components not fully engaged during assembly.	Re-check assembly of the seat, backrest, and armrests to ensure they are securely clicked into place.

5. PRODUCT SPECIFICATIONS

- **Model:** RS839
- **Brand:** MOPEDIA
- **Manufacturer:** Moretti S.p.A.
- **Material:** Aluminum frame, Polyethylene seat and backrest
- **Color:** White
- **Features:** Height adjustable, Armrests, Backrest, Suction cup rubber tips
- **Assembly:** Tool-free
- *(Note: Specific dimensions and weight capacity are not provided in the product data. Please refer to product packaging for detailed measurements and maximum weight limit.)*

6. WARRANTY INFORMATION

MOPEDIA products are manufactured to high-quality standards. For specific warranty details regarding your RS839 Shower Seat, please refer to the warranty card included with your purchase or visit the official MOPEDIA website. Keep your proof of purchase for any warranty claims.

7. CUSTOMER SUPPORT

If you have any questions, require assistance with assembly, or need replacement parts, please contact MOPEDIA customer support. Contact information can typically be found on the product packaging or the manufacturer's official website.

Manufacturer: Moretti S.p.A.