

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Graham-Field](#) /
- › [Lumex Bariatric Imperial Offset Cane User Manual](#)

Graham-Field 6334A-1



Lumex Bariatric Imperial Offset Cane

MODEL: 6334A-1

1. Product Overview

The Lumex Bariatric Imperial Offset Cane is an adjustable aluminum walking stick designed to provide mobility assistance for individuals requiring additional support. Its offset handle design helps center the user's weight over the cane shaft for improved stability.





Image: The Lumex Bariatric Imperial Offset Cane, featuring its black finish and offset handle design.

2. Safety Information

- Always consult with a healthcare professional to determine if a cane is appropriate for your mobility needs and to receive proper training on its use.
- Ensure all adjustment mechanisms are securely locked before each use.
- Regularly inspect the cane tip for wear. Replace worn tips immediately to prevent slipping.
- Do not exceed the maximum weight capacity specified for this cane.
- Use caution on uneven or slippery surfaces.

3. Setup and Adjustment

The Lumex Bariatric Imperial Offset Cane is designed for easy height adjustment to ensure proper fit for the user.

Adjusting Cane Height:

1. Loosen the locking collar (if present) near the adjustment pin.
2. Press the push-button pin on the lower shaft of the cane.
3. Slide the lower shaft up or down to the desired height. Ensure the push-button pin fully protrudes through the corresponding hole, indicating it is securely locked.
4. Tighten the locking collar (if present) to further secure the adjustment.

Proper Height: When standing upright with arms relaxed at your sides, the top of the cane handle should align with the crease in your wrist. When holding the cane, your elbow should have a slight bend (approximately 15-20 degrees).

4. Operating Instructions

Using your Lumex Bariatric Imperial Offset Cane correctly is crucial for stability and safety.

Walking with the Cane:

- Hold the cane in the hand opposite to your weaker leg.
- Move the cane forward simultaneously with your weaker leg.
- Step forward with your stronger leg, passing the cane.
- Maintain a steady, natural gait.

Using Stairs:

- **Going Up:** Lead with your stronger leg, then bring the cane and weaker leg up to the same step.
- **Going Down:** Place the cane on the lower step, then bring your weaker leg down, followed by your stronger leg.
- Always use a handrail if available.

5. Maintenance and Care

Regular maintenance will extend the life of your cane and ensure its safe operation.

- **Cleaning:** Wipe the cane with a damp cloth and mild soap. Dry thoroughly. Do not use abrasive cleaners.
- **Inspection:** Periodically check all parts for wear, damage, or looseness. Pay close attention to the cane tip, adjustment pin, and handle.
- **Storage:** Store the cane in a dry place away from extreme temperatures.

6. Troubleshooting

Problem	Possible Cause	Solution
Cane feels unstable or wobbly.	Adjustment pin not fully engaged; locking collar loose.	Ensure adjustment pin is fully engaged and locking collar is tightened.
Cane tip slips easily.	Worn or damaged cane tip.	Replace the cane tip immediately with a new, compatible tip.

Problem	Possible Cause	Solution
Difficulty adjusting height.	Debris in adjustment holes; pin stuck.	Clean adjustment holes. If pin is stuck, apply a small amount of silicone spray and try again. Do not force.

7. Specifications

Model Number: 6334A-1

Brand: Graham-Field (Lumex)

Material: Aluminum

Type: Offset Cane

Package Dimensions: 31.9 x 8.4 x 2.5 inches

Item Weight: 3.5 Pounds

Country of Origin: China

8. Warranty and Support

Graham-Field products are manufactured to the highest standards and are backed by a commitment to quality. For specific warranty information regarding your Lumex Bariatric Imperial Offset Cane, please refer to the documentation included with your purchase or visit the official Graham-Field website.

For technical support, replacement parts, or any inquiries regarding your product, please contact Graham-Field customer service. Contact information can typically be found on the manufacturer's website or on the product packaging.



Image: The Graham-Field family of brands, including Lumex.

