

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [NRS Healthcare](#) /
- › [NRS Healthcare Easy Rise Cushion N98712 User Manual](#)

NRS Healthcare N98712

NRS Healthcare Easy Rise Cushion N98712 User Manual

Model: N98712

1. INTRODUCTION

The NRS Healthcare Easy Rise Cushion N98712 is designed to assist individuals in transitioning from a seated to a standing position with greater ease. This cushion utilizes a specially layered foam construction that provides a gentle lifting action as the user leans forward, effectively increasing the seat height by 10 cm (4 inches). Its lightweight and portable design makes it suitable for use in various locations, offering a practical alternative to permanent furniture risers.



Figure 1: Top-down view of the NRS Healthcare Easy Rise Cushion.

2. SETUP AND PLACEMENT

The Easy Rise Cushion requires minimal setup. Follow these guidelines for optimal use:

1. **Unpacking:** Remove the cushion from its packaging.
2. **Placement:** Position the cushion centrally on the seat of your chosen chair. Ensure the chair is stable and on a level surface. The cushion should sit flat and securely without slipping.
3. **Stability Check:** Before use, gently press down on the cushion to confirm it is stable and does not shift.

excessively on the chair.



Figure 2: The Easy Rise Cushion correctly positioned on an armchair.

3. OPERATING INSTRUCTIONS

The Easy Rise Cushion is designed to assist with standing. Follow these steps for effective operation:

1. **Sitting Down:** Carefully sit onto the center of the cushion. Ensure you are fully supported and balanced.
2. **Preparing to Stand:** When ready to stand, position your feet firmly on the floor.
3. **Leaning Forward:** Lean your upper body forward. As you shift your weight forward, the layered foam within the cushion will compress and then gently push upwards, providing a natural spring action to assist your rise.
4. **Standing Up:** Continue to push up with your legs as the cushion assists, transitioning smoothly to a standing position.

Note: The cushion works by utilizing your body weight and forward momentum. It is not an automatic lifting device but an aid to reduce the effort required to stand.

4. MAINTENANCE AND CARE

To ensure the longevity and hygiene of your NRS Healthcare Easy Rise Cushion, follow these maintenance guidelines:

- **Cleaning:** The cushion features a black leatherette finish, which is easy to clean. Wipe the surface with a damp cloth and mild soap solution. Do not use abrasive cleaners or harsh chemicals.
- **Drying:** Allow the cushion to air dry completely before reuse. Do not machine wash or tumble dry.
- **Storage:** Store the cushion in a cool, dry place away from direct sunlight and extreme temperatures when not in use.
- **Inspection:** Periodically inspect the cushion for any signs of wear, tears, or damage to the cover or foam. Discontinue use if significant damage is observed.



Figure 3: Close-up of the cushion's zipper, indicating removable cover for cleaning.

5. TROUBLESHOOTING

If you encounter issues with your Easy Rise Cushion, consider the following:

- **Cushion feels unstable:** Ensure the cushion is placed centrally and securely on a flat, stable chair seat. Avoid using on chairs with uneven or very soft seating surfaces.
- **Cushion does not provide enough lift:** Ensure you are leaning forward sufficiently when attempting to stand. The layered foam requires a forward weight shift to activate its lifting mechanism. Verify the cushion is not excessively compressed or damaged.

- **Discomfort during use:** Check that the cushion is positioned correctly and that your posture is aligned. If discomfort persists, consult with a healthcare professional to ensure the product is suitable for your specific needs.

6. SPECIFICATIONS

Feature	Detail
Model Number	N98712
Brand	NRS Healthcare
Dimensions (L x W x H)	46 x 46 x 11 cm (approx. 18.1 x 18.1 x 4.3 inches)
Weight	900 g (approx. 1.98 lbs)
Material	Foam core with black leatherette cover
Color	Black
Seat Height Increase	10 cm (4 inches)
Special Features	Lightweight, Portable, Easy to Clean
Recommended Age Range	Adult

7. WARRANTY AND SUPPORT

For information regarding warranty, returns, or technical support, please refer to the purchase documentation or contact NRS Healthcare directly through their official website or customer service channels. Keep your proof of purchase for any warranty claims.