

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Beurer](#) /
- › [Beurer Bluetooth Digital Fingertip Pulse Oximeter PO60 User Manual](#)

Beurer PO60

Beurer Bluetooth Digital Fingertip Pulse Oximeter PO60 - User Manual

Brand: Beurer | Model: PO60

INTRODUCTION

The Beurer Bluetooth Digital Fingertip Pulse Oximeter, Model PO60, is designed for the precise measurement of arterial oxygen saturation (SpO2) and heart rate (pulse). This portable device is suitable for use at home or on the go, offering a convenient way to monitor key health metrics. Its compact design and Bluetooth connectivity allow for easy data transfer to the Beurer HealthCoach app for storage and evaluation of measured values.

WHAT'S INCLUDED

- 1 x Beurer Bluetooth Digital Fingertip Pulse Oximeter (PO60)
- 1 x Lanyard
- 1 x Protective Storage Bag
- 2 x AAA Batteries (included)
- 1 x Instruction Manual (PDF version available via link below)

SETUP GUIDE

1. Battery Installation

The device requires 2 AAA batteries, which are included. To install:

1. Locate the battery compartment cover on the back of the oximeter.
2. Slide the cover open.
3. Insert the 2 AAA batteries, ensuring correct polarity (+/-) as indicated inside the compartment.

4. Close the battery compartment cover securely.



Image: The Beurer PO60 Pulse Oximeter, ready for use after battery installation.

OPERATING INSTRUCTIONS

To obtain accurate measurements:

1. Ensure your finger is clean and free of nail polish or artificial nails.
2. Open the clamp of the pulse oximeter.
3. Insert one of your fingers (preferably the index finger) into the opening until it reaches the stop. Ensure the finger is placed correctly.
4. Release the clamp.
5. Press the power button to turn on the device.
6. Keep your hand still during the measurement. The device will display your SpO2 and pulse rate readings within a few seconds.
7. The device features an automatic switch-off function when no finger is detected.

Determining arterial oxygen saturation (SpO₂) and Heart rate (pulse frequency)

 german
engineering



Image: Proper finger placement for accurate SpO₂ and heart rate measurement.

Display Features

- **Color LED Display:** Provides clear readings of oxygen saturation and pulse.
- **Automatic Screen Rotation:** The display orientation adjusts automatically for easy viewing from any direction.
- **Real-time Pulse Curve Display:** Visual representation of your pulse.

EASY-TO-READ COLOR DISPLAY



 german
engineering

Image: The easy-to-read color display with automatic rotation for various viewing angles.

Bluetooth Connectivity

The PO60 oximeter connects to your smartphone via Bluetooth through the free Beurer HealthCoach app. This allows for wireless data transfer, enabling you to store and evaluate up to 100 oxygen and pulse measurements directly on your device.

- SMART BLUETOOTH -

Transfer measurements to your
smartphone



Image: Seamless data transfer from the oximeter to your smartphone via Bluetooth.

Official Product Video

Your browser does not support the video tag.

Video: This video provides a brief overview of the Beurer Bluetooth Pulse Oximeter PO60, highlighting its smart Bluetooth connectivity for wireless data transfer, color display with four display formats, suitability for individuals with heart failure, chronic obstructive pulmonary diseases, and bronchial asthma, and its ability to save up to 100 measurements. It also shows the device being stored securely in its belt pouch.

KEY FEATURES

- **Accurate Measurements:** German-engineered for accuracy, providing real-time pulse and oxygen data.
- **Portable Design:** Lightweight and easy to use, ideal for monitoring during exercise or at high altitudes.
- **Bluetooth Connectivity:** Seamlessly connects to the Beurer HealthCoach app for data transfer and tracking.
- **Memory Function:** Stores up to 100 oxygen and pulse measurements directly on the device.

- **Adjustable Display:** Bright color LED display with automatic screen rotation for clear viewing.
- **Included Accessories:** Comes with a lanyard and protective storage bag for convenience.



Image: Key features of the Beurer PO60, including blood oxygen saturation, pulse rate measurement, one-button operation, finger opening, and lanyard loop hole.

SMALL AND LIGHT

– for use at home and
on the move



Incl. retaining strap and belt pouch

 **german
engineering**

Image: The compact and lightweight design of the oximeter, suitable for home and travel, shown with its included belt pouch.

SUITABLE FOR SPORTS AT HIGH ALTITUDES

 german
engineering



Image: The oximeter is suitable for various activities, including sports at high altitudes like mountaineering, skiing, and aviation sports.

TECHNICAL SPECIFICATIONS

Specification	Value
Product Dimensions	3 x 2 x 2 inches
Item Weight	5.93 ounces
Model Number	PO60
Batteries	2 AAA batteries required (included)
Measuring Range (SpO2)	70% - 100%

Measuring Range (Pulse Rate)	30 - 250 bpm
Color	Silver
Manufacturer	Beurer North America LP

PRODUCT MEASUREMENTS



Image: Detailed product measurements for the Beurer PO60 Pulse Oximeter.

MAINTENANCE AND CARE

- **Cleaning:** Wipe the device with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the oximeter in its protective storage bag in a cool, dry place away from direct sunlight and extreme temperatures.
- **Battery Replacement:** Replace batteries when the low battery indicator appears on the display. Remove batteries if the device will not be used for an extended period.

TROUBLESHOOTING

If you encounter issues with your Beurer PO60 Pulse Oximeter, please refer to the full user manual for detailed troubleshooting steps. Common issues may include incorrect readings due to improper finger placement, low battery, or movement during measurement. Ensure the device is clean and batteries are correctly installed.

WARRANTY AND SUPPORT

The Beurer PO60 Pulse Oximeter comes with a limited lifetime warranty. For warranty claims, technical support, or any product-related inquiries, please contact Beurer North America LP. Their contact information is:

*1 Oakwood Boulevard, Suite 255,
Hollywood, FL 33020, USA.*

ADDITIONAL RESOURCES

For a comprehensive guide and detailed instructions, please download the official User Manual (PDF):

[Download User Manual
\(PDF\)](#)