

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [McKesson](#) /

› [McKesson Lightweight Aluminum Rollator Walker with Seat and Wheels Instruction Manual](#)

McKesson B0785N3TYX

McKesson Lightweight Aluminum Rollator Walker with Seat and Wheels

Model: B0785N3TYX

Instruction Manual

INTRODUCTION

This manual provides essential information for the safe and effective use of your McKesson Lightweight Aluminum Rollator Walker. Please read all instructions carefully before assembly and operation. This 4-wheel rollator features a durable aluminum frame, an oversized padded seat and backrest, and a convenient storage pouch. It is designed to support users up to 300 pounds and is suitable for both indoor and outdoor use.

SAFETY INFORMATION

- Always ensure all parts are securely fastened before use.
- Do not exceed the maximum weight capacity of 300 pounds.
- Engage the brakes when sitting on the seat or when the rollator is stationary.
- Use caution on uneven surfaces or inclines.
- Regularly inspect the rollator for wear and tear, especially the wheels and braking system.
- This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

PARTS LIST

Your McKesson Rollator Walker includes the following main components:

- Aluminum Frame
- Adjustable Handles with Loop-Lock Brakes
- Padded Seat
- Hinged Padded Backrest
- 6-inch Non-Marring Caster Wheels (4)
- Under-seat Storage Pouch

SETUP AND ASSEMBLY

The McKesson Rollator Walker is designed for easy setup. Follow these steps:

1. **Unfold the Rollator:** Carefully unfold the rollator frame until it locks into place. Ensure all locking mechanisms are engaged.
2. **Adjust Handle Height:** Locate the adjustment levers on the handles. Press the lever and slide the handle up or down to your desired height. Ensure the handles are at a comfortable height, allowing for a slight bend in your elbows when standing upright. Release the lever to lock the handles in position.
3. **Position the Backrest:** The hinged padded backrest can be folded up, down, or removed as needed for comfort or storage.
4. **Check Brakes:** Ensure the deluxe loop-lock braking system is functioning correctly. Squeeze the hand brakes to slow down or stop. Push the handle into the down position to lock the brakes for a stationary position.
5. **Verify Stability:** Before initial use, gently push down on the handles and seat to confirm the rollator is stable and all components are securely attached.



Image: Fully assembled McKesson Rollator Walker, showcasing its design and components.

OPERATING INSTRUCTIONS

Walking with the Rollator

1. Stand upright behind the rollator, holding both handles firmly.
2. Push the rollator forward a short distance.
3. Step forward with one foot, then the other, maintaining a comfortable pace.
4. Keep your body close to the rollator for optimal balance and support.

Provides support and stability



Image: The rollator in use, demonstrating its role in providing support and stability.

Using the Brakes

- **To Slow Down or Stop:** Squeeze the hand brake levers located on the handles.
- **To Lock the Brakes:** Push the hand brake levers down until they click into the locked position. This will prevent the rollator from moving. Always lock the brakes before sitting on the seat.
- **To Release the Brakes:** Pull the hand brake levers up to unlock the wheels.

Deluxe braking system



Image: Detail of the deluxe loop-lock braking system for secure stopping and parking.

Using the Seat and Storage

- **Sitting:** Ensure the rollator is on a flat, stable surface. Lock both brakes before attempting to sit. Sit slowly and carefully on the padded seat.
- **Storage:** The under-seat pouch provides convenient storage for personal items. Lift the seat to access the pouch.

Convenient storage and seating



Image: The rollator's seat and accessible storage pouch for personal belongings.

Folding for Storage and Transport

The rollator folds down for easy storage and transport. Consult the product's quick-start guide or manufacturer's instructions for specific folding mechanisms, as these can vary by model.

MAINTENANCE

- **Cleaning:** Wipe down the frame and seat with a damp cloth and mild soap. Do not use abrasive cleaners.
- **Wheel Inspection:** Periodically check the 6-inch non-marring caster wheels for debris or damage. Ensure they swivel freely.
- **Brake Adjustment:** If the brakes feel loose or do not engage properly, they may need adjustment. Consult a qualified technician or the manufacturer's support for assistance.
- **Fasteners:** Ensure all screws, bolts, and other fasteners are tight. Do not overtighten.

TROUBLESHOOTING

- **Rollator pulls to one side:** Check for uneven handle height adjustment or debris in the wheels. Ensure wheels are properly aligned.

- **Brakes not engaging:** Inspect the brake cables for damage or looseness. The brake mechanism may require adjustment by a professional.
- **Difficulty folding:** Ensure the seat is fully lifted and any locking pins are disengaged before attempting to fold.

SPECIFICATIONS

Feature	Detail
Brand	McKesson
Model Number	B0785N3TYX
Color	Red
Material	Aluminum, Metal
Weight Capacity	300 lbs
Item Weight	15.4 Pounds
Product Dimensions	20"D x 24.5"W x 37"H
Handle Height Adjustment	32 to 37 inches
Wheel Size	6-inch non-marring casters



Image: Detailed dimensions of the rollator for reference.

OFFICIAL PRODUCT VIDEO

Your browser does not support the video tag.

Video: An overview of the McKesson Rollator Walker, demonstrating its features and ease of use.

WARRANTY AND SUPPORT

For warranty information, product support, or to purchase replacement parts, please contact McKesson customer service directly. Refer to the packaging or the official McKesson website for the most current contact details.

© 2023 McKesson. All rights reserved.

	<p>McKesson Lightweight Rollator User Manual and Assembly Instructions</p> <p>Comprehensive guide for the McKesson Lightweight Rollator, including assembly instructions, safety precautions, loop lock operation, and hand brake system adjustment. Features model numbers 146-R726SL, 146-R726GR, 146-R726BK, 146-R726BL, and 126-R726RD.</p>
	<p>McKesson ErgoCart Cleaning and Disinfection Instructions</p> <p>Detailed cleaning and disinfection instructions for McKesson ErgoCart storage carts and mobile workstations, including recommended cleaning agents and processes.</p>
	<p>McKesson Refrigerator/Freezer Thermometer TRIPLE DISPLAY SINGLE PROBE MCK821RFV</p> <p>Detailed information and operating instructions for the McKesson Refrigerator/Freezer Thermometer, Triple Display, Single Probe (Model MCK821RFV). Features include high accuracy, min/max memory, alarm settings, and NIST traceability.</p>
	<p>McKesson Disposable Vaginal Specula: Instructions for Use and Product Information</p> <p>Comprehensive guide to McKesson's disposable vaginal specula, covering usage instructions, cautions, indications, sterilization, and product specifications for medical professionals.</p>
	<p>McKesson Wireless Refrigerator/Freezer Thermometer MCK895WRFT - Triple Display, Single Probe</p> <p>User manual and specifications for the McKesson Wireless Refrigerator/Freezer Thermometer, model MCK895WRFT. Details features like triple display, single probe, large digit, wireless connectivity, alarm settings, and error codes.</p>
	<p>McKesson Refrigerator/Freezer Thermometer TRIPLE DISPLAY DUAL PROBE MCK821RFV2 User Manual</p> <p>User manual and specifications for the McKesson Refrigerator/Freezer Thermometer TRIPLE DISPLAY DUAL PROBE (Model MCK821RFV2). Features include high accuracy, min/max memory, low/high alarm settings, and NIST traceability.</p>