

CANMOV Power Lift

CANMOV Power Lift Electric Recliner Chair User Manual

Model: Power Lift

IMPORTANT SAFETY INSTRUCTIONS

Please read and understand all instructions before using this product. Failure to follow these instructions may result in injury or product damage.

- Consult your doctor before using this product.
- Discontinue use and consult your doctor if you experience persistent pain in a muscle or joint after prolonged use, as this could indicate a more serious condition.
- Ensure all connections are secure and the chair is properly assembled before operation.
- Keep children and pets away from the moving parts of the chair.
- Do not exceed the maximum weight capacity of 300 pounds.

PACKAGE CONTENTS

Your CANMOV Power Lift Electric Recliner Chair package includes the following components:

- Main Chair Unit (with integrated power mechanism)
- Two Armrests
- Backrest
- Power Supply Unit and Cables
- Wired Remote Control for Lift/Recline
- Wired Remote Control for Massage and Heating Functions (if applicable)
- Instruction Manual

SETUP AND ASSEMBLY

Assembly of your CANMOV Power Lift Electric Recliner Chair is designed to be straightforward and typically takes less than 10 minutes. No additional tools are required.

1. **Unpack Components:** Carefully remove all components from their packaging. The chair typically arrives in three main boxes: one for the main unit, one for the two armrests, and one for the backrest.
2. **Locate Power Supplies:** The power supply unit and cables are often found in a bag labeled "spare parts." Ensure these are not discarded.
3. **Attach Armrests:** Slide the armrests onto the main chair unit until they securely lock into place.
4. **Attach Backrest:** Slide the backrest onto the main chair unit until it securely lock into place.
5. **Connect Power Cables:** Connect the power cables from the chair to the power supply unit. Refer to the included diagram in the printed manual for specific connection points.
6. **Plug into Power Outlet:** Plug the power supply unit into a standard electrical outlet.

PRODUCT DETAILS



Side Pockets



Sturdy Metal
Frame & Motor



Figure 1: Illustration of the sturdy metal frame and motor mechanism, highlighting the robust construction of the chair.

Your browser does not support the video tag.

Video 1: Demonstrates the assembly process and basic functionality of the CANMOV Power Lift Recliner Chair.

OPERATING INSTRUCTIONS

Your CANMOV Power Lift Electric Recliner Chair offers multiple positions and features for enhanced comfort and ease of use.

Power Lift and Recline Functions

The chair's height and recline positions are controlled by a side button located on the armrest.

- **To Recline:** Press and hold the appropriate button to smoothly adjust the chair to your desired reclining position. The chair can recline to various angles suitable for reading, watching TV, or napping.
- **To Lift/Stand:** Press and hold the lift button to gently raise the chair, assisting you in standing up from a seated position.



Figure 2: Close-up of the two built-in buttons on the chair's side, controlling the lift and recline functions, and an integrated USB port.



Figure 3: Visual guide demonstrating the chair's versatility for different activities and assistance with standing.

Massage and Heating Functions (if applicable)

If your model includes massage and heating features, these are controlled by a separate remote. Please refer to the specific remote control instructions provided in your package for detailed operation.

Your browser does not support the video tag.

Video 2: Shows the power lift and recline functions of the CANMOV Oversized Wide Seat Lift Recliner Chair.

MAINTENANCE

To ensure the longevity and appearance of your CANMOV Power Lift Electric Recliner Chair, follow these simple maintenance guidelines:

- **Cleaning:** Wipe clean with a soft, damp cloth. Avoid harsh chemicals or abrasive cleaners that may damage the upholstery.
- **Regular Inspection:** Periodically check all connections and moving parts for any signs of wear or looseness. Tighten as necessary.
- **Spill Management:** Blot spills immediately with a clean, dry cloth. Do not rub.

TROUBLESHOOTING

If you encounter any issues with your CANMOV Power Lift Electric Recliner Chair, please refer to the following common troubleshooting steps:

- **Chair Not Responding:**
 - Ensure the power cord is securely plugged into both the chair and a functional wall outlet.
 - Check all cable connections underneath the chair to ensure they are firmly attached.
 - Verify that the power supply unit is receiving power (some units may have an indicator light).
- **Remote Control Not Working:**
 - Confirm the remote control cable is securely connected to the chair.
 - If the remote has batteries (for massage/heat functions), ensure they are correctly inserted and not depleted.
- **Unusual Noises During Operation:**
 - Check for any obstructions around the chair's moving parts or mechanism.

- Ensure the chair is on a level surface.

If the issue persists after attempting these steps, please contact CANMOV customer service for further assistance.

SPECIFICATIONS

Feature	Detail
Brand	CANMOV
Model Name	Power Lift
Color	Chocolate
Product Dimensions	39"D x 37"W x 40"H
Seat Height	20.4 Inches
Seat Depth	22 inches
Maximum Weight Recommendation	300 Pounds
Upholstery Fabric Type	Suede
Special Features	Reclining, Cushion Availability, Durable, Arm Rest, Upholstered, Side Pocket
Reclining Position Count	3
Product Care Instructions	Wipe Clean

Product Size



Figure 4: Detailed product dimensions for the CANMOV Power Lift Electric Recliner Chair.

WARRANTY AND SUPPORT

CANMOV is committed to providing high-quality products and excellent customer service.

- **Warranty:** The product typically comes with a manufacturer's warranty. Please refer to the warranty information included in your product packaging for specific terms and duration.
- **Customer Service:** For any questions, concerns, or assistance with your CANMOV Power Lift Electric Recliner Chair, please contact CANMOV customer service. We aim to respond to all inquiries within 24 hours.
- **Protection Plans:** Extended protection plans may be available for purchase separately to cover your product beyond the manufacturer's warranty.