

[Manuals.plus](#) /

> [CANMOV](#) /

> CANMOV Power Lift Recliner Chair Instruction Manual (Model: Power Lift)

## CANMOV Power Lift

# CANMOV Power Lift Recliner Chair Instruction Manual

Model: Power Lift

Brand: CANMOV

## 1. INTRODUCTION

---

This manual provides essential information for the safe and effective use of your CANMOV Power Lift Recliner Chair. This chair is designed to offer enhanced comfort and assistance, particularly for individuals who may benefit from a powered lift and recline mechanism. Its robust construction and overstuffed design ensure a comfortable and supportive experience.



*Image 1: The CANMOV Power Lift Recliner Chair in its upright position.*

## 2. SAFETY INFORMATION

---

- **Consult a Doctor:** For your health, consult your doctor before using this product.
- **Discontinue Use if Pain Occurs:** If you experience pain in a muscle or joint after using the chair for a prolonged period, discontinue use and consult your doctor. Persistent pain could indicate a more serious condition.
- **Power Source:** Ensure the chair is connected to a stable power source. Avoid overloading electrical outlets.
- **Clear Surroundings:** Keep the area around the chair clear of obstructions to allow for safe reclining and lifting movements.
- **Children and Pets:** Keep children and pets away from the moving parts of the chair to prevent injury.
- **Weight Limit:** Do not exceed the maximum weight capacity of 300 pounds.

## 3. PACKAGE CONTENTS

---

Upon receiving your CANMOV Power Lift Recliner Chair, please verify that all components are present and

undamaged. The chair typically arrives in multiple boxes for easier transport and assembly.

- Main Chair Unit (includes base and motor)
- Backrest
- Two (2) Armrests
- Power Supply Unit and Cables
- Remote Control for Lift/Recline
- Remote Control for Massage/Heat (if applicable to your model)
- Instruction Manual

**Note:** The power supply and cables are often found in a bag labeled 'Spare Parts' within the packaging. Do not discard this bag, as these components are essential for the chair's operation.

## 4. SETUP AND ASSEMBLY

---

Assembly of the CANMOV Power Lift Recliner Chair is designed to be straightforward and typically takes less than 10 minutes, requiring no additional tools beyond what is provided or your hands. Follow these general steps:

1. **Unpack Components:** Carefully remove all parts from their packaging.
2. **Attach Armrests:** Slide the armrests onto the main chair unit until they securely lock into place.
3. **Attach Backrest:** Align the backrest with the main chair unit and slide it down until it clicks and locks securely.
4. **Connect Power Cables:** Locate the power supply unit and all necessary cables. Connect them according to the instructions provided in the included manual. Ensure all connections are firm.
5. **Plug In:** Plug the power supply into a suitable electrical outlet.



Image 2: Illustrates the simple process of attaching the backrest to the chair base.

## Assembly Demonstration Video

Video 1: A demonstration of the assembly process and basic operation of the CANMOV Power Lift Recliner Chair.

## 5. OPERATING INSTRUCTIONS

Your CANMOV Power Lift Recliner Chair features easy-to-use remote controls for its various functions.

### 5.1. Lift and Recline Functions

The primary remote control manages the chair's lift and recline movements. It typically features two main buttons:

- **Up Button:** Press and hold to gently lift the chair, assisting the user to a standing position.
- **Down Button:** Press and hold to recline the chair backward and raise the footrest. Continue holding to return the chair to an upright sitting position.

## 2 Built-in Buttons

Adjust the angle you want at the click of a button



Image 3: Illustrates the two built-in buttons for adjusting the chair's angle, showing positions for reading, watching TV, napping, and standing up.

### 5.2. Massage and Heat Functions (if applicable)

Some models include a separate remote for massage and heat functions. Refer to your specific model's manual for detailed instructions on these features.

### Operation Demonstration Video

Video 2: Demonstrates the lift and recline functions of the CANMOV Power Lift Recliner Chair, highlighting its use for various users.

## 6. MAINTENANCE

- **Cleaning:** To clean the upholstery, simply wipe it clean with a damp cloth. Avoid harsh chemicals or abrasive cleaners.
- **Regular Checks:** Periodically inspect all connections and moving parts to ensure they are secure and functioning correctly.

## 7. TROUBLESHOOTING

- **Chair Not Responding:**

- a. Ensure the power cord is securely plugged into both the chair and the wall outlet.
- b. Check that all internal connections (found in the 'Spare Parts' bag during assembly) are properly connected.
- c. Verify that the wall outlet is functional by plugging in another device.

- **Unusual Noises:** If the chair makes unusual noises during operation, discontinue use and inspect the mechanism for any obstructions or loose parts. Contact customer service if the issue persists.

## 8. SPECIFICATIONS

Feature	Specification
Brand	CANMOV
Model Name	Power Lift
Color	Light Tan
Product Dimensions (D x W x H)	39"D x 39"W x 42"H
Seat Height	20.4 Inches
Seat Depth	22.5 inches
Maximum Weight Recommendation	300 Pounds
Upholstery Fabric Type	Suede
Fill Material	Foam
Special Features	Reclining, Cushion Availability, Durable, Arm Rest, Side Pocket
Item Weight	109 pounds
UPC	651354988065

# DIMENSIONS



Image 4: Detailed dimensions of the recliner chair.



*Image 5: Highlights the convenient side pocket and the robust metal frame and motor mechanism.*

## 9. WARRANTY AND SUPPORT

CANMOV is committed to providing high-quality products and customer satisfaction. For any inquiries, assistance, or to report issues, please contact CANMOV customer service. They are available to provide helpful, prompt, and friendly support.

Additional protection plans may be available for purchase, offering extended coverage beyond the manufacturer's warranty. Please refer to your purchase documentation for details on available protection plans and their terms.