



Manuals.plus /

› TRG /

› TRG CoFlex TLC Calamine 2-Layer Compression Bandage System User Manual

TRG CoFlex TLC Calamine 2-Layer Compression Bandage System

TRG CoFlex TLC Calamine 2-Layer Compression Bandage System User Manual

MODEL: CoFLEX TLC CALAMINE 2-LAYER COMPRESSION BANDAGE SYSTEM

Introduction

The TRG CoFlex TLC Calamine 2-Layer Compression Bandage System is designed to provide sustained compression and support for various medical applications. This system consists of two distinct layers: a medicated foam comfort roll impregnated with calamine and a non-latex moisture-resistant cohesive compression bandage. It is engineered for ease of application and effective performance, delivering 25-30 mmHg compression.

Product Components

- **Layer 1: Medicated Foam Comfort Roll**
This layer is impregnated with calamine, measuring 4 inches x 6 yards. It provides comfort and soothing properties to the skin.
- **Layer 2: Cohesive Compression Bandage**
This non-latex, moisture-resistant bandage measures 4 inches x 7 yards. It adheres directly to the first layer, providing the necessary compression and securing the system in place.
- **Nylon Stocking**
A nylon stocking is included to be worn over the entire system, facilitating easy movement and protecting the bandage.

Setup and Application

Follow these steps for proper application of the CoFlex TLC Calamine 2-Layer Compression Bandage System:

1. **Unpack the System:** Carefully remove all components from the packaging.
2. **Apply Layer 1 (Calamine Foam Roll):** Begin by applying the medicated foam comfort roll (Layer 1) to the affected area. Ensure it is applied smoothly and evenly, covering the intended region. This layer is 4 inches x 6 yards.
3. **Apply Layer 2 (Cohesive Compression Bandage):** Directly over Layer 1, apply the non-latex moisture-resistant cohesive compression bandage (Layer 2). This bandage is designed to stick securely to the first layer. Apply with consistent tension to achieve the desired 25-30 mmHg compression. This layer is 4 inches x 7 yards. The system is designed to apply cleanly, similar to a dry foam first layer, without the mess associated with traditional Unna boots.
4. **Wear Nylon Stocking:** Once both layers are applied, gently pull the provided nylon stocking over the entire system. This helps to keep the bandage in place and allows for easier movement.



This image displays the packaging and the two distinct layers of the CoFlex TLC Calamine Compression Bandage System: the medicated foam comfort roll and the cohesive compression bandage.

Maintenance and Care

Proper care of the compression bandage system is crucial for effectiveness and hygiene:

- **Consult a Healthcare Professional:** Always follow the specific instructions provided by your healthcare professional regarding the duration of wear, frequency of changes, and any specific care requirements.
- **Keep Dry:** Keep the bandage system dry. If it becomes wet, consult your healthcare provider for guidance on replacement or drying procedures.
- **Do Not Reuse:** This product is designed for single use or for use over a specified period as directed by a healthcare professional. Do not wash or attempt to reuse the bandages.
- **Disposal:** Dispose of used bandages according to local regulations for medical waste.

Troubleshooting

If you experience any issues while using the compression bandage system, consider the following:

- **Discomfort or Pain:** If you experience increased pain, numbness, tingling, or excessive tightness, immediately loosen or remove the bandage and consult your healthcare professional.
- **Skin Irritation:** Should you notice any skin redness, itching, rash, or blistering under the bandage, remove it and seek medical advice.
- **Bandage Slippage:** Ensure the bandage is applied with consistent tension and that the cohesive

layer is properly adhered. If slippage persists, reapply or consult a professional.

- **Swelling Increase:** If swelling increases or new swelling develops, contact your healthcare provider immediately.

Specifications

Feature	Detail
Product Name	CoFlex TLC Calamine 2-Layer Compression Bandage System
Brand	TRG
Compression Level	25-30 mmHg
Layer 1 Material	Medicated foam comfort roll impregnated with calamine
Layer 1 Dimensions	4 inches x 6 yards
Layer 2 Material	Non-latex moisture-resistant cohesive compression bandage
Layer 2 Dimensions	4 inches x 7 yards
Sterility	Non-Sterile
Product Style	Cohesive, Self-Adherent
Package Dimensions	6.54 x 4.96 x 3.9 inches
Item Weight	0.38 Kilograms (13.44 ounces)

Safety Information

Please read and adhere to the following safety guidelines:

- **Consult a Healthcare Professional:** Always consult with a doctor or qualified healthcare professional before using this product, especially if you have pre-existing medical conditions, circulatory issues, or open wounds.
- **External Use Only:** This product is for external application only.
- **Not for Open Wounds:** Do not apply directly to open wounds or broken skin without specific medical advice and appropriate wound dressings underneath.
- **Allergies:** If you have known allergies to latex or any components of the bandage, do not use this product. The cohesive bandage is non-latex, but always check for other sensitivities.
- **Monitor Skin:** Regularly inspect the skin under and around the bandage for any signs of irritation, discoloration, swelling, or discomfort.
- **Keep Out of Reach of Children:** Store this product safely away from children.

Warranty and Support

Specific warranty information for the TRG CoFlex TLC Calamine 2-Layer Compression Bandage System is not provided in this manual. For details regarding product warranty, returns, or technical support, please contact the manufacturer, TRG, directly through their official customer service channels. Keep your purchase receipt for any warranty claims.

