



[Manuals.plus](#) /

› [Cleanse Right](#) /

› Cleanse Right Electronic Ear Wax Removal Device Instruction Manual

## Cleanse Right Electronic Ear Wax Removal Device

# Cleanse Right Electronic Ear Wax Removal Device Instruction Manual

Model: Electronic Ear Wax Removal Device

## INTRODUCTION

This manual provides detailed instructions for the safe and effective use of your Cleanse Right Electronic Ear Wax Removal Device. This device is designed to assist in the removal of earwax blockages, which can contribute to hearing loss, earaches, tinnitus, and dizziness.



Image: The Cleanse Right Electronic Ear Wax Removal Device kit, including the main unit, ear drops, and reusable tips.

## IMPORTANT SAFETY INFORMATION

### WARNING:

- Always consult your physician before using the Cleanse Right system to ensure you are healthy enough for this product.
- Do not use if you have ear infections, a perforated eardrum, ear tubes, or have recently had ear surgery.
- Avoid using cotton swabs or ear candling, as these methods can worsen impaction and cause long-term damage.
- Ensure the device tips are securely attached before use to prevent detachment inside the ear canal.

For additional safety information regarding earwax removal tips, please watch the following video:

Your browser does not support the video tag.

Video: Important information regarding the safety of earwax removal tips, highlighting potential dangers of improper tools.

## WHAT'S INCLUDED

---

Your Cleanse Right Electronic Ear Wax Removal Device kit includes the following components:

- Electronic Ear Wax Removal Device (main unit)
- Charging cord
- USA Made Ear Drops (15 mL)
- USA Made, Reusable, Dishwasher Friendly Tips with Splash Dish (3 units)
- Detailed Instructions for use (this manual)

For a visual guide on unboxing and identifying components, please refer to the video below:

Your browser does not support the video tag.

Video: An overview of the Cleanse Right Electronic Ear Wax Remover kit contents and initial setup.

---

## SETUP AND PREPARATION

### 1. Charge the Device

Before first use, ensure the device is fully charged. Connect the charging cord to the device and plug the USB adapter into a wall outlet. Allow at least 10 minutes for initial charging.

### 2. Fill the Water Basin

Open the water basin cap at the bottom of the device and fill it with warm water. Ensure the cap is securely closed until it snaps into place.

### 3. Attach the Neck and Tip

Push the wider end of the neck component into the top of the device body. It should fit snugly but allow rotation. Then, gently attach one of the reusable tips to the top of the neck. The integrated splash dish on the tip prevents over-insertion into the ear canal.

# Ear Wax Removal Kit Assembly Instructions



Gently squeeze 5-10 drops into the ear to soften the ear wax. Wait 3-5 minutes.



Ensure the battery is charged. Attach a clean, reusable tip to the Cleanse Right Electronic Ear Wax Remover and fill the water basin at the base of the device with water.



Insert the clean, reusable tip into the ear and spray for 5-10 seconds to remove the blockage. NOTE: The Cleanse Right Electronic Ear Wax Remover can be used in the shower as it's waterproof.



If blockage is not dislodged after first spraying, do another 5-10 seconds or until blockage is removed.



When finished, clean the reusable tip by putting into the dishwasher or soaking in warm, soapy water for 5 minutes and then let air dry to be used again.

Image: Visual guide for assembling the Cleanse Right Electronic Ear Wax Removal Device, showing how to attach the tip and fill the water basin.

## OPERATING INSTRUCTIONS

### 1. Pre-treatment with Ear Drops

For optimal results, pre-treat the affected ear with the Cleanse Right USA Made Ear Drops. Gently squeeze 5-10 drops into the ear while lying on your side. Allow the drops to settle for 3-5 minutes. You may feel a fizzing sensation, which indicates the drops are working to soften the earwax. After the waiting period, turn your head over a sink to allow the solution and softened wax to drain out.

Your browser does not support the video tag.

Video: Instructions on how to use Cleanse Right Ear Wax Drops for softening earwax prior to removal.

## 2. Using the Electronic Ear Wax Remover

Stand over a sink or use the device in the shower, as it is waterproof. Insert the clean, reusable tip into your ear. The integrated splash dish will prevent it from going in too deeply.

The device features 3 separate speeds and a pulse mode. To test the pressure, point the device into the sink and cycle through the settings:

1. Press the power button once for the low setting.
2. Press again for the high setting.
3. Press a third time for pulse mode.
4. Press a fourth time to turn the device off.

Once you are comfortable with the pressure, insert the tip into your ear and spray for 5-10 seconds to remove the blockage. If the blockage is not dislodged after the first spraying, repeat the 5-10 second spray until the blockage is removed.



The image shows a white, ergonomic electronic ear wax remover with a clear water reservoir at the base and a long, thin handle with a splash guard at the tip. The device is shown in a vertical orientation.

# Blue EchoCare

## 3 MODES

The Electronic Ear Wax Remover has 3 Modes, including Pulse Mode to help remove as much Ear Wax as possible

-  SAFE
-  PAINLESS
-  EFFECTIVE
-  PORTABLE



Three reusable ear wax tips are shown at the bottom of the graphic, arranged in a row and slightly overlapping. Each tip is white and has a small, circular opening at the end.

Image: The Cleanse Right Electronic Ear Wax Remover highlighting its three modes and key features.

## MAINTENANCE AND CLEANING

---

After each use, disconnect the tip from the device. The reusable tips are dishwasher friendly. Alternatively, you can wash them in warm, soapy water. Allow the tips to air dry completely before storing or reusing. Each tip can be washed up to 100 times before needing replacement.

Drain any excess water from the device's water basin after use.

## TROUBLESHOOTING

---

### Dizziness During Use

You may experience some dizziness while using the device. This is normal. If this happens, sit down and continue using the device with a Cleanse Right catch basin (not included in this kit) instead of over a sink or in the shower, to manage any discomfort.

### Stubborn Wax Blockage

If the earwax blockage is still present after initial attempts, repeat the pre-treatment with ear drops and the irrigation process as many times as necessary to clear the blockage. Ensure you are using warm water, as cold water can cause dizziness.

### Tip Detachment

Always ensure the reusable tip is securely attached to the device's neck before use. The patent-pending design with a splash dish is intended to prevent over-insertion and secure attachment.

## PRODUCT SPECIFICATIONS

---

- **Manufacturer:** Blue Echo Care
- **ASIN:** B07DHT2SS6
- **Package Dimensions:** 11.42 x 4.8 x 3.11 inches
- **Item Weight:** 4.5 ounces
- **Waterproof:** Yes
- **Speeds:** 3 (Low, High, Pulse)
- **Tips:** Reusable, Dishwasher Friendly, USA Made
- **Ear Drops:** USA Made, Carbamide Peroxide 6.5%

## WARRANTY AND CUSTOMER SUPPORT

---

Cleanse Right is a brand of Blue Echo Care, an American company committed to customer satisfaction. If you have any questions, concerns, or require assistance with your Cleanse Right Electronic Ear Wax Removal Device, please reach out to our American Customer Service Team. We are available 24/7 to provide fast and responsive support.

For contact information, please refer to the packaging or visit the official Cleanse Right website.