

[Manuals.plus](#) /

› [SOUNDLINK](#) /

› SOUNDLINK Hearing Aid Drying Capsule Desiccant User Manual

SOUNDLINK Drying Capsule Desiccant



Hearing Aid Drying Capsule Desiccant User Manual

Brand: SOUNDLINK | Model: Drying Capsule Desiccant

INTRODUCTION

The SOUNDLINK Hearing Aid Drying Capsule Desiccant is designed to absorb daily moisture from hearing aids and earmolds, helping to prolong their life and maintain optimal performance. Regular use of a drying capsule is crucial for the care of your hearing devices.

Your browser does not support the video tag.

Video: A 360-degree view of the drying capsules, providing a comprehensive visual overview of the product packaging and individual capsules.

PRODUCT FEATURES

- **Color:** Orange
- **Weight:** 30 g per capsule
- **Dimensions:** Approximately 63mm diameter x 19mm height
- Designed for use with a drying jar or case for hearing aids and earmolds.
- Effectively absorbs daily moisture.
- Helps to prolong the lifespan of hearing aids.

Image: Dimensions of the SOUNDLINK drying capsule, showing its diameter (64mm) and height (19mm), and weight (30g).

SETUP

Before first use, prepare the drying capsule by removing its protective seal.

1. Carefully tear apart the aluminum foil paper from the capsule. This exposes the desiccant beads.
2. Ensure the capsule is intact and the desiccant beads are visible. The beads are typically orange when fresh.



Image: A hand peeling back the aluminum foil seal from a SOUNDLINK drying capsule, preparing it for use.

Your browser does not support the video tag.

Video: This video demonstrates how to unpackage the drying capsules from their blister packs, showing the individual capsules ready for use.

OPERATING INSTRUCTIONS

To effectively dry your hearing aids, follow these steps:

1. Place the prepared drying capsule into a suitable drying jar or case designed for hearing aids.
2. Ensure your hearing aids are clean and free from debris before placing them into the drying jar. Open the battery doors of your hearing aids to allow for better air circulation and moisture absorption.
3. Place your hearing aids into the drying jar, ensuring they do not directly touch the desiccant capsule. Close the lid of the drying jar securely.
4. Leave your hearing aids in the drying jar overnight, or for at least 6-8 hours, to allow the desiccant to absorb moisture.

Your browser does not support the video tag.

Video: This video demonstrates the process of placing a 3g drying capsule into a hearing aid drying box, showing how it fits and is secured.

MAINTENANCE

The drying capsule contains desiccant beads that change color as they absorb moisture, indicating when replacement is needed. When the orange beads turn a pale yellow or white, the capsule has absorbed its maximum moisture and should be replaced with a new one.



Image: A top-down view of a fresh, orange SOUNDLINK drying capsule, illustrating the initial color of the desiccant beads.

TROUBLESHOOTING

- **Hearing aids still damp:** Ensure the drying capsule is fresh (orange color) and the drying jar lid is securely closed. Also, verify that the hearing aid battery doors are open during the drying process.
- **Capsule not changing color:** If the capsule remains orange after prolonged use, it might not be exposed to sufficient moisture, or the drying jar might not be sealed properly.

SPECIFICATIONS

Brand	SOUNDLINK
Model	Drying Capsule Desiccant

Color	Orange
Weight	30 g/pc
Dimensions	63mm (diameter) x 19mm (height)
ASIN	B07GBS42KR
Date First Available	August 10, 2018

IMPORTANT SAFETY INFORMATION

WARNING: DESICCANT. DO NOT EAT. Keep out of reach of children and pets. If ingested, seek immediate medical attention. The desiccant beads are not toxic but can pose a choking hazard.

WARRANTY AND SUPPORT

For any questions, concerns, or support regarding your SOUNDLINK Hearing Aid Drying Capsule Desiccant, please contact the manufacturer directly through their official website or the retailer where the product was purchased. Please refer to your purchase documentation for specific warranty details.

© 2024 SOUNDLINK. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

