

Homecraft 49332

Homecraft Corner Shower Stool Instruction Manual

Model: 49332
Brand: Homecraft

PRODUCT OVERVIEW

The Homecraft Corner Shower Stool is designed to provide a comfortable and safe seating solution for individuals requiring assistance during showering or bathing. Its unique quadrant design allows it to fit snugly into shower and bath corners, maximizing space. This adjustable stool is suitable for elderly, handicapped, and disabled users, as well as those recovering from injuries, offering stability and support.



Image: The Homecraft Corner Shower Stool, showcasing its quadrant seat and adjustable legs with non-slip feet.

SAFETY INFORMATION

Please read and understand all safety instructions before using the Homecraft Corner Shower Stool. Failure to do so may result in injury.

- **Weight Capacity:** Do not exceed the maximum weight capacity of 265 pounds (120 kg).
- **Stability Check:** Before each use, ensure all three legs are firmly on the ground and the stool is stable. Do not use if the stool wobbles or feels unstable.
- **Leg Adjustment:** Ensure the push-pins are fully engaged in the adjustment holes after changing leg height. Uneven leg heights can cause instability.
- **Non-Slip Feet:** Verify that the rubber feet are clean and free of soap residue or debris to maintain maximum slip resistance.
- **Surface:** Use the stool only on flat, even surfaces within the shower or bath.

- **Inspection:** Regularly inspect the stool for any signs of wear, damage, or corrosion. Discontinue use immediately if any damage is found.
- **Children:** This product is not a toy. Keep out of reach of children when not in use by an adult requiring assistance.

PACKAGE CONTENTS

Upon opening the package, please verify that all components are present and undamaged:

- 1 x Homecraft Corner Shower Stool (fully assembled)

Note: The Homecraft Corner Shower Stool typically comes pre-assembled, requiring only height adjustment.

SETUP

Assembly and Height Adjustment

The Homecraft Corner Shower Stool is designed for minimal assembly, primarily focusing on adjusting the leg height to suit the user's needs.

1. **Unpack:** Carefully remove the shower stool from its packaging.
2. **Determine Desired Height:** Sit on a stable chair of a comfortable height. Measure the distance from the floor to the back of your knee. This measurement can serve as a guide for the stool's seat height. The stool is adjustable from 17 inches to 21 inches.
3. **Adjust Legs:**
 - Locate the push-pins on each leg of the stool.
 - Press the push-pin inward and slide the lower leg section up or down to the desired height.
 - Ensure the push-pin fully snaps into one of the adjustment holes.
 - Repeat this process for all three legs, ensuring they are adjusted to the *same height setting* to maintain stability.
4. **Verify Stability:** Before placing the stool in the shower, gently press down on the seat and check for any wobbling. All three legs must be evenly set and firmly on the ground.



Image: The Homecraft Corner Shower Stool demonstrating its adjustable leg mechanism, allowing users to customize the height for optimal comfort and safety.



Image: A close-up view of the push-pin mechanism on the stool's leg, which is used for secure height adjustment.

OPERATING INSTRUCTIONS

Using your Homecraft Corner Shower Stool safely and effectively:

1. **Placement:** Position the stool firmly in a corner of your shower or bath. Ensure all three non-slip rubber feet are in full contact with the shower floor.
2. **Entering the Shower:** Carefully enter the shower area. Use grab bars if available for additional support.
3. **Sitting Down:** Turn around and back up to the stool until your legs touch the seat. Slowly lower yourself onto the seat, using your arms for support if needed. Ensure you are centered and balanced on the seat.
4. **During Showering:** Remain seated while showering. Avoid leaning excessively or standing on the stool.
5. **Exiting the Shower:** Carefully push up from the seat, using your legs and arms for support. Once standing, carefully exit the shower area.

Always exercise caution when entering and exiting the shower, especially on wet surfaces.

MAINTENANCE

Proper maintenance will extend the life of your Homecraft Corner Shower Stool and ensure continued safe operation.

- **Cleaning:** Clean the stool regularly with a mild soap and water solution. Rinse thoroughly and wipe dry to prevent water spots and soap scum buildup. Avoid abrasive cleaners or harsh chemicals, as these can damage the vinyl seat and frame finish.
- **Inspection:** Periodically inspect the stool for:
 - Loose or damaged push-pins on the legs.
 - Cracks, tears, or excessive wear on the vinyl seat pad.
 - Damage or hardening of the rubber non-slip feet.
 - Signs of rust or corrosion on the frame.
- **Storage:** When not in use for extended periods, store the stool in a dry, cool place away from direct sunlight.

If any part of the stool appears damaged or worn, discontinue use immediately and contact customer support for assistance.

TROUBLESHOOTING

Below are common issues and their potential solutions:

Problem	Possible Cause	Solution
Stool wobbles or is unstable.	Legs are not adjusted to the same height. Push-pins are not fully engaged. Uneven shower floor surface. Damaged rubber feet.	Ensure all three legs are set to the identical height setting. Press push-pins firmly until they snap into place. Reposition the stool on a flatter section of the shower floor. Inspect feet for damage; replace if necessary.
Difficulty adjusting leg height.	Push-pin is stuck or corroded. Debris inside the leg tubing.	Apply a small amount of silicone spray to the push-pin and gently work it. Clean out any debris from the leg tubing.
Stool slips on shower floor.	Rubber feet are dirty or have soap residue. Rubber feet are worn or hardened.	Clean the rubber feet thoroughly with soap and water. Replace worn or hardened rubber feet.

If the problem persists after attempting these solutions, please contact Homecraft customer support.

SPECIFICATIONS

Feature	Detail
Model Number	49332

Feature	Detail
Brand	Homecraft
Seat Dimensions (W x D)	19 inches x 15 inches (48.3 cm x 38.1 cm)
Adjustable Height Range	17 inches - 21 inches (43.2 cm - 53.3 cm)
Weight Capacity	265 pounds (120 kg)
Seat Material	Vinyl covered pad
Legs/Frame Material	Coated metal (e.g., aluminum or steel) with rubber non-slip feet

WARRANTY INFORMATION

Homecraft products are manufactured to high standards and are subject to strict quality control. For specific warranty details regarding your Homecraft Corner Shower Stool, please refer to the warranty card included with your purchase or visit the official Homecraft website. Generally, warranty covers manufacturing defects under normal use.

Please retain your proof of purchase for any warranty claims.

CUSTOMER SUPPORT

If you have any questions, require assistance with assembly, or need to troubleshoot an issue not covered in this manual, please contact Homecraft customer support.

- **Online Support:** Visit the official Homecraft website for FAQs, product registration, and contact forms.
- **Contact Information:** Refer to the packaging or the Homecraft website for the most current customer service phone numbers and email addresses.

When contacting support, please have your model number (49332) and purchase date ready.

© 2024 Homecraft. All rights reserved.

This manual is for informational purposes only and is subject to change without notice.