

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [FCS](#) /

> [FCS Ear Hook Instruction Manual for Siemens Signia XTM XP P6 BTE Hearing Devices](#)

FCS Signia XTM XP P6

FCS Ear Hook Instruction Manual for Siemens Signia XTM XP P6 BTE Hearing Devices

INTRODUCTION

This manual provides detailed instructions for the installation, use, and maintenance of your FCS Replacement Ear Hook for Siemens Signia XTM XP P6 BTE hearing devices. This accessory is designed to ensure a secure fit and optimal performance of your hearing aid.

PRODUCT FEATURES

- **Perfect Fit:** Designed to seamlessly integrate with your hearing amplifier or hearing aid, ensuring a secure and comfortable fit.
- **Howling Prevention:** Effectively reduces feedback and prevents annoying howling sounds, enhancing your hearing experience.
- **Easy Installation and Removal:** Allows for quick and effortless attachment and detachment, simplifying maintenance and adjustments.
- **Durable ABS Material:** Crafted from high-quality Acrylonitrile Butadiene Styrene (ABS) material for long-term reliability and durability.
- **Comfortable and Secure:** Ergonomic design provides a snug and stable fit, keeping your hearing aid in place during active movements.



Image: Key features of the FCS Ear Hook, including howling prevention and easy installation.

COMPATIBILITY

The FCS Ear Hook is compatible with the following hearing aid models:



Image: Detailed compatibility chart for various hearing aid brands and models.

Compatible Hearing Aid Models

Brand	Models
Signia	Fast P, Fun P, Fun SP, Run P, Run SP, Intuis 2 SP, Intuis 3 P/M/SP, Prompt P/SP
Rexton	Arena P1, Arena P3, Arena HP3, Targa P5, Targa HP5, Targa P5A, Targa HP5A
A&M	STF P T1, STF P T3, STF XP T3, XTM P A4, XTM XP A4, XTM P P4, XTM XP P4, XTM XP P6, XTM P P6
Audio Service	Volta P B, Volta P T, Volta HP T
Vesuvio	STF P (V-P), STF XP (V-XP)

SETUP AND INSTALLATION

Follow these steps to correctly install the FCS Ear Hook onto your hearing device:

- 1. Identify the Connection Point:** Locate the existing ear hook or the connection point on your BTE (Behind-The-Ear) hearing aid where the ear hook attaches.
- 2. Remove Old Ear Hook (if applicable):** Gently twist or pull the old ear hook to detach it from the hearing aid. Be careful not to damage the hearing aid.
- 3. Align New Ear Hook:** Take your new FCS Ear Hook and align its threaded or snap-on end with the corresponding connection point on your hearing aid.
- 4. Attach Securely:** Push or twist the FCS Ear Hook firmly until it is securely attached. Ensure there is no looseness, as this can affect sound quality and retention.
- 5. Check Fit:** Once installed, gently tug on the ear hook to confirm it is firmly in place.



Image: Visual guide demonstrating the easy installation process of the ear hook.



Image: Illustration of the ear hook correctly fitted to a hearing aid and worn on the ear.

OPERATING THE EAR HOOK

The FCS Ear Hook is a passive component designed to securely hold your BTE hearing aid in place and direct sound into the ear canal. Once installed, no further operational steps are required for the ear hook itself. Its primary function is to provide a stable and comfortable fit, preventing the hearing aid from dislodging and minimizing feedback.



Image: The ear hook providing a secure and comfortable fit for the hearing aid on the ear.

MAINTENANCE

Proper maintenance ensures the longevity and hygiene of your FCS Ear Hook:

- Regular Cleaning:** Clean the ear hook regularly with a soft, damp cloth. You may use a mild soap solution if

necessary, but ensure it is thoroughly rinsed and dried.

- **Avoid Harsh Chemicals:** Do not use abrasive cleaners, solvents, or harsh chemicals, as these can damage the ABS material.
- **Inspect for Wear:** Periodically inspect the ear hook for any signs of wear, cracks, or damage. Replace the ear hook if it shows significant signs of deterioration to maintain optimal performance and hygiene.
- **Storage:** Store the ear hook in a clean, dry place away from direct sunlight and extreme temperatures when not in use.

TROUBLESHOOTING

- **Ear Hook Feels Loose:** Ensure the ear hook is fully twisted or pushed onto the hearing aid's connection point. If it continues to feel loose, inspect both the ear hook and the hearing aid's connection for damage. A worn ear hook may need replacement.
- **Increased Howling/Feedback:** A loose or improperly seated ear hook can contribute to feedback. Re-seat the ear hook firmly. Also, ensure your ear mold or dome is properly inserted into your ear canal.
- **Discomfort:** If the ear hook causes discomfort, ensure it is correctly positioned. While the FCS Ear Hook is designed for comfort, individual ear shapes vary. If discomfort persists, consult with your hearing aid professional.

SPECIFICATIONS

- **Material:** Acrylonitrile Butadiene Styrene (ABS)
- **Color:** Transparent
- **Item Weight:** Approximately 0.705 ounces
- **Package Dimensions:** Approximately 3.15 x 1.97 x 1.18 inches
- **Manufacturer:** FCS Accessories
- **Model Number:** Signia XTM XP P6



Image: The durable ABS material used in the construction of the ear hooks.

WHAT'S IN THE BOX

- 1 x FCS Ear Hook



Image: The FCS Ear Hook as packaged.

SUPPORT

For further assistance or inquiries, please contact FCS Accessories customer support through the retailer's platform or visit the official FCS store.

